

# How to define a bio protocol in the LCD User Interface



<b>Intro message.....</b>	<b>3</b>
<b>Main page.....</b>	<b>4</b>
<b>Define Bio Protocol.....</b>	<b>5</b>
<b>Intro Denaturing.....</b>	<b>6</b>
<b>Denaturing.....</b>	<b>6</b>
<b>Annealing.....</b>	<b>7</b>
<b>Extension.....</b>	<b>7</b>
<b>Final Extension.....</b>	<b>8</b>
<b>Main Page temp, cycles, stage.....</b>	<b>9</b>

- Intro message
- Main page
- Main menu
  - Define Bio Protocol
    - § Intro Denaturing
    - § Denaturing
    - § Annealing
    - § Extension
    - § Final Extension
  - Load Bio Protocol
    - § Dog.txt
    - § Cat.txt
    - § Horse.txt
    - § Tulip.txt
  - Cooldown
  - Settings
    - § PID calibration
    - § Buzzer On/Off

This is the introduction message after the activation of buzzer (x2 beep) then the user is ready to explore the menu.

## Intro message




1<sup>st</sup> line: It appears the stage of each bio protocol (Intro Denaturing, Annealing, Extension, Final Extension)

2<sup>nd</sup> line: It appears the temperature of LID heater, left side the current temperature, right side the set temperature which in our case is 100°C

3<sup>rd</sup> line: It appears in the left side the temperature for the peltier module (current temperature and set temperature). Also the completion of cycles will demonstrate in that line on the right side.

4<sup>th</sup> line: The progress of the whole PCR experiment will appear in %

## Main page




By clicking the push button a menu will appear in the LCD. Rotating the encoder the ' > ' will move and indicate in which sub-menu will select.

## Main menu




## Define Bio Protocol




The first stage is Intro Denaturing, the user should define the Temperature, the duration in seconds (min 10 seconds, max 300 seconds. There is a choice to go 'Back' in the main menu or select 'Next' to continue with the definition of the bio protocol.

## Intro Denaturing



The second stage is Denaturing, the user should define the Temperature in Celsius, the cycles and the duration in seconds. There is a choice to go 'Back' in the main menu or select 'Next' to continue with the definition of the bio protocol.


## Denaturing






The third stage is Annealing, the user should define the Temperature in Celsius, the cycles and the duration in seconds. There is a choice to go 'Back' in the main menu or select 'Next' to continue with the definition of the bio protocol.

## Annealing



The fourth stage is Extension, the user should define the Temperature in Celsius, the cycles and the duration in seconds. There is a choice to go 'Back' in the main menu or select 'Next' to continue with the definition of the bio protocol.

## Extension




The fifth stage is Final Extension. the user should define the Temperature, the duration in seconds (min 10 seconds, max 300 seconds. There is a choice to go 'Back' in the main menu or select 'Set' to finalize the definition of the bio protocol.

## Final Extension



After successfully defining all the stages then a message 'Bio Protocol defined' appeared and the thermocycler will start to procedure all the commands.



## Main Page temp, cycles, stage

