

| Fortress Rh Phenotype +K (BLOOD GROUPING GEL CARDS) |                       |
|-----------------------------------------------------|-----------------------|
| BGGCRHPK-12                                         | BGGCRHPK-24           |
| 12 FORTRESS GEL CARDS                               | 24 FORTRESS GEL CARDS |
| STORE AT 2-8°C                                      |                       |
| FOR IN-VITRO DIAGNOSTIC USE ONLY                    |                       |

#### INTENDED USE

This test allows the detection of clinically important Rh antigens and K antigen in gel technique.

#### INTRODUCTION

All Rh antibodies should be considered to have the potential to cause hemolytic transfusion reactions and hemolytic disease of the fetus and newborn (HDFN). Anti-c is clinically the most important Rh antibody after anti-D and often causes severe HDFN, whereas anti-C, -E, and -e rarely cause HDFN, and when they do the disease is generally mild.

Antibodies of the Kell system should be considered potentially clinically significant, both from the point of view of causing severe HDFN and HTRs.

#### PRINCIPLE OF TEST

Fortress Rh Phenotype +K Test allow the detection of clinically important Rh antigens and K antigen on a single gel card. The microtubes contain gel matrix. Seven microtubes contain monoclonal antibodies and last microtubes (Ctl) don't contain antibody in the gel. Ctl microtubes are used for control of the test.

The red blood cells will agglutinate in the presence of antibody directed towards the antigen. Agglutinated cells forming a red button on the surface of the gel or agglutinates dispersed in the gel is a positive result and indicates the presence of the corresponding antibody.

#### KIT CONTENTS

12 or 24 Fortress Rh Phenotype +K Test gel card.

1 gel card contains 8x1 microtubes

|   |   |   |   |   |      |      |     |
|---|---|---|---|---|------|------|-----|
| C | c | E | e | K | DVI- | DVI+ | Ctl |
|---|---|---|---|---|------|------|-----|

Ready to use. Contains Sodium azide as preservative.

Each microtube of the card contains buffered gel medium with preservative and mixed different reagent.

- microtube C : monoclonal anti-C IgM (clone P3X25513G8)
- microtube c : monoclonal anti-c IgM (clone H48)
- microtube E : monoclonal anti-E IgM (cloneDEM1)
- microtube e : monoclonal anti-e IgM (cloneP3GD512, MS-63)
- microtube K : polyclonal anti-K
- microtube D<sup>VI</sup>- : monoclonal anti-D IgM (clone BS 225)

This anti-D monoclonal reagent detects weak D and partial variants of the D antigen. But DVI variant of the D antigen will not be detected with this monoclonal reagent.

- microtube D<sup>VI</sup>+ : monoclonal anti-D IgM (clones LDM1 / ESD1M)

This anti-D monoclonal reagent detects weak D and partial variants of the D antigen, including the DVI variant. The ESD1M antibody will directly detect RhDVI red cells

- microtube Ctl: buffered solution without antibodies (control microtube)

#### MATERIAL REQUIRED BUT NOT PROVIDED

- 10 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Test tubes.
- Fortress LISS Diluent
- Centrifuge for Blood Grouping Gel System.

#### STORAGE CONDITIONS AND STABILITY

- Store the gel cards at 2-8 °C in upright position.
- Do not freeze.
- Avoid exposure of the cards to any heat source, direct air-conditioning sources or direct sunlight.
- Do not use reagents over the expiration date.
- All the components of the kit are stable until the expiration date on the label when stored 2-8 °C.
- Do not remove foil seal until ready to use. Once opened, the gel may begin to dry out which could affect test results.
- Use the card within 2 hours after removing the aluminum foil. Otherwise, don't use the card.
- Stability of the opened gel card is 2 hours.

#### PRECAUTIONS AND WARNINGS

- This reagent kit is for in vitro diagnosis only.
- This reagent kit is for professional use only.
- All human sourced material ( monoclonal antibodies) are tested free from HBsAg, Anti-HCV and Anti-HIV. However, all materials of human origin should be considered as potentially infectious and it is recommended that all kit materials and test specimens to be handled using established good laboratory practice.
- All kit components and specimens should be regarded as potential hazards to health. It should be used and discarded according to your own laboratory's safety procedures.

#### SAMPLE COLLECTION

Fresh red blood cells are preferred for testing. Use the red blood cells collected with anticoagulants for determination of the antigens of the ABO/Rh system. If necessary, samples stored at 2-8 °C can be used up to 48 hours. Red blood cells collected in ACD, CPD, SAGM and PAGGSM can also be used until the expiry date indicated on the label of the bag at 2-8 °C.

Fibrin residues may interfere with the reaction pattern.

#### PREPARATION OF BLOOD SAMPLE

##### 5 % red blood cell suspension :

Do not use haemolysed, cloudy or contaminated samples or those containing clots

- 1- Bring Fortress LISS Diluent to reach room temperature before use
- 2- Dispense 0,5 mL of Fortress LISS Diluent into a clean tube.
- 3- Add 50 µL of whole blood or 25 µL of packed cells to the diluent and mix gently.

##### 0.8 % red blood cell suspension :

Do not use haemolysed, cloudy or contaminated samples or those containing clots.


- 1- Bring Fortress Diluent 1 to reach room temperature before use
- 2- Dispense 1 mL of Fortress Diluent 1 into a clean tube.
- 3- Add 20 µL of whole blood or 10 µL of packed cells to the diluent and mix gently.

#### TEST PROCEDURE

Allow the test card and cell reagent to reach room temperature before use.

- 1- Label the appropriate gel card with patient's name or identification number. Remove the aluminium foil carefully.
- 2- Pipette 10-12 µL (5 %) or 50 µL (0.8 %) of the patient's red cell suspension to the microtubes 1– 8.
- 3- Centrifuge the gel cards for 10 minutes in the gel card centrifuge.
- 4- Read the results.

#### Interpretation of the results:



|          |                                                                                                              |
|----------|--------------------------------------------------------------------------------------------------------------|
| Negative | Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column. |
| +/-      | Scarce small-sized agglutinations in the lower half of the column                                            |
| +1       | Some small-sized agglutinations in the column                                                                |
| +2       | Small or medium-sized agglutinations throughout the column                                                   |
| +3       | Most agglutinated red blood cells remain in the upper half of the gel column                                 |
| +4       | Agglutinated red blood cells form a band at the top of the gel column                                        |

The microtube Ctl must be negative. If it is positive, invalidate the test.

A positive reaction (+ to +++) indicates presence of the corresponding antigen.

#### Stability of the results:

it is recommended an immediate reading of the results after centrifuging the cards. Do not leave processed cards in horizontal position. If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2-8°C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.



Harmful

**LIMITATIONS**

- The Aluminum foil on the matrix gel card should be removed carefully and gently.
- Do not use matrix gel cards, which show signs of drying.
- The gel cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use
- Do not touch the pipet tip inside the gel card. Carryover problem may occur.
- Do not use haemolysed, cloudy or contaminated samples or with clot presence.
- Use of suspension solutions other than Fortress LISS Diluent may modify the reactions.
- Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test.
- Transfused patients or those subjected to bone marrow transplant may present images of double population.
- Patients with high-potency cold antibodies may coat the red blood cells completely so spontaneous agglutination may occur.

**PERFORMANCE CHARACTERISTICS**

**Diagnostic sensitivity and specificity of**

**Fortress Rh Phenotype +K Test:**

The diagnostic sensitivity and specificity of the antibodies present in Fortress Rh Phenotype +K Test for the determination of the antigens have been studied in a representative number of positive and negative samples. In the study using 1004 blood samples the sensitivity and specificity Fortress Rh Phenotype +K Test was as follows.

|        | Number of samples | Sensitivity | Specificity |
|--------|-------------------|-------------|-------------|
| Anti-C | 1004              | % 100       | % 100       |
| Anti-c | 1004              | % 100       | % 100       |
| Anti-E | 1004              | % 100       | % 100       |
| Anti-e | 1004              | % 100       | % 100       |
| Anti K | 1004              | % 100       | % 100       |

**Precision:**

Inter-assay and intra-assay reproducibility tests; no false positive or false negative results were obtained. Variability between agglutination view in positive samples were 1 agglutination rating or less in all assays

**BIBLIOGRAPHY**

- 1-Marion E. Reid Christine Lomas-Francis THE BLOOD GROUP ANTIGEN FactsBook 2004, Elsevier Ltd.
- 2-Geoff Daniels, Imelda Bromilow, Essential Guide to Blood Groups , Blackwell Publications 2010
- 3-Harvey G. Klein MD, David J. Anstee ,Mollison’s Blood Transfusion in Clinical Medicine 11TH EDITION Blackwell Publishing Ltd
- 4-Michael F. Murphy , Derwood H. Pamphilon , Practical Transfusion Medicine Third Edition 2009 Blackwell Publishing Ltd
- 5-Neil D. Avent and Marion E. Reid , The Rh blood group system: Blood 2000 95: 375-387
- 6- Christopher D. Hillyer ... BLOOD BANKING AND TRANSFUSION MEDICINE, Second Edition, Churchill Livingstone 2007
- 7-Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109 –113.

**Label Symbols**



Batch code



Expiry date



Storage temperature



Consult instructions for use



In vitro diagnostic medical device



Product code



Manufacturer