

ULTRASOFT

MIZO

BY DENTIGROUP

ADULT



THREE KEY FEATURES

The essentials of a gentler everyday clean.

01

DENSELY TUFTED BRISTLES

Gentle gum line care

Densely tufted bristles help clean comfortably along the gum line while remaining gentle on sensitive gums.

02

0.01 MM TAPERED BRISTLES

Ultra-soft cleaning

Fine tapered bristle tips provide a soft, gentle clean while helping reach around the teeth and along the gum line.

03

COMFORT GRIP HANDLE

Steady everyday control

Designed for a secure and comfortable hold, giving you better control throughout your daily brushing routine.


HOW TO BRUSH

Simple daily brushing guidance.




Brush for two minutes

Brush gently, twice a day.




Angle at the gum line

Position the bristles at approximately 45° towards the gum line.



Use gentle circular movements

Use small circular motions rather than forceful scrubbing.



Spit, then do not rinse immediately

Allow the toothpaste to remain on the teeth a little longer.

MIZO by DentiGroup | UltraSoft Adult Toothbrush

Replace your toothbrush every three to six months, or earlier if the bristles become worn or frayed.