





info@matrixcellbio.com



+82(02)970-9915



www.matrixcellbio.com
Seoul National Univ. of Sci. & Tech.,




The 2nd startup incubation center 131, Gongneung-ro 232, Nowon-gu, Seoul 01811, Republic of Korea


## APPLICATION BIOINK(GEL) ALIQUOTER

## Experimental data and results

ALIQUATED GEL




Aliquoting precision up to 96%




Linearity up to 98%

## CELL VIABILITY AND DRUG RELEASE



Viability up to 88%



## ALIQUOTED GEL IMAGES

Fig. A




Fig. B

Total weight (g) without syringe	Height
15kg	450 mm