

Ankura™

AAA Stent Graft System




Excellent Fixation, Advanced Conformability



Lifetech
SCIENTIFIC

| CE 0344



Unique designed **barbs** with excellent durability for better suparenal fixation



Tip capture mechanism for precise deployment



Spiral woven connecting bar for great flexibility



New **markers** for better visibility and indication



Over **10** years clinical experiences

Stable Support, Minimized Endoleak


- 6 barbs of the bare stent avoid migration effectively
- e-PTFE dual membrane provides low permeability and excellent biocompatibility
- Self-expanding Nitinol stent offers stable radial force
- Smaller waves on the main body improve the flexibility

Easy to Use, Simple to Assemble


- 3 piece components meet various clinical needs
- 7 radiopaque markers facilitate accurate positioning
- Hydrophilic coating sheath provides easier access




Ankura™ AAA Stent Graft Cases



Tortuous neck



Tortuous neck



Short neck




Huge aneurysm



Conical aortic neck

Product Specification



- D0 Proximal diameter
- D1 Ipsilateral/ Contralateral limb diameter
- L1 Main body covered length (Ipsilateral limb length included)
- L2 Main body covered length
- L3 Contralateral limb length

Main Body Specification

Code	D0 (mm)	D1 (mm)	L1 (mm)	L2 (mm)	L3 (mm)	Delivery System (Fr.)
AAA2010B120	20	10	120	50	30	21
AAA2210B120	22	10	120	50	30	21
AAA2412B120	24	12	120	50	30	21
AAA2612B120	26	12	120	50	30	21
AAA2812B120	28	12	120	50	30	23
AAA3012B120	30	12	120	50	30	23
AAA3212B120	32	12	120	50	30	23
AAA3412B120	34	12	120	50	30	23



Iliac Extension Product Specification

Code	Proximal Diameter (D1) (mm)	Distal Diameter (D2) (mm)	Length (L1) (mm)	Delivery System (Fr.)
cuff1212C040	12	12	40	18
cuff1212C060	12	12	60	18
cuff1212C080	12	12	80	18
cuff1212C100	12	12	100	18
cuff1212C120	12	12	120	18
cuff1410C040	14	10	40	18
cuff1410C060	14	10	60	18
cuff1410C080	14	10	80	18
cuff1410C100	14	10	100	18
cuff1410C120	14	10	120	18
cuff1412C040	14	12	40	18
cuff1412C060	14	12	60	18
cuff1412C080	14	12	80	18
cuff1412C100	14	12	100	18
cuff1412C120	14	12	120	18
cuff1414C040	14	14	40	18
cuff1414C060	14	14	60	18
cuff1414C080	14	14	80	18
cuff1414C100	14	14	100	18
cuff1414C120	14	14	120	18
cuff1416C040	14	16	40	18
cuff1416C060	14	16	60	18
cuff1416C080	14	16	80	18
cuff1416C100	14	16	100	18
cuff1416C120	14	16	120	18
cuff1418C040	14	18	40	18
cuff1418C060	14	18	60	18
cuff1418C080	14	18	80	18
cuff1418C100	14	18	100	18
cuff1418C120	14	18	120	18
cuff1420C040	14	20	40	18
cuff1420C060	14	20	60	18
cuff1420C080	14	20	80	18
cuff1420C100	14	20	100	18
cuff1420C120	14	20	120	18



Iliac Extension Product Specification

Code	Proximal Diameter (D1) (mm)	Distal Diameter (D2) (mm)	Length (L1) (mm)	Delivery System (Fr.)
cuff1422C040	14	22	40	18
cuff1422C060	14	22	60	18
cuff1422C080	14	22	80	18
cuff1422C100	14	22	100	18
cuff1422C120	14	22	120	18
cuff1616C040	16	16	40	18
cuff1616C060	16	16	60	18
cuff1616C080	16	16	80	18
cuff1616C100	16	16	100	18
cuff1616C120	16	16	120	18
cuff1618C040	16	18	40	18
cuff1618C060	16	18	60	18
cuff1618C080	16	18	80	18
cuff1618C100	16	18	100	18
cuff1618C120	16	18	120	18
cuff1620C040	16	20	40	18
cuff1620C060	16	20	60	18
cuff1620C080	16	20	80	18
cuff1620C100	16	20	100	18
cuff1620C120	16	20	120	18
cuff1622C040	16	22	40	18
cuff1622C060	16	22	60	18
cuff1622C080	16	22	80	18
cuff1622C100	16	22	100	18
cuff1622C120	16	22	120	18

Lifetech Scientific (Shenzhen) Co., Ltd.

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park,
Nanshan District, Shenzhen 518057, P. R. China
Tel: +86-755-86026250 Fax: +86-755-86026251
Email: marketing@lifetechmed.com
Web: www.lifetechmed.com

EU authorized representative:
LifeTech Scientific (Europe) Coöperatief U.A.
Oudemolenstraat 60, 6416CB Heerlen, the Netherlands