


Order Information

CE
0344

CeraFlex™ ASD


Occluder Product Specification



Code	[D]Waist Diameter (mm)	[D1]Right Disc Diameter (mm)	[D2]Left Disc Diameter (mm)	[L]Waist Length (mm)	Matching Sheath Size SteerEase™ (Fr.)
LT-ASDf-06	6	14	18	4	SFA8F-f
LT-ASDf-08	8	16	20	4	SFA8F-f
LT-ASDf-10	10	18	22	4	SFA8F-f
LT-ASDf-12	12	22	26	4	SFA9F-f
LT-ASDf-14	14	24	28	4	SFA10F-f
LT-ASDf-16	16	26	30	4	SFA10F-f
LT-ASDf-18	18	28	32	4	SFA10F-f
LT-ASDf-20	20	30	34	4	SFA12F-f
LT-ASDf-22	22	32	36	4	SFA12F-f
LT-ASDf-24	24	34	38	4	SFA12F-f
LT-ASDf-26	26	36	40	4	SFA12F-f
LT-ASDf-28	28	38	42	4	SFA12F-f
LT-ASDf-30	30	40	44	4	SFA14F-f
LT-ASDf-32	32	42	46	4	SFA14F-f

CeraFlex™ MF-ASD


Occluder Product Specification



Code	[D1]Right Disc Diameter (mm)	[D2]Left Disc Diameter (mm)	[L]Waist Length (mm)	Matching Sheath Size SteerEase™ (Fr.)
LT-ASDf-MF-1818	18	18	3	SFA9F-f
LT-ASDf-MF-2518	25	18	3	SFA10F-f
LT-ASDf-MF-2525	25	25	3	SFA10F-f
LT-ASDf-MF-3030	30	30	3	SFA12F-f
LT-ASDf-MF-3525	35	25	3	SFA14F-f
LT-ASDf-MF-3535	35	35	3	SFA14F-f
LT-ASDf-MF-4040	40	40	3	SFA14F-f

CeraFlex™ PDA


Occluder Product Specification



Code	[D]Retention Skirt (mm)	[D1]Proximal Diameter (mm)	[D2]Distal Diameter (mm)	[L]Waist Length (mm)	Matching Sheath Size SteerEase™ (Fr.)
LT-PDAf-0406	10	4	6	7	SFP6F-f
LT-PDAf-0608	12	6	8	7	SFP7F-f
LT-PDAf-0810	14	8	10	7	SFP7F-f
LT-PDAf-1012	16	10	12	7	SFP8F-f
LT-PDAf-1214	20	12	14	7	SFP9F-f
LT-PDAf-1416	22	14	16	8	SFP9F-f
LT-PDAf-1618	24	16	18	8	SFP10F-f
LT-PDAf-1820	26	18	20	9	SFP12F-f
LT-PDAf-2022	28	20	22	9	SFP12F-f
LT-PDAf-2224	30	22	24	10	SFP14F-f

CeraFlex™ PFO



Occluder Product Specification



Code	[D1]Right Disc Diameter (mm)	[D2]Left Disc Diameter (mm)	[L]Waist Length (mm)	Matching Sheath Size SteerEase™ (Fr.)
LT-PFOf-1818	18	18	3	SFA9F-f
LT-PFOf-2518	25	18	3	SFA10F-f
LT-PFOf-3030	30	30	3	SFA12F-f
LT-PFOf-3525	35	25	3	SFA14F-f

CeraFlex™

Innovative pre-loaded occluder device



Lifetech Scientific

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: 86-755-86026250

Fax: 86-755-86026251

Web: www.lifetechmed.com

Lifetech Scientific India Private Limited

Add: #3, III Floor, 1st Main, BDA Layout Kodihalli, HAL 2nd Stage, Bangalore, India 560008

Tel: +91-80-41245252 / 41245129

Fax: +91-80-25200364

Lifetech Scientific (Europe) Coöperatief U.A.

Add: Kruisdonk 64, 6222 PH Maastricht, Netherlands




Version NO.: 2014-EN-A-1.0



CE
0344

Innovative pre-loaded delivery device

Easy and convenient to use



Maximum flexibility—Pivot 360°

Accurate positioning during procedure.

Proprietary Titanium nitride coating technology

Accelerates endothelialization;
Prevents 93% of nickel release compared with traditional uncoated occluders.

Minimizing the amount of material on the left atrial disc

Avoid thrombotic complications.



IF design award nomination

CeraFlex™ streamline handle design to fits perfectly
in our hands.

