


NXT ▶ MED

Official German distribution Partner of  **APT**Medical

INTERVENTIONAL CARDIOLOGY





March[™]
Guiding Catheter

PREMIUM PERFORMANCE YOU CAN EXPECT

March™ guiding catheter combines superior torqueability and back-up support due to braided structure and polymer material. Together with soft tip and large lumen design to give you more space for your procedures.

Optimal Back-up Support

- Polymer material with braided design provides excellent kink resistance, back-up support and torque response




Soft and Atraumatic Tip

- Decreased layers of the tip to make sure the softness and safety

Large Lumen Design

6F: 0.071"

- Flat wire braided structure for a bigger lumen
- Provides smooth and compatible pathway for multiple and various devices



MULTI-SEGMENT DESIGN

Torque/Push Segment

- Stiff segment ensures 1:1 torque response promotes quicker catheter positioning with less whipping

Anti-kink Segment

- A gradually decreasing stiffness between torque/push and co-axial segments to preventing kink when torque
- Retain the Shape

Co-Axial Segment

- Co-Axial segments provide easier entrance of devices into the coronary artery
- Retain the shape

Radiopaque Tip

- Radiopaque soft tip ensures atraumatic navigation and clear visibility

	5 Fr	6 Fr	7 Fr
I.D.(Inch/mm)	0.058" / 1.47mm	0.071" / 1.80mm	0.081" / 2.06mm
O.D.(Inch/mm)	0.069" / 1.75mm	0.0815" / 2.07mm	0.093" / 2.36mm

ORDERING INFORMATION

Curve	Tip	5F (100cm)	6F (100cm)	6.5F (100cm)	7F (100cm)	8F (100cm)
Extra Back-up Left	ABU3.0	53510096	53610094	53210094	53710094	53810094
Extra Back-up Left	ABU3.25	53510098	53610096	53210096	-	-
Extra Back-up Left	ABU3.5	535100100	53610098	53210098	53710096	53810096
Extra Back-up Left	ABU3.75	535100102	536100100	532100100	53710098	53810098
Extra Back-up Left	ABU4.0	535100104	536100102	532100102	537100100	538100100
Extra Back-up Left	XB3.0	53510004	53610001	53210004	53710004	53810004
Extra Back-up Left	XB3.5	53510008	53610002	53210008	53710008	53810008
Extra Back-up Left	XB4.0	53510012	53610003	53210012	53710012	53810012
Judkins Left	JL3.0	53510016	53610022	53210016	53710016	53810016
Judkins Left	JL3.5	53510020	53610004	53210020	53710020	53810020
Judkins Left	JL4.0	53510024	53610005	53210024	53710024	53810024
Judkins Left	JL4.5	53510028	53610006	53210028	53710028	53810028
Judkins Left	JL5.0	53510032	53610007	53210032	53710032	53810032
Judkins Right	JR3.0	53510036	53610023	53210036	53710036	53810036
Judkins Right	JR3.5	53510040	53610008	53210040	53710040	53810040
Judkins Right	JR4.0	53510044	53610009	53210044	53710044	53810044
Judkins Right	JR4.5	53510048	53610010	53210048	53710048	53810048
Judkins Right	JR5.0	53510052	53610011	53210052	53710052	53810052
Amplatz Left	AL1	53510056	53610012	53210056	53710056	53810056
Amplatz Left	AL2	53510060	53610013	53210060	53710060	53810060
Amplatz Left	AL3	53510064	-	53210064	53710064	53810064
Amplatz Right	AR	53510068	53610015	53210068	53710068	53810068
Amplatz Right	AR1	53510072	53610016	53210072	53710072	53810072
Amplatz Right	AR2	53510076	53610017	53210076	53710076	53810076
Multipurpose	MPA	53510084	53610019	53210084	53710084	53810084
Multipurpose	MPB	53510088	53610021	53210088	53710088	53810088

Curve (with side holes)	Tip	5F (100cm)	6F (100cm)	6.5F (100cm)	7F (100cm)	8F (100cm)
Extra Back-up Left	ABU3.0	53510095	53610093	53210093	53710093	53810099
Extra Back-up Left	ABU3.25	53510097	53610095	53210095	-	-
Extra Back-up Left	ABU3.5	53510099	53610097	53210097	53710095	53810093
Extra Back-up Left	ABU3.75	535100101	53610099	53210099	53710097	53810095
Extra Back-up Left	ABU4.0	535100103	536100101	532100101	53710099	53810097
Extra Back-up Left	XB3.0	53510002	53610026	53210002	53710002	53810002
Extra Back-up Left	XB3.5	53510006	53610028	53210006	53710006	53810006
Extra Back-up Left	XB4.0	53510010	53610031	53210010	53710010	53810010
Judkins Left	JL3.0	53510014	53610034	53210014	53710014	53810014
Judkins Left	JL3.5	53510018	53610037	53210018	53710018	53810018
Judkins Left	JL4.0	53510022	53610040	53210022	53710022	53810022
Judkins Left	JL4.5	53510026	53610043	53210026	53710026	53810026
Judkins Left	JL5.0	53510030	53610046	53210030	53710030	53810030
Judkins Right	JR3.0	53510034	53610049	53210034	53710034	53810034
Judkins Right	JR3.5	53510038	53610052	53210038	53710038	53810038
Judkins Right	JR4.0	53510042	53610055	53210042	53710042	53810042
Judkins Right	JR4.5	53510046	53610058	53210046	53710046	53810046
Judkins Right	JR5.0	53510050	53610061	53210050	53710050	53810050
Amplatz Left	AL1	53510054	53610064	53210054	53710054	53810054
Amplatz Left	AL2	53510058	53610067	53210058	53710058	53810058
Amplatz Left	AL3	53510062	53610070	53210062	53710062	53810062
Amplatz Right	AR	53510066	53610073	53210066	53710066	53810066
Amplatz Right	AR1	53510070	53610076	53210070	53710070	53810070
Amplatz Right	AR2	53510074	53610079	53210074	53710074	53810074
Multipurpose	MPA	53510082	53610085	53210082	53710082	53810082
Multipurpose	MPB	53510086	53610088	53210086	53710086	53810086

NXT MED

Official Distribution Partner

Address: Fritz-Vomfelde-Straße 6
40547 Düsseldorf, Germany
Email: hello@nxt-med.com
Tel: +49 211 58678600

@NXT MED



www.nxt-med.com

APT Medical

Manufacturer

Address: Room 601, Building B, Tongfang Information Port,
No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China
Email: international@aptmed.com **Tel:** +86 755-83591013
PN: ENG-6.5F Set-210X285X4P-20230922

@APT Medical



www.aptmed.com

CE 1639