

VASCULAR GRAFTS

VASCULAR
polythese®
polyarch
polybranch

Woven Polyester




Value Life

polythese® polyarch - polybranch WOVEN POLYESTER...

... intended for **replacement or bypass of arteries** during open-surgical vascular repair procedures.

The **polythese polyarch & polybranch** prosthesis are especially indicated for replacement or bypass of **the ascending aorta, aortic arch and thoracic aorta**.

Adapted for **high pressure area!**



1 OR 4 BRANCHES !



4 branches corresponding to the brachiocephalic, common left carotid and subclavian arteries.
They can also have a 4th branch intended for introduction of extracorporeal circulation.

polythese polyarch

Impregnated with type I bovine collagen on both side

1 branch with a small diameter intended for introduction of extracorporeal circulation.



polythese polybranch

MAKE THE DIFFERENCE

1, 3 or 4 branches depending your need

No formaldehyde nor glutaraldehyde included during manufacturing

REFERENCES

poly**these**[®] polyarch
(3 branches)



Reference	Body diameter (mm)	Body length (cm)	Branches diameter (brachiocephalic/common carotid/subclavian) (mm)	Branches len- gth (cm)
-----------	-----------------------	---------------------	--	------------------------------

poly**these**[®] polyarch
(4 branches)



Reference	Body diameter (mm)	Body length (cm)	Branches diameter (brachiocephalic/common carotid/subclavian/ECC) (mm)	Branches len- gth (cm)
-----------	-----------------------	---------------------	--	------------------------------



Reference	Body diameter (mm)	Body length (cm)	Branch diameter (mm)	Branch length (cm)
PTB1608	16	40	8	30
PTB1808	18			
PTB2008	20			
PTB2208	22			
PTB2408	24			
PTB2608	26			
PTB2808	28			
PTB3008	30			
PTB3208	32			
PTB3408	34			
PTB3608	36			
PTB3808	38			
PTB1610	16		10	
PTB1810	18			
PTB2410	24			
PTB2610	26			
PTB2810	28			
PTB3010	30			
PTB3210	32			
PTB3410	34			
PTB3610	36			
PTB3810	38			

FOR FURTHER INFORMATION, PLEASE CONTACT: questions@vygon.com

The specification shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

VYGON – 5, RUE ADELIN – 95440 ECOUEN - FRANCE

RECEPTION: +33 (0)1 39 92 63 63 – SERVICE CLIENTS FRANCE: +33 (0)1 39 92 63 81

EXPORT CUSTOMER SERVICE: +33 (0)1 39 92 64 15

www.vygon.com