


VASCULAR GRAFTS

VASCULAR
polythese®
IC/ICT

Woven Polyester



Value Life

polythese® IC/ICT
WOVEN POLYESTER...


... intended for **replacement or bypass of arteries** during open-surgical vascular repair procedures.

The **polythese IC-ICT** prosthesis is especially indicated in **thoracic, abdominal and peripheral repair surgery**.

Adapted for **high pressure area!**



DESIGNED FOR A **BETTER POSITIONING**...



MAKE THE DIFFERENCE

Two black guide lines, single and double lines opposite each other over the entire length of the prosthesis

► To avoid the “twist” of the prosthesis during the surgery.

No formaldehyde nor glutaraldehyde included during manufacturing

REFERENCES

polythese® IC straight tube

		Minimum usable lengths (cm)				
		15 cm	40 cm	60 cm	70 cm	100 cm
Nominal relaxed and pressurized internal diameters (mm)	6 mm	1201 0602	-	1201 0600	1201 0603	1201 0601
	7 mm	1201 0702	-	1201 0700	1201 0703	1201 0701
	8 mm	1201 0802	-	1201 0800	1201 0803	1201 0801
	10 mm	1201 1002	-	1201 1000	1201 1003	1201 1001
	12 mm	1201 1202	-	1201 1200	1201 1203	-
	14 mm	1201 1402	1201 1400	1201 1401	-	-
	16 mm	1201 1602	1201 1600	1201 1601	-	-
	18 mm	1201 1802	1201 1800	1201 1801	-	-
	20 mm	1201 2002	1201 2000	1201 2001	-	-
	22 mm	1201 2202	1201 2200	1201 2201	-	-
	24 mm	1201 2402	1201 2400	1201 2401	-	-
	26 mm	1201 2602	1201 2600	1201 2601	-	-
	28 mm	1202 2802	1201 2800	1201 2801	-	-

polythese® IC bifurcated



Reference	Main tube nominal relaxed and pressurized inner diameter (mm)	Branches nominal relaxed and pressurized inner diameter (mm)	Minimum usable overall length (cm)
1202 1206	12	6	45
1202 1407	14	7	
1202 1608	16	8	
1202 1809	18	9	
1202 2010	20	10	
1202 2211	22	11	
1202 2412	24	12	

polythese® ICT straight thoracic tube



		Minimum usable lengths (cm)		
		15 cm	30 cm	60 cm
Nominal relaxed and pressurized internal diameters (mm)	30 mm	1301 3002	1301 3000	1301 3001
	32 mm	1301 3202	1301 3200	1301 3201
	34 mm	1301 3402	1301 3400	1301 3401
	36 mm	1301 3602	1301 3600	1301 3601
	38 mm	1301 3802	1301 3800	1301 3801

04/2023 - DB BU CS 18 450 E

FOR FURTHER INFORMATION, PLEASE CONTACT: questions@vygon.com

The specification shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

VYCON – 5, RUE ADELIN – 95440 ECOUEN - FRANCE

RECEPTION: +33 (0)1 39 92 63 63 – SERVICE CLIENTS FRANCE: +33 (0)1 39 92 63 81

EXPORT CUSTOMER SERVICE: +33 (0)1 39 92 64 15

www.vygon.com

