

Sistema de Cardioplégia desechable Unique®

Apirógeno. Compuesto por un intercambiador de calor, manómetro y circuito de líneas con filtros para soluciones cristaloides.

Permite la infusión de solución cardioplégica, además de promover la mezcla de sangre oxigenada con una solución preformulada (4:1).



| Código interno | Descripción | Presentación |
|----------------|-----------------------------------|--------------|
| OX-CTD41 | Brazo soporte para sistema Unique | Pieza |
| OX-EDR17 | Sistema de Cardioplégia Unique | Pieza |



UNIQUE

Cardioplegia System

Special Features

- Flexibility and safety for a wide range of use
- Effective performance in a compact low priming design
- Pressure manometer and temperature probe port
- High heat exchanger performance
- Sterilization : ETO

Specifications

| | |
|--|----------------------|
| Heat exchanger effective area | 0.02 m ² |
| Heat Exchanger priming volume | 30 mL |
| Equipment priming volume | 105 mL |
| System total volume | 135 mL |
| Crystalloid solution filter pore diameter | 0.8 microns |
| Filtering area crystalloid solution filter | 14.5 cm ² |

Materials

| | |
|----------------|--------------------|
| Heat exchanger | Anodized Aluminium |
| Central Valve | Polyethylene |