

HTContainer

s t o r e

Honesty and T rust



Sterilization Container Systems

Surgical Trays & Wire Baskets



HTContainer
store



TABLE OF CONTENTS

04-10 **Introduction**

10-19 **Sterilization Containers**

20-28 **Perforated Trays & Wire
Baskets**

30-34 **Components**



Type and Size Variations

The Sterilization Containers are manufacturing arange of sizes, offering specific solutions to meet the requirements of healthcare professionals.

Designed with user-friendly features that improve convenience in packaging and storage processes, optimizing workflow efficiency for CSSD.



Standard Height Sizes

(mm)

100

130

150


200

260

There are 5 different height options available in standard container types

Specific Needs Particular Solutions

The Sterilization Containers are typically available in three standard sizes: full (1/1), medium (3/4), and half (1/2) which meet flexibility and compatibility with a wide range of surgical instrument sets and medical equipment. Additionally, to meet specific needs, we offer various types of containers, providing particular solutions.





Colors



The containers offer six different color choices.



Ease of Identification

The containers feature color-coded lids, offering a quick and suitable method for categorization by color, which is easy to identify.





STANDARD TYPE CONTAINER SYSTEMS

Model 1

10

Full Size (1/1)

Model 6

16

Flat Size

Model 2

11

Quarter Size (3/4)

Model 7

17

Mini Size

Model 3

12

Half Size (1/2)

Model 8

18

Dental Size

Model 4

13

Long Body

Model 9

19

1/2 Dental Size

Model 5

14

Large Body



Lid Perforated / Base Closed



Lid Safety Lid / Base Perforated



Lid & Base Perforated



Lid Safety Lid / Base Closed



STANDARD TYPE

MODEL 1

Full Size (1/1)



Base
Closed

Outside Dimensions (mm)

Outside Dimensions (mm)	Code
585X280X100	C1-B-10
585x280x135	C1-B-13
585x280x150	C1-B-15
585x280x200	C1-B-20
585X280X260	C1-B-26

585x280 (mm)

Code



Perforated Lid

●	C1-LF-01
●	C1-LF-02
●	C1-LF-03
●	C1-LF-04
●	C1-LF-05
●	C1-LF-06

Outside Dimensions (mm)

Outside Dimensions (mm)	Code
585X280X100	C1-BP-10
585x280x135	C1-BP-13
585x280x150	C1-BP-15
585x280x200	C1-BP-20
585X280X260	C1-BP-26



Base
Perforated *

585x280 (mm)

Code




Safety Lid

●	C1-LS-01
●	C1-LS-02
●	C1-LS-03
●	C1-LS-04
●	C1-LS-05
●	C1-LS-06

WIRE BASKETS

Outside Dimensions (mm)

Code




Full Size (1/1)

540x255x30	B1-BW-03
540x255x50	B1-BW-05
540x255x70	B1-BW-07
540x255x100	B1-BW-10

PERFORATED BASKETS

Outside Dimensions (mm)

Code



Full Size (1/1)

540x255x30	B1-BP-03
540x255x50	B1-BP-05
540x255x70	B1-BP-07
540x255x100	B1-BP-10



STANDARD TYPE

MODEL 2

Quarter Size (3/4)



Base
Closed


Outside Dimensions (mm)

Code

465X280X100	C2-B-10
465x280x135	C2-B-13
465x280x150	C2-B-15
465x280x200	C2-B-20

465x280 (mm)

Code



Lid
Perforated

●	C2-LF-01
●	C2-LF-02
●	C2-LF-03
●	C2-LF-04
●	C2-LF-05
●	C2-LF-06



Base
Perforated *

Outside Dimensions (mm)

Code

465X280X100	C2-BP-10
465x280x135	C2-BP-13
465x280x150	C2-BP-15
465x280x200	C2-BP-20

465x280 (mm)

Code



Safety
Lid

●	C2-LS-01
●	C2-LS-02
●	C2-LS-03
●	C2-LS-04
●	C2-LS-05
●	C2-LS-06

WIRE BASKETS




Quarter Size (3/4)

Outside Dimensions (mm)

Code

405x255x30	B2-BW-03
405x255x50	B2-BW-05
405x255x70	B2-BW-07
405x255x100	B2-BW-10

PERFORATED BASKETS



Quarter Size (3/4)

Outside Dimensions (mm)

Code

405x255x30	B2-BP-03
405x255x50	B2-BP-05
405x255x70	B2-BP-07
405x255x100	B2-BP-10

* Product code does not include the filters,
for visual purposes only.



STANDARD TYPE

MODEL 3

Half Size (1/2)



Base
Closed

Outside Dimensions (mm)

Code

295X280X100	C3-B-10
295x280x135	C3-B-13
295x280x150	C3-B-15
295x280x200	C3-B-20
295X280X260	C3-B-26



Lid
Perforated

295x280 (mm)

Code

	C3-LF-01
	C3-LF-02
	C3-LF-03
	C3-LF-04
	C3-LF-05
	C3-LF-06



Base
Perforated

Outside Dimensions (mm)

Code

295X280X100	C3-BP-10
295x280x135	C3-BP-13
295x280x150	C3-BP-15
295x280x200	C3-BP-20
295X280X260	C3-BP-26



Safety
Lid

295x280 (mm)

Code


	C3-LS-01
	C3-LS-02
	C3-LS-03
	C3-LS-04
	C3-LS-05
	C3-LS-06

WIRE BASKETS

Outside Dimensions (mm)

Code

245x255x30	B3-BW-03
245x255x50	B3-BW-05
245x255x70	B3-BW-07
245x255x100	B3-BW-10




Half Size (1/2)

PERFORATED BASKETS

Outside Dimensions (mm)

Code

245x255x30	B3-BP-03
245x255x50	B3-BP-05
245x255x70	B3-BP-07
245x255x100	B3-BP-10



Half Size (1/2)



STANDARD TYPE

MODEL 4

Long Body

Outside Dimensions (mm)

Code

710x280x200

C4-B-20



Base
Closed

710x280 (mm)

Code

- C4-LF-01
- C4-LF-02
- C4-LF-03
- C4-LF-04
- C4-LF-05
- C4-LF-06



Lid
Perforated

Outside Dimensions (mm)

Code

710x280x200

C4-BP-20



Base
Perforated *

WIRE BASKETS

Outside Dimensions (mm)

Code

655x255x30

B4-BW-03

655x255x50

B4-BW-05

655x255x70

B4-BW-07

655x255x100

B4-BW-10



Long Body

PERFORATED BASKETS

Outside Dimensions (mm)

Code

655x255x30

B4-BP-03

655x255x50


B4-BP-05

655x255x70

B4-BP-07

655x255x100

B4-BP-10



Long Body

* Product code does not include the filters,
for visual purposes only.



STANDARD TYPE

MODEL 5

Wide Body



Base
Closed

Outside Dimensions (mm)

Code

600x400x100

C5-B-10

600x400x200

C5-B-20



Lid
Perforated

600x400 (mm)

Code

- C5-LF-01
- C5-LF-02
- C5-LF-03
- C5-LF-04
- C5-LF-05
- C5-LF-06



Base
Perforated *

Outside Dimensions (mm)

Code

600x400x100

C5-BP-10

600x400x200

C5-BP-20

WIRE BASKETS

Outside Dimensions (mm)

Code



Wide Body

530x330x30

B5-BW-03

530x330x50

B5-BW-05

530x330x70

B5-BW-07

530x330x100

B5-BW-10

PERFORATED BASKETS

Outside Dimensions (mm)

Code



Wide Body

530x330x30

B5-BP-03

530x330x50

B5-BP-05

530x330x70


B5-BP-07

530x330x100

B5-BP-10

HTContainer

s t o r e



STANDARD TYPE

MODEL 6

Flat Size



Base
Closed



Base
Perforated *

Outside Dimensions (mm)	Code
295x280x50	C6-B-05
295x280x70	C6-B-07

Outside Dimensions (mm)	Code
295x280x70	C6-BP-07



Lid
Perforated

295x280 (mm)	Code
	C6-LF-01
	C6-LF-02
	C6-LF-03
	C6-LF-04
	C6-LF-05
	C6-LF-06

WIRE BASKETS


Outside Dimensions (mm)	Code
245x255x30	B6-BW-03
245x255x50	B6-BW-05



Half Size (1/2)

PERFORATED BASKETS

Outside Dimensions (mm)	Code
245x255x30	B6-BP-03
245x255x50	B6-BP-05




Half Size (1/2)

STANDARD TYPE

MODEL 7

Mini Size



Base
Closed

Outside Dimensions (mm)

Code

310x130x40

C7-B-04

310x130x70

C7-B-07

310x130x100

C7-B-10



Base
Perforated *

Outside Dimensions (mm)

Code

310x130x40


C7-BP-04

310x130x70

C7-BP-07

310x130x100

C7-BP-10



Lid
Perforated

310x140 (mm)

Code

● C7-LF-01

● C7-LF-02

● C7-LF-03

● C7-LF-04

● C7-LF-05

● C7-LF-06

WIRE BASKETS

Outside Dimensions (mm)

Code

270x125x50

B7-BW-05

270x125x70

B7-BW-07



Mini Size

PERFORATED BASKETS

Outside Dimensions (mm)

Code

270x125x50

B7-BP-05

270x125x70

B7-BP-07



Mini Size

* Product code does not include the filters,
for visual purposes only.

STANDARD TYPE

MODEL 8

Dental Size



Base
Closed



Base
Perforated *

Outside Dimensions (mm)

Code

310x190x40	C8-B-04
310x190x60	C8-B-06
310x190x80	C8-B-08
310x190x100	C8-B-10
310x190x130	C8-B-13

310x190x40	C8-BP-04
310x190x60	C8-BP-06
310x190x80	C8-BP-08
310x190x100	C8-BP-10
310x190x130	C8-BP-13

310x190 (mm)

Code

●	C8-LF-01
●	C8-LF-02
●	C8-LF-03
●	C8-LF-04
●	C8-LF-05
●	C8-LF-06



Lid
Perforated

WIRE BASKETS

Outside Dimensions (mm)

Code

270x170x50	B8-BW-05
270x170x70	B8-BW-07
270x170x100	B8-BW-10



Dental Size

PERFORATED BASKETS

Outside Dimensions (mm)

Code

270x170x50	B8-BP-05
270x170x70	B8-BP-07
270x170x100	B8-BP-10



Dental Size

STANDARD TYPE

MODEL 9

1/2 Dental Size



Base
Closed

Outside Dimensions (mm)


210x155x50

210x155x60

Code

C9-B-05

C9-B-06



Lid
Perforated

210x155 (mm)

Code

- C9-LF-01
- C9-LF-02
- C9-LF-03
- C9-LF-04
- C9-LF-05
- C9-LF-06

WIRE BASKETS



1/1 Dental Size

Outside Dimensions (mm)

190x140x30

Code

B9-BW-03

PERFORATED BASKETS



1/2 Dental Size


Outside Dimensions (mm)

190x140x30

Code

B9-BP-03

* Product code does not include the filters,
for visual purposes only.



PERFORATED BASKETS

Standard Model


34

Optic, Mini & Dental Model


35

Micro, Small, Large Model

36-37



Mesh Dimensions 4x4 mm
Range of Mesh 1 mm



Suitable for effortless adaptation of the components



Easy handling and carrying even when the lid is closed

PERFORATED BASKETS

Standard Model



Types

Outside Dimensions (mm)

 Polymer Foot

Model 1

Full Size

B1-BP-03
540x255x30

Model 2

Middle Size

B2-BP-03
405x255x30

Model 3

Half Size

B3-BP-03
245x255x30

Model 4

Long Size

B4-BP-03
655x255x30

Model 5

Wide Size

B5-BP-03
530x330x30

Model 6

Flat Size

B6-BP-03
245x250x30



Base (With Feet)

B1-BP-05
540x255x50

B2-BP-05
405x255x50

B3-BP-05
245x255x50

B4-BP-05
655x255x50

B5-BP-05
530x330x50

B6-BP-05
245x250x50

B1-BP-07
540x255x70

B2-BP-07
405x255x70

B3-BP-07
245x255x70

B4-BP-07
655x255x70

B5-BP-07
655x252x70

B1-BP-10
540x255x100

B2-BP-10
405x255x100

B3-BP-10
245x255x100

B4-BP-10
655x255x100

B5-BP-10
530x330x100



Lid

B1-BP-00
540x255

B2-BP-00
405x255

B3-BP-00
245x255

B4-BP-00
655x255

B5-BP-00
530x330

B6-BP-00
245x250

Model 7

Mini Size

Model 8

Dental Size

Model 9

1/2 Dental Size

Model 48

DIN Size

B9-BP-03
190x140x30

B48-BP-03
485x255x30



Base (With Feet)

B7-BP-05
270x125x50

B8-BP-05
270x170x50

B48-BP-05
485x255x50

B7-BP-07
270x125x70

B8-BP-07
270x170x70

B48-BP-07
485x255x70

B8-BP-10
270x170x100

B48-BP-10
485x255x100



Lid

B7-BP-00
270x125

B8-BP-00
270x170

B48-BP-00
485x255

PERFORATED BASKETS

Optic Model



Types

Outside Dimensions (mm)

Single Scope



Base & Lid

B27-BLP-04S
272x83x44

B45-BLP-04S
452x83x44

Double Scope




Base & Lid

B27-BLP-04D
272x83x44

B45-BLP-04D
452x83x44

STANDARD

Optic Wire Baskets



Types

Outside Dimensions (mm)

Single Scope



Base & Lid

B29-BLW-05S1
290x80x55

B46-BLW-05S1
460x80x55

B67-BLW-05S1
670x80x55

Double Scope

B29-BLW-05D1
290x80x55

B46-BLW-05D1
460x80x55

B67-BLW-05D1
670x80x55



Silicone Holder
(Type 1)

Single Scope



Base & Lid


B29-BLW-05S2
290x80x55

B46-BLW-05S2
460x80x55


Double Scope

B29-BLW-05D2
290x80x55

B46-BLW-05D2
460x80x55



Silicone Holder
(Type 2)



PERFORATED TRAYS WITH HINGES


Micro& SmallTrays

Code	Outside Dimensions (mm)	Lock System	• Embossed lock system • 2x2mm perforated bottom and lid
B50-BLP-10	40x40x25	Embossed lock system	
B50-BLP-11	60x60x25	Embossed lock system	
B50-BLP-12	60x60x40	Embossed lock system	
B50-BLP-13	80x80x25	Embossed lock system	
B50-BLP-14	80x80x40	Embossed lock system	
B50-BLP-15	80x40x40	Embossed lock system	

Code	Outside Dimensions (mm)	Lock System	• Sliding flat lock system • 2x2mm perforated bottom and lid
B51-BLP-10	100x30x25	Push lock system	
B51-BLP-11	160x30x25	Push lock system	
B51-BLP-12	200x30x25	Push lock system	
B51-BLP-13	240x40x40	Push lock system	
B51-BLP-14	300x40x40	Push lock system	
B51-BLP-15	450x40x40	Push lock system	
B51-BLP-16	115x90x45	Sliding flat lock system	
B51-BLP-17	115x115x25	Sliding flat lock system	
B51-BLP-18	145x150x50	Sliding flat lock system	
B51-BLP-19	170x135x35	Sliding flat lock system	
B51-BLP-20	210x145x25	Sliding flat lock system	
B51-BLP-21	210x145x45	Sliding flat lock system	
B51-BLP-22	240x240x45	Sliding flat lock system	
B51-BLP-23	275x175x35	Sliding flat lock system	
B51-BLP-24	280x170x60	Sliding flat lock system	

PERFORATED TRAYS WITH HINGES

Micro& SmallTrays



Small Trays


Code	Outside Dimensions (mm)	Lock System	• Sliding flat lock system • 4x4mm perforated bottom and lid
B52-BLP-10	80x80x40	Push lock system	
B52-BLP-11	100x80x30	Sliding flat lock system	
B52-BLP-12	130x115x40	Sliding flat lock system	
B52-BLP-13	130x115x60	Sliding flat lock system	
B52-BLP-14	195x70x40	Sliding flat lock system	
B52-BLP-15	225x130x40	Sliding flat lock system	
B52-BLP-16	250x115x25	Sliding flat lock system	
B52-BLP-17	250x115x40	Sliding flat lock system	
B52-BLP-18	250x115x60	Sliding flat lock system	
B52-BLP-19	480x115x50	Sliding flat lock system	
B52-BLP-20	270x123x45	Sliding flat lock system	
B52-BLP-21	270x123x65	Sliding flat lock system	
B52-BLP-22	270x170x50	Sliding flat lock system	
B52-BLP-23	270x170x100	Sliding flat lock system	



*Dimensions: width x length x height (with tolerances up to ± 5 mm)


HTContainer

store



STANDARD

Wire Baskets



Types

Outside Dimensions (mm)

Model 1

Full Size

B1-BW-03
540x255x30

Model 2

Middle Size

B2-BW-03
405x255x30

Model 3

Half Size

B3-BW-03
245x255x30

Model 4

Long Size

B4-BW-03
655x255x30

Model 5

Wide Size

B5-BW-03
530x330x30

Model 6

Flat Size

B6-BW-03
245x250x30



Base (With Feet)

B1-BW-05
540x255x50

B2-BW-05
405x255x50

B3-BW-05
245x255x50

B4-BW-05
655x255x50

B5-BW-05
530x330x50

B6-BW-05
245x250x50

B1-BW-07
540x255x70

B2-BW-07
405x255x70

B3-BW-07
245x255x70

B4-BW-07
655x255x70

B5-BW-07
655x252x70

B1-BW-10
540x255x100

B2-BW-10
405x255x100

B3-BW-10
245x255x100

B4-BW-10
655x255x100

B5-BW-10
530x330x100



Lid

B1-BW-00
540x255

B2-BW-00
405x255

B3-BW-00
245x255

B4-BW-00
655x255

B5-BW-00
530x330

B6-BW-00
245x250

Model 7

Mini Size

Model 8

Dental Size

Model 9

1/2 Dental Size

Model 48

DIN Size

B9-BW-03
190x140x30

B48-BW-03
485x255x30



Base (With Feet)

B7-BW-05
270x125x50

B8-BW-05
270x170x50

B48-BW-05
485x255x50


B7-BW-07
270x125x70

B8-BW-07
270x170x70

B48-BW-07
485x255x70

B8-BW-10
270x170x100

B48-BW-10
485x255x100



Lid

B7-BW-00
270x125

B8-BW-00
270x170

B48-BW-00
485x255


Product Details

Base & Lid
Plain Weave



STANDARD

Fine Mesh Baskets



Base & Lid

B53-BLW-10
80x80x40

B53-BLW-11
145x150x50

B53-BLW-12
220x140x50


B53-BLW-13
275x175x35

B53-BLW-14
280x170x60



Base & Lid
Fine Wire Mesh

Surgical Laundry Baskets




Base

B55-BW-15
605x405x150

B55-BW-30
605x405x300

HTContainer

s t o r e




COMPONENTS

Container Components

54-55

Perforated Tray & Wire Basket Components

56-61



COMPONENTS

Container Components



Filter Retainer
(For Standard Type)

A1-FR-01



Filter Retainer
(For Mini&DentalContainers)

A1-FR-02



Filter Retainer
(For 1/2 Dental Container)

A1-FR-03



Security Seal

A1-S-03
100 Pcs.



Security Seal

A1-S-04
100 Pcs.



Security Seal
With Indicator

A1-SI-03
100 Pcs.



Security Seal
With Indicator

A1-SI-04
100 Pcs.

COMPONENTS

Container Components



Paper Filter
(For 1/1, 3/4, 1/2 Size)

A1-FP-01
100 Pcs.



Paper Filter
(For Mini&DentalContainers)

A1-FP-02
100 Pcs.



Paper Filter
(For 1/2 DentalContainer)

A1-FP-03
100 Pcs.



PTFE Filter
(For 1/1, 3/4, 1/2 Size)

A1-FT-01



PTFE Filter
(For Mini&DentalContainers)

A1-FT-02



PTFE Filter
(For 1/2 DentalContainer)

A1-FT-03



Identification Labels

- A1-LA-01
- A1-LA-02
- A1-LA-03
- A1-LA-04
- A1-LA-05
- A1-LA-06



Paper Labels
with indicators (Big)

A1-LP-01
100 Pcs.





Paper Labels
with indicators (Small)

A1-LP-02
100 Pcs.

COMPONENTS

Silicone Mats

Code	Dimensions mm
A1-SM-01	500x235
A1-SM-02	375x235
A1-SM-03	225x235
A1-SM-07	270X125



COMPONENTS

Silicone Holders



F1-FR-01

17x25x40 mm



F1-FR-02

17x100x40 mm



F1-FR-03

17x225x40 mm



F2-FR-01

17x25x40 mm



F2-FR-02

17x100x40 mm



F2-FR-03

17x225x40 mm



F3-FR-01

17x25x37mm

Ø12mm



F3-FR-04

17x25x37 mm

Ø14mm



F3-FR-02

17x25x30 mm

Ø12mm



F3-FR-05

17x25x30 mm

Ø14mm



F3-FR-03

17x24x21mm



Ø5mm





F3-FR-06

17x24x21 mm

Ø10mm



Type	Packing		Steam sterilization procedure acc. To DIN EN 868-8:2009-09				
	Description	Symbolic Delivery and Packaging System	Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure
With filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. To DIN EN 868-8:2009-09		+	+	+	+	+
			0	0	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated to DIN EN 868-8:2009-09		-	-	+	+	+
			-	-	+	+	0
	Textiles in sterilization container with filter in bottom, lid not perforated to DIN EN 868-8:2009-09		0	+	+	+	+
			-	0	+	+	+
	Instruments in wire basket acc. To DIN 58952 part III, in sterilization container with filter in lid and bottom acc. To DIN EN 868-8:2009-09		+	+	+	+	+
	Instruments in wire basket acc. To DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. To DIN EN 868-8:2009-09		-	+	+	+	+
	Instruments in wire basket acc. To DIN 58952 part III, in sterilization container with filter in bottom acc. To DIN EN 868-8:2009-09		0	+	+	+	+
							
With valve (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. To DIN EN 868-8:2009-09		-	0	+	+	0
			-	-	+	+	0
	Textiles in sterilization container with valves in lid, no valves on the bottom		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and bottom		-	+	+	+	+
	Instruments in wire basket in sterilization container with valve in lid no valves on the bottom		-	+	+	+	+
							



Honesty and Trust

HT Container Store Co. Inc.

Esenevler Mahallesi 310.Sokak No:43/6
55200 Atakum, Samsun, Türkiye
Tel. +90 537 445 16 42
htcontainerstore.com

Sterilization Container Systems

Edition - 01 2024 Version

Copyright © 2025, HT Container Store Co., Inc. All rights reserved. All visual and textual contents in this catalog are protected by copyrights and other rights owned by HT Container Store Co., Inc. We reserve the right to modify our products at any time.

htcontainerstore.com