

Lifeguard

Lifeguard

Guarding Lives Innovatively...



Lifeguard-ECG Holter

3 Channel ECG Holter



Features:

- Recording: 3 Channel ECG
- Compact/Lightweight (75*55*16mm, 50g (Battery not included))
- Sampling Rate : Up to 1024 Hz/16Bit (Sampling rate is up to 1024 samples/channel/second 16 bit resolution)
- Display : 1.3 inches OLED screen, 128*64 dot matrix
- Real-time wave forms display (Provide real-time OLED display of ECG waveform of all channels)
- Recording time: Up to 7 days
- Powered by 1 AAA alkaline batteries, up to 168 hours ECG recording
- 3 channel digital recording, automatic selection via patient cable
- Storage medium: Removable micro SD Card, ≥8 GB
- Data Transfer: Hi-Speed USB 2.0, SD card-reader
- Support multi-languages, easy and friendly to operate
- Pacemaker detection: Yes
- Sleep Management: yes

Technical Specifications

General

Weight	: 50g±5g (battery not Included)
Dimensions	: 70*55*16mm(±2mm)
Battery	: 1 AAA alkaline battery
Recording time	: 1-7 days optional
Storage medium	: Removable micro SD Card, ±8GB
Data Transfer	: Hi-Speed USB 2.0, SD Card-reader
Display	: 1.3 inches OLED screen 128*64 dot matrix

Marketed By

Lifeguard

Medical Systems Pvt. Ltd.

An ISO 13485:2016 Certified Co.

www.lifeguard-med.com

SCO-123, Honey View Market, Sector – 124, Mohali – 140301, Punjab, INDIA

Mobile: +91-7696257100

Mail ID: sales@lifeguard-med.com

Parameters

No. Of ECG Channels	: 3 Channels
A/d Resolution	: 16 bit
Sampling Rate	: 1024 samples/channel/second
Input Dynamic range	: $\pm 5\text{mV}$
Input Impedance	: $\geq 10\text{M}\Omega$
CMRR	: $\geq 80\text{dB}$
Offset voltage range	: $\pm 300\text{mV}$
Pulse width	: 0.1-2ms
Pulse amplitude	: 2-200mV
Frequency Response	: +0.4dB, -3.0dB (0.05-60Hz)
Event mark	: Yes

Accessories

Accessory	Quantity	Note
7 Lead Patient Cable (IEC)	1 Piece	Standard
Micro SD Card (8 GB)	1 Piece	Standard
AAA Battery	2 Piece	Standard
Carrying Bag & Strap	1 Piece	Standard
Recorder Operation Manual	1 Piece	Standard
Software Operation Manual	1 Piece	Standard

Marketed By

Lifeguard

Medical Systems Pvt. Ltd.

An ISO 13485:2016 Certified Co.

www.lifeguard-med.com

SCO-123, Honey View Market, Sector – 124, Mohali – 140301, Punjab, INDIA

Mobile: +91-7696257100

Mail ID: sales@lifeguard-med.com