

Lifeguard

Lifeguard ECG Patch Holter is a comfortable, medical-grade wearable solution for continuous and reliable cardiac rhythm monitoring.

INTRODUCING



Lifeguard Medical Systems Pvt Ltd is the Manufacturer of Medical Equipments based in Chandigarh, India. Lifeguard medical Systems is ISO Certified 13485:2016 with Manufacturing License. We are also one of India's Fastest growing companies in Healthcare Sector. Dedicated towards bringing latest in Hospital and home care health technology to our customers at the most competent prices. Our aim is to give our customers full value for their money by bringing to them only the most accurate and reliable health equipments.

SCO 123, Honey View Market, Sector 124, Mohali - 140301 (Punjab)
INDIA

Mobile : 7696257100

Email @ : sales@lifeguard-med.com

www.lifeguard-med.com

Lifeguard

GUARDING LIVES INNOVATIVELY!!!



ISO 13485: 2016 Certified



Range of Products

Lifeguard Mfg Licence No. MFG/MD/2023/000533



Lifeguard - 12 Channel ECG Machine

Specifications

- Recording: 12 Channel ECG
- Full keyboard, 7" touch screen operation
- Sampling rate up to 64000Hz, pacemaker detection
- Save up to 1000 ECG records internally (expandable with SD card and USB flash-disk)
- Advanced ECG analysis: VCG/TVCG/VLP/HRV
- Support external USB printer, Keyboard, Mouse, Bar-code scanner etc.
- Data transfer to PC directly by Ethernet, USB or Wi-Fi (optional)



Lifeguard - 3 Channel ECG Machine

Specifications

- 7" color touch screen
- Automatic measurement, calculation, analysis, waveform freezing.
- ² Thermal printer with 80mm print paper, synchronization print.
- ² Built-in rechargeable li-ion battery (11.1V/4000mAh), AC/DC power conversion.
- Adapt to 100-240V, 50/60Hz AC power supply
- Support USB flash disk stored reports can be saved, opened, analyzed in PC via ECG playback software.
- The machine can store over 400 ECG reports in its built-in flash.



Lifeguard - 3 Channel Animal ECG Machine

Specifications

- 7" Color Touch Screen
- Inbuilt Thermal Printer 80mm.
- The machine can store over 400 ECG reports in its built-in flash.
- Built-in Rechargeable Li-ion Battery.
- Simultaneously acquisition of 12 lead.



Lifeguard - 3 Channel ECG Holter

Specifications

- Recording: 3 Channel ECG.
- Compact/Lightweight (75*55*16mm, 50g (Battery not included))
- Sampling Rate : Up to 1024 Hz/16Bit.
- Display : 1.3 inches OLED screen, 128*64 dot matrix.
- Real-time wave forms display.
- Recording time: Up to 7 days.
- Powered by 1 AAA alkaline batteries, up to 168 hours ECG recording.
- Storage medium: Removable micro SD Card.
- Data Transfer: Hi-Speed USB 2.0, SD card-reader.
- Pacemaker detection: Yes.



Lifeguard - 12 Channel ECG Holter

Specifications

- Recording: 12 Channel ECG.
- Compact/Lightweight (75*55*16mm, 50g (Battery not included))
- Sampling Rate : Up to 1024 Hz/16Bit.
- Display : 1.3 inches OLED screen, 128*64 dot matrix.
- Real-time wave forms display.
- Recording time: Up to 4 days.
- Powered by 1 AAA alkaline batteries, up to 72/96 hours ECG recording.
- Storage medium: Removable micro SD Card.
- Data Transfer: Hi-Speed USB 2.0, SD card-reader.
- Pacemaker detection: Yes.



Lifeguard - SpiroLife

Specifications

- USB Compatible Handle
- Lung Age Calculation
- Graph and Grid color selection
- PDF Report Generation
- 3 Manuever Result Display
- Auto Interpretation in Result