

Lifeguard

GUARDING LIVES INNOVATIVELY!!!



Since 2010....



Lifeguard Mfg Licence No. MFG/MD/2023/000533

ISO 13485: 2016 Certified


 sales@lifeguard-med.com

 +91-7696257100 | +91-9316606410 | +91-7814900743

 www.lifeguard-med.com

12. CHANNEL ECG MACHINE

- 12 Channel ECG Recording.
- Full keyboard, 7" touch screen operation.
- Sampling rate up to 64000Hz, pacemaker detection.
- Save up to 1000 ECG records internally (expandable with SD card and USB flash-disk).
- Advanced ECG analysis VCG/TVCG/VLP/HRV.
- Support external USB printer, Keyboard, Mouse, Bar-code scanner (optional).
- Data transfer to PC directly by Ethernet, USB or Wi-Fi (optional).
- Dual battery option.



3 CHANNEL ANIMAL ECG MACHINE

- 7" Color Touch Screen.
- Inbuilt Thermal Printer 80mm.
- The machine can store over 400 ECG reports in its built-in flash.
- Built-in Rechargeable Li-ion Battery.
- Simultaneously acquisition of 12 lead.

3 CHANNEL ECG MACHINE

- 7" color touch screen.
- Automatic measurement, calculation, analysis, waveform freezing.
- Built-in Thermal Printer (80mm).
- Built-in rechargeable li-ion battery, AC/DC power conversion.
- Adapt to 100-240V, 50/60Hz AC power supply.
- Support USB flash disk stored reports can be saved, opened, analyzed in PC via ECG playback software.
- The machine can store over 400 ECG reports in its built-in flash.



3 CHANNEL ECG HOLTER

- 3 Channel ECG Recording.
- Compact & Lightweight- 75*55*16mm, 50g (Battery not included).
- Sampling Rate: Up to 1024 Hz/16Bit.
- 1.3 inches OLED screen, 128*64 dot matrix.
- Recording time: Up to 7 days.
- Powered by 1 AAA alkaline batteries.




12 CHANNEL ECG HOLTER

- 12 Channel ECG Recording.
- Compact & Lightweight- 75*55*16mm, 50g (Battery not included).
- Sampling Rate: Up to 1024 Hz/16Bit.
- 1.3 inches OLED screen, 128*64 dot matrix.
- Recording time: Up to 4 days.
- Powered by 1 AAA alkaline batteries.

SINGLE CHANNEL ECG PATCH HOLTER

- Single Channel ECG Recording.
- Ultra-Lightweight (22 Grams).
- Recording Time up to 11 days.
- Pacemaker Detection.
- Water & splash Resistant.
- Accurate AF Screening.
- Long Battery Life.
- Biocompatible Material & Disposable Patch.



3 CHANNEL ECG PATCH HOLTER

- 3 Channel ECG Recording.
- Ultra-Lightweight (22 Grams).
- Recording Time up to 07 days.
- Pacemaker Detection.
- Water & splash Resistant.
- Accurate AF Screening.
- Long Battery Life.
- Biocompatible Material & Disposable Patch.

12 CHANNEL ECG PATCH-HOLTER

- 12 Channel ECG Recording.
- Ultra-Lightweight (22 Grams).
- Recording Time up to 05 days.
- Hi-Speed USB Data Transfer.
- Pacemaker Detection.
- Water & splash Resistant.
- Accurate AF Screening.
- Long Battery Life.
- Biocompatible Material & Disposable Patch.



SPIROMETER

- USB Compatible Handle.
- Lung Age Calculation.
- Graph and Grid color selection.
- PDF Report Generation.
- 3 Manuever Result Display.
- Auto Interpretation in Result.

Scan Here
Sample ECG Holter &
Spirometer Report



Single-Channel



3 - Channel



12 - Channel



Spirometer

PC ECG MACHINE

- Portable USB-powered 12-lead ECG device.
- Compact, lightweight, and user-friendly design.
- Real-time cardiac monitoring.
- Easy connectivity with PC/Laptop and mobile devices via OTG support.
- No separate battery or charging required.
- Real-time ECG waveform display and analysis.
- Quick sharing via WhatsApp, email, and other platforms.
- Ground isolation technology for enhanced patient safety.
- Ideal for hospitals, clinics, and diagnostic centers.



ENQUIRE NOW

📍 SCO 123, Honey View Market, Sector 124,
Mohali-140301 Punjab, India

🌐 www.lifeguard-med.com

✉ sales@lifeguard-med.com

☎ +91-7696257100 | +91-9316606410 | +91-7814900743

📷 <https://www.instagram.com/lifeguard962>