



SampleRight™

Sponge Samplers

For Environmental Monitoring



Resists tearing & crumbling
associated with cellulose sponge

Increased rigidity
makes it easier to handle & manipulate

Polyurethane sponge
medical grade, biocide-free

Allows for vigorous scrubbing
of surfaces to disrupt biofilm

SampleRight™

Upgrade to a innovative polyurethane sponge
for easier handling and better sampling

The SampleRight™ Sponge Sampler with a high-density polyurethane sponge provides another innovative solution to solve everyday surface sampling challenges. The sponge is easy to manipulate when moving it in and out of the sample bag and easy to hold while collecting a surface sample. The polyurethane material resists tearing and flaking, making it ideal for rough surfaces. Each unit is hydrated with a 10 ml of buffer, providing sufficient volume of liquid to perform multiply quantitative tests for microorganisms from a single sample. Select the dry version and add your choice of liquid for a customized option.



Polyurethane advantages

- Resists crumbling and tearing
- Biocide free material
- Better recovery of organisms from surfaces
- Better release of organisms from the sponge
- More consistent solution volume for testing

Choose the Right Device

Match the sponge to the application

Available dry or pre-hydrated with various collection solutions

NO gloves



WITH gloves



HiCap™ Neutralizing Broth

An advance sample collection solution designed specifically for food safety, that improves the recovery of microorganism from surfaces and unlike other collection solutions, it maintains the viability of stressed and injured cells for [up to 72 hours after collection](#).



Online Resource Library

Learn how to improve the effectiveness of your environmental monitoring programs. Visit our online Education Library to download white papers and studies, or view recorded presentations detailing the efficacy and advantages of HiCap Neutralizing Broth.

worldbioproducts.com/education

Beyond Basic Sample Collection

We design tools to support an **Environmental Monitoring Program** that meet today's unique food safety challenges



EZ Reach™ Sponges

Environmental monitoring requires collecting samples from a wide variety of surfaces, from processing equipment and conveyors, to drains and expansive floor or wall surfaces. The EZ Reach™ sponge sampler family offers a range of sizes and functions to fit most sampling needs. All EZ Reach sponges combine a high purity polyurethane sponge with our patented twist-off handle.

US Patent 8,127,627 B2



Dilutions & Enrichments

Improve laboratory efficiencies with prepared diluents and media. FlipRight™ prefilled vials are available in ready-to-use 9 ml and 10 ml fill volumes for serial dilutions or hydrating sponges. NovaLock™ bottles come with ready-to-use 90 ml or 99 ml fill volumes for convenient sample preparation.

US Patent 8590 735 B2



Dividable Grid Containers

Safely and efficiently move samples around the laboratory. Reconfigurable dividers provide clear separation of samples, allowing for easy handling and processing while minimizing the risk of cross contamination.

www.worldbioproducts.com
orders@worldbioproducts.net
PO Box 947, Bothell, WA 98041-0947
(877) 260-6441

WBMKT_BROSREN2102

