

# DATASHEET CLEAN FLASH OPA

Rapid and Effective Disinfection for  
High-Risk Medical Instruments and  
Endoscopes

## CHEMICAL COMPOSITION

Ortho-phthalaldehyde%0,55, Corrosion inhibitors and auxiliary agents.

## PRODUCT SPECIFICATIONS

Clean Flash OPA is ready-to-use and disinfects and sterilize in 5 minutes and usable for maximum 14 days without loss of efficiency. It contains corrosion inhibitor. It can be easily checked with test strips. It is a ready-to use, colourless, odour-free and nontoxic disinfectant for disinfection of medical equipment, flexible and rigid endoscopes and it is phenol-free and formaldehyde-free.

## AREAS OF USE

Clean Flash Opa Solution is very suitable for disinfection of any type of surgical equipment, heatsensitiveanaesthetic accessories (ambu mask, etc.) and all endoscopes such as Karl Storz, Olympus, Pentax, Fujinon, flexible and rigid) and for disinfection of all kinds of surgical equipment and laparoscopy equipment. It does not cause any harm on the endoscope covers such as rust, corrosion, softening, colour change, etc.

## ANTI-MICROBIAL CHARACTERISTICS

Bactericidal, Fungicidal, Sporocidal, Virucidal, Tuberculocidal, Micobactericidal.



+970562600007



www.activemed.ps



info@activemed.ps

