

DATASHEET

ULTRASOUND GEL

Ensuring Clear and Uninterrupted Images
for Ultrasound and Medical Applications

PRODUCT SPECIFICATIONS

In all kinds of ultrasonography and doppler applications, it cuts the air between the skin and the probe and ensures that the ultrasonography waves come to the device screen more clearly and uninterruptedly.

AREAS OF USE

Used in many areas like Ultrasound, ECG, Cardiac stress test, ECG electrodes, Epilation applications, etc. It does not contain oil and oily substances. It does not contain formal- dehyde and salt. Easy to clean.

STORAGE CONDITIONS

Store in a cool place between 5 - 30 C and away from sunlight.

