

VivaChek



BioSieve™

Multi-Drug Urine Test Cup


- FDA 510K Cleared
- CLIA Waived
- OTC Approved

BioSieve™

Multi-Drug Urine Test Cup

Test Item

CLIA-waived	Abbreviation	Cut-off (ng/mL)
Amphetamine	AMP	500, 1000
Secobarbital	BAR	300
Buprenorphine	BUP	10
Oxazepam	BZO	300
Cocaine	COC	150, 300
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	EDDP	300
Methylenedioxyamphetamine	MDMA	500
Methamphetamine	MET	500, 1000
Morphine	MOP	300
Methadone	MTD	300
Opiate	OPI	2000
Oxycodone	OXY	100
Phencyclidine	PCP	25
Propoxyphene	PPX	300
Nortriptyline	TCA	1000
Marijuana	THC	50



Specification

Product Catalog	Popular Products					Other Options
Number of Panel	5	10	12	14	15	2-15
Test Items	AMP, COC, MET, MOP/OPI, THC	AMP, BAR, BZO, COC, MET, MOP/OPI, MTD, OXY, PCP, THC	AMP, BAR, BZO, COC, MET, MDMA, MOP/OPI, MTD, OXY, PCP, TCA, THC	AMP, BUP, BAR, BZO, COC, MET, MDMA, MOP/OPI, MTD, OXY, PCP, PPX, TCA, THC	AMP, BUP, BAR, BZO, COC, EDDP, MET, MDMA, MOP/OPI, MTD, OXY, PCP, PPX, TCA, THC	Any combination of 2-15 items

VivaChek

Multi-Drug Urine Test Cup




- All-In-One Cup, screen for up to 15 drugs
- Convenient & Reliable
- Results within 5 minutes




VivaChek

BioSieve™

Test Procedure



Test Results



VivaChek
VivaChek Diagnostics, Inc.
1401 21st St. STE R, Sacramento, CA
<https://vivachekdx.com>

Number: 1626008201
Effective Date: 2023-09-02

VivaChek