



**BioSieve™** / **VivaChek**

**BioSieve™**

**Multi-Drug Test**

# About VivaChek



The Global leader in Blood Glucose Meter Innovation and Manufacturing is creating its base in the U.S with launch of VivaChek USA based in Sacramento, CA. VivaChek specializes in OEM, ODM and Custom Manufacturing of various Drugs of Abuse Tests in various standard and tailor-made configurations.

We aspire to become a trusted partner in manufacturing for the U.S to cater the Professional and OTC market for DOA products with strategically placed warehouses in the U.S.

## Milestones



○ 2023  
VivaChek Diagnostics Inc. USA Established.  
Products Available in 120+ countries,  
150+ patents

○ 2022  
DOA urine cup of 15 Drugs: 510(k) cleared,  
CLIA Waived & OTC approved

○ 2021  
Revenue reached 180+ million USD

○ 2019  
Blood Ketone+ BGM dual system  
510(k) cleared

○ 2017  
1<sup>st</sup> BGM 510(k) cleared

○ 2013  
VivaChek established



# Urine Drug Test

## 510 (k) Cleared, CLIA Waived & OTC Approved

- Multi-Drug Urine Test Cup/Panel
- Marijuana Single Dip Panel/Strip

- ✓ Any combination up to 15 drugs
- ✓ Convenient & Reliable
- ✓ Results within 5 minutes
- ✓ Long expiration dating

## Product Catalogue


Urine Drug Test - 510 (k) Cleared, CLIA Waived & OTC Approved — 4

- Multi-Drug Urine Test Cup — 5
- Multi-Drug Urine Test Panel — 6
- Marijuana Single Dip Panel/Strip — 7

Urine Drug Test – FUO — 8

Oral Fluid Drug Test – FUO — 9

Accessories & Related Products — 10



# 01 Multi-Drug Urine Test Cup

- 510 (k) Cleared, CLIA Waived & OTC Approved

## Products Information

Drug Test	Cut-off (ng/mL)
2-ethylidene-1,5dimethyl-3,3-diphenylpyrrolidine (EDDP)	300
Amphetamine(AMP)	500,1000
Buprenorphine(BUP)	10
Cocaine (COC)	150, 300
Marijuana (THC)	50
Methylenedioxyamphetamine (MDMA)	500
Methadone(MTD)	300
Methamphetamine(MET)	500, 1000
Morphine(MOP)	300
Morphine(OPI)	2000
Nortriptyline (TCA)	1000
Oxazepam(BZO)	300
Opiate (OPI)	2000
Oxycodone (OXY)	100
Phencyclidine (PCP)	25
Propoxyphene(PPX)	300
Secobarbital (BAR)	300




# 02 Multi-Drug Urine Test Panel

- 510 (k) Cleared, CLIA Waived & OTC Approved

## Products Information

Drug Test	Cut-off (ng/mL)
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300
Amphetamine (AMP)	500, 1000
Barbiturates (BAR)	300
Buprenorphine (BUP)	10
Cocaine (COC)	150, 300
Ecstasy (MDMA)	500
Marijuana (THC)	50
Methadone (MTD)	300
Methamphetamine (MET)	500, 1000
Morphine (MOP/OPI )	300, 2000
Nortriptyline (TCA)	1000
Oxazepam (BZO)	300
Oxycodone (OXY )	100
Phencyclidine (PCP)	25
Propoxyphene (PPX)	300





# 03 Marijuana Single Dip Panel/Strip

## Product Information

### Drug Test

Marijuana (THC)

### Cut-off (ng/mL)

20 - 510 (k) Cleared, CLIA Moderate

50 - 510 (k) Cleared, CLIA Waived & OTC Approved



# 04 Urine Drug Test - FUO




## Products Information

Drug Test	Cut-off (ng/mL)
6-MAM	10
7-Aminoclonazepam (7-ACL)	100, 200, 300
AB-Pinaca (ABP)	10
Acetaminophen (ACE)	5000
Alprazolam (ALP)	100
Caffeine (CAF)	1000
Cathine (CAT)	150
Carfentanyl (CFYL)	500
Clonazepam (CLO)	150, 400
Cannabinol (CNB)	500
Cotinine (COT)	10, 50, 100, 200
Diazepam (DIA)	300
Ethyl Glucuronide (EtG)	500, 1000
Fentanyl (FYL)	10, 20
Gabapentin (GAB)	2000

γ-Hydroxybutyric acid (GHB)	10 µg/mL
Synthetic Marijuana (K2)	25, 30, 50
Ketamine (KET)	100, 300, 500, 1000
LSD	20
Methcathinone (MCAT)	500
Tenamfetamine (MDA)	500
MDPV	1000
Mephedrone HCl (MEP)	100
Methylphenidate (MPD)	300, 1000
Meperidine (MPRD)	100
Methaqualone (MQL)	300
Pregabalin (PGB)	50000
Tramadol (TML)	100,300
Tropicamide (TRO)	350
Trazodone (TZD)	200
Urine Alcohol	0.02%
Xylazine (XYL)	500
Zolpidem (ZOL)	50
α-Pyrrolidinovalerophenone (α-PVP)	300, 500, 1000

# 05 Oral Fluid Drug Test - FOUO



## Products Information

Drug Test	Cut-off (ng/mL)		
AB-Pinaca (ABP)	10	Ketamine (KET)	1000
Amphetamine (AMP)	300, 500	MDPV	1000
Benzodiazepines (BZO)	100, 200	Methamphetamine (MET)	300, 500
Cannabinoids (THC)	25, 50	Methadone (MTD)	200, 300
Clonazepam (CLO)	100	Methaqualone (MQL)	300
Cocaine (COC)	100, 150	Methylphenidate (MPD)	300
Cotinine (COT)	200	Monoacetylmorphine (6-MAM)	10
EDDP	100	Morphine (MOP/OPI)	100
Ethyl Glucuronide (EtG)	500	Pregabalin (PGB)	1000
Fentanyl (FYL)	200	Synthetic Cannabinoid (K2)	50
		Tramadol (TML)	100, 200, 300

# 06 Accessories & Related Products



DOA FIA Test Kit \*

\* Coming soon

# 07 FIA Testing System



BioSieve™ FIA Reader

## Products Information

Analyzer	Specification
Size	124.5 × 72.5 × 40 mm
Weight	165 g
Screen	1.54" LCD
Connectivity	Bluetooth
Power	Battery & AC/DC adapter
Printer	External bluetooth printer (optional)
Memory	999 results