



FQbeauty Group
Forever innovation

HK FT BEAUTY EQUIPMENT CO.,LIMITED
(Guangzhou Fuqiang Beauty Equipment Co.,Limited)

Add:No. 103, Huasheng South Road, Donghua
Industrial Park, Renhe Town, Baiyun District,
Guangzhou, Guangdong, China (Mainland)

Experienced manufacturer

Our team has specialized in beauty industrial for more than 16 years.

Main product:

Diode laser machine, Fractional co2 laser machine, cryolipolysis machine, nd yag laser machine, 980nm laser machine, hifu machine, opt shr machine, cavitation rf machine, Emsculpture



Our team



Production room



OEM/ODM service

We provide a complete solution from software, hardware, body, structure design for the customized clients need OEM/ODM CE Certificate

With us your business is safe
and your money is safe.

G22(FQ808)

Feature:

Screen size: 10.4 inch HD display

Air outlet: Water cooling plus air cooling

Spot size: 15mm*40mm

Voltage: AC110V-220V 50-60Hz

Instrument size: 130cm*50cm*40cm

Handle: two

Material: sheet metal


Operating position: A B C D E

Frequency: 1-10/hz

Energy: 1-50/jm

Application:

Painless hair removal combines skin rejuvenation help skin smooth and light



Noninvasive 532

532 wavelength is aimed at light color and color, red mouth, yellow, green and light coffee. Light blue and red are difficult to wash and can only be diluted, mainly washing tattoos.



532

All very light colors and colors choose 532nm wavelength to wash light color tattoos, light red eyebrows and light blue hair.



Black faced doll

Black faced dolls have thick cuticle, dark yellow complexion, acne marks, large pores and strong oil secretion.



LP1064/532

The buckle reset corresponds to 1064 wavelength. For dark dark black, dark gray, dark blue, that is, freckle spots and nevus: fade the epidermis, pull out the buckle, corresponding to anti yellowing and rejuvenation.



Noninvasive 1064

1064 wavelength is for non-invasive eyebrow washing (dark black, dark gray, dark blue, dark coffee).




Diode / Pico laser machine (FQ808P)

Feature:

Display: 15.6" colorful touch screen
 Laser wavelength: 755nm+808nm+1064nm
 Spot size: 12*40mm
 Pico: 5 cartridges, (532/755/1064nm)
 Energy: 1-2000
 Frequency: 1-10Hz
 Power supply: 220V, 50Hz or 110V, 60Hz
 Packing size: 65.5*59*131.5cm
 Gross Weight: 68Kg

Application:

Full body Hair remove, skin rejuvenation



3 in 1 all skin type can be used



- 755 nm for white skin (fine, golden hair removal);
- 808 nm for yellow neutral skin (brown hair removal);
- 1064 nm for black skin (black hair removal)



808 semiconductor hair removal(FQ094)

Feature:

Spot size of hand-piece: 12*20mm

Voltage: 110V-220V / 50-60Hz

Energy: 1-30j / cm


Frequency: 1hz-10hz

Screen size: 10.4 Inch Touch Screen

Instrument size: 40 * 50 * 130cm

Application:

Noninvasive eyebrow washing
Freezing point depilation
Photon rejuvenation



3 wavelength DIODE LASER MACHINE(FQ808E)

Feature:

Screen size: 15 inch HD display
 Air outlet: Water cooling plus air cooling
 Spot size: 15mm*40mm
 Voltage: Ac110v-220v50-60hz
 Instrument size: 133cm*37cm*38cm
 Handle: two
 Material: sheet metal
 Operating position: S M L XL
 Skin color : I II III IV V VI

Application:



Fast hair removal, safe and painless, longer lasting, Three simultaneous working effects take effect at the same time



808 nm



1064+755+808nm



OUR PRODUCT APPEARANCE



808 Semiconductor Laser Hair Removal Instrument (FQ066)

Feature:

Screen size: 15.6 inch HD display

Air outlet: Water cooling plus air cooling

Spot size: 15mm*40mm

Voltage: AC110V-220V 50-60Hz

Package Size: 130*57*60

Weight: 88kg

Handle: two

Material: ABS

Application:

Hair removal, photon rejuvenation, non-invasive eyebrow washing, tattoo washing, fade pigment, make black faced dolls, point moles



NON-INVASIVE ENERGY GATHERING HAND TOOLS



808 HANDLE



ONE MACHINE TO DEAL WITH MULTIPLE SKIN PROBLEMS



LP1064/532

The built-in test corresponds to 1064 wavelength. For dark skin (black, dark gray, dark blue, etc.), freckle spots and nevus, fade the epidermis, pull out the built-in, corresponding to anti-yellowing and rejuvenation.



Noninvasive 1064

1064 wavelength is for non-invasive eyebrow washing (dark black, dark gray, dark blue, dark coffee).



Noninvasive 532

532 wavelength is aimed at light color and color, red mouth, yellow, green and light coffee. Light blue and red are difficult to wash and can only be diluted, mainly washing tattoos.



532

All very light colors and colors choose 532nm wavelength to wash light color tattoos, light red eyebrows and light blue hair.



Black faced doll

Black faced dolls have thick cuticle, dark yellow complexion, acne marks, large pores and strong oil secretion.



Photon rejuvenation



arm hair removal



Method for removing spots



eyebrow tattoo



tattoo removal



black face doll

2 in 1 Binary Star 808 diode laser hair removal machine(FQ096)

Feature:

Screen size: 10 inch HD display
 Air outlet: Water cooling plus air cooling
 Spot size: 15mm*40mm
 Voltage: AC110V-220V 50-60Hz
 Package Size: 124*73*50 weight: 90kg
 Handle: two
 Material: ABS

Application:

Hair removal, photon rejuvenation, non-invasive eyebrow washing, tattoo washing, fade pigment, make black faced dolls



Non-invasive energy gathering hand tools
 (Standard configuration 5 heads)




808nm



SKIN PROBLEMS

One machine to deal with multiple skin problems



LP1064/532

The bubble reset corresponds to 1064 wavelength.
 For dark dark black, dark gray, dark blue, that is, freckle spots and nevus: fade the epidermis, pull out the bubble, corresponding to anti-yellowing and rejuvenation




Noninvasive 1064

1064 wavelength is for non-invasive eyebrow washing (dark black, dark gray, dark blue, dark coffee).



Noninvasive 532

532 wavelength is aimed at light color and color, red mouth, yellow, green and light coffee. Light blue and red are difficult to wash and can only be diluted, mainly washing tattoos.



Diode laser with pico laser machine (FQG3000)

Feature:

Display: 10.4" colorful touch screen

Diode laser: 808nm

Wavelength: 532/1064/1320nm

Spot size: 15*15mm

Cooling way: Water circulation system + fan cooling

Net weight: 50kg

Packing size: 125*72*64cm

Working time: More than 12 hours

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, tattoo remove, freckle remove, spot remove



Before

After



Black Baby tip



Freckle



Wash eyebrows



Tattoo removal



Before

After

Diode laser with pico laser machine (FQ808G)

Feature:

Display: 8.4" colorful touch screen

Diode laser: 808nm

Wavelength: 755/808/1064nm

Spot size: 15*15mm

Cooling way: Water circulation system + fan cooling

Net weight: 45kg

Packing size: 123*72*64cm

Working time: More than 12 hours

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, tattoo remove, freckle remove, spot remove

808NM DIODE LASER PICO LASER

2in1 beauty machine

- 01 Fast cooling
- 02 Short focus
- 03 Strong power
- 04 Effective treatment



Before

After

Before

After

808 semiconductor hair removal(FQ063)

Feature:


Voltage: ac10v-220v50-60hz
 screen: (10 inch color touch screen)
 wavelength non-invasive
 5321064nmlp1064 / 532nm
 black doll: 532nms, point size
 10X10mm, 1515mm, 15x20
 frequency: 01-10hz,
 working time 10ms-500ms
 energy density: 01-135j / m2

Application:

Painless hair removal combineskin
 rejuvenation help skinsmooth and light



The non-invasive and normal modes can be adjusted.



808 semiconductor hair removal(FQ065)

Feature:

Voltage: ac10v-220v50-60hz
 screen: (10 inch color touch screen)
 wavelength non-invasive
 5321064nmlp1064 / 532nm
 black doll: 532nms, point size
 10X10mm, 1515mm, 15x20
 frequency: 01-10hz,
 working time 10ms-500ms
 energy density: 01-135j / m2

Application:

Painless hair removal combineskin
 rejuvenation help skinsmooth and light



HAIR REMOVAL EFFECT

808nm for all skin types



HAIR REMOVAL HANDLE

Painless/Comfortable/Safe/Permanent



808 semiconductor hair removal(FQ093)

Feature:

Screen size: 10 inch HD display
 Air outlet: Water cooling plus air cooling
 Language: Chinese, English, Japanese
 Voltage: AC110V-220V 50-60Hz
 Instrument size: 135cm*55cm*55cm
 Handle: four
 Material: sheet metal
 Package size: 131*62*60
 Weight: 40kg

Application:

Painless hair removal combines skin rejuvenation help skin smooth and light



Vertical 810nm diode laser machine (FQ808E)

Feature:

Display:8.4" colorful touch screen

Wavelength:808nm

Spot size:12*12mm

Shoot:5million

Laser power:2000w

Pulse width:10-400ms

Frequency:1-10hz

Cooling way:Water circulation +Air

Working time:More than 12 hours

Power supply: 220V, 50Hz or 110V, 60Hz

Gross weight:50KGS

Packing size:124*50*61cm

Application:

Full body Hair remove,skin rejuvenation



Portab; e Diode laser (FQ808P)

Feature:

Display:10.4" colorful touch screen

Laser wavelength:755nm+808nm+1064nm

Spot size:12*40mm

Frequency:1-10Hz

Power supply: 220V, 50Hz or 110V, 60Hz

Packing size:130*52*62cm

Gross Weight:50Kg

Application:

Full body Hair remove,skin rejuvenation



808 semiconductor hair removal (FQ082)

Feature:

Screen size: 10 inch HD display
 Air outlet: Water cooling plus air cooling
 Language: Chinese, English
 Voltage: AC110V-220V 50-60Hz
 Package Size: 67*58*50
 weight: 40kg
 Handle: four
 Material: sheet metal

Application:


Painless hair removal combines skin
 rejuvenation help skin smooth and light



face



Hand

Key Switches and
Emergency SwitchesHandle and
rotatable screen

Bikini



Legs



Extractor fan



Power cord interface



808 semiconductor hair removal(FQ808)

Feature:

voltage: 110V-220V / 50-60hz

808nm laser flow: 001-030j / cm

Energy: 40-100%


working time: 10-100msf

frequency: 1-10hz

scope of application: legs, hands,
armpits and chest

Application:

Painless hair removal combineskin
rejuvenation help skinsmooth and light



OPT SHR Nd Yag Laser Rf Machine (FQA31-4)

Feature:

Display: 8.4" colorful touch screen
 Wavelength: 530/590/640nm
 Nd yag laser: 532/1064/black doll
 Nd yag laser frequency: 1-10
 Energy: 1-50J/cm²
 Elight Discharge time: 1-4s, Elight plus: 1-9ms
 Elight Pulse repetition rate: 1-6 P/s
 OPT SHR plus: 1-8ms
 OPT SHR frequency: 1-5s
 Power: 3000W
 Package size: 56*64*131CM
 Continuous working time: 10 hours
 Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy,
 Vascular therapy, Acne therapy



8 Available language

中文 中文繁体

English Español

Русский Português

Việt 한국어



Multifunction 360 magneto hair removal machine (FQA31-6)

Feature:

Elight Display: 8.4" colorful touch screen
 Power supply: 220V, 50Hz or 110V, 60Hz
 Paching size: 64*56*131cm
 Energy: 1-50j
 Spot size: 15*40mm²
 First plug width: 01-09ms
 Second plug width: 01-09ms
 First plug delay: 5-90ms
 Pulse number: 1-6P/s
 Plus width: 1-8ms

Application:

Hair remove, skin rejuvenation, Pigmentation therapy,
 Vascular therapy, Acne therapy, tattoo remove, skin
 tightening



Multifunction 360 magneto hair removal machine (FQA31-2)

Feature:

Elight Display: 8.4" colorful touch screen
 Laser Display: 4.3" colorful touch screen
 Handle: 4 handle (HR, SR, RF, Nd Yag laser)
 Nd Yag laser: 532/1064/black doll
 RF handle: 3 different size probe
 Laser energy: 500-2000mj; Laser frequency: 1-10HZ
 SHR Energy: 1-50J/cm²; SHR frequency: 1-6Hz
 Elight Discharge time: 1-4s, Elight plus: 1-9ms
 Elight Pulse repetition rate: 1-6 P/s
 SHR plus: 1-8ms; Power: 3000W
 Packing size: 131*56*64cm; GW: 65KG
 Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy,
 Vascular therapy, Acne therapy, tattoo remove, skin
 tightening



8 Available language

中文 中文繁体

English Español

Русский Português

Việt 한국어





Multifunction 360 magneto hair removal machine (FQA31-5)

Feature:

Display: 10" colorful touch screen
 Handle: 4 handle (HR, SR, RF, Nd Yag laser)
 Nd Yag laser: 532/1064/black doll
 RF handle: 3 different size probe
 Laser energy: 500-2000mj; Laser frequency: 1-10HZ
 SHR Energy: 1-50J/cm²; SHR frequency: 1-6Hz
 Elight Discharge time: 1-4s, Elight plus: 1-9ms
 Elight Pulse repetition rate: 1-6 P/s
 SHR plus: 1-8ms; Power: 3000W
 Packing size: 131*56*64cm; GW: 65KG
 Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy,
 Vascular therapy, Acne therapy, tattoo remove, skin
 tightening



Pico laser tattoo laser hair remove machine (FQA32-6)

Feature:

Display: 10" colorful touch screen

Handle: 2 hand (HR+Pico laser)

Nd yag laser: 532/1064/black doll

Size: 15*40mm²

Energy: 1-50J/cm²

Plus 1-9ms (adjustable)

Fast model plus: 1-9ms (adjustable)

Frequency: 1-10Hz

Power: 3000W

64*56*13154*74*111cm

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy, Vascular therapy, Acne therapy, tattoo remove, skin tightening



Black doll



532nm



1064nm



Support the design of dual-head output, bald automatically select the bald independent count design, the operation of the handle at the same time with the light and the wave of two light filter device, three kinds of filter cover a variety of light therapy range



before



after



before



after



before



after



before



after



before



after



before



after

Nd yag laser tattoo laser hair remove machine (FQA32-3)

Feature:

Display: 8" colorful touch screen

Handle: 2 hand (HR+nd yag laser)

Nd yag laser: 532/1064/black doll

Size: 15*40mm²

Energy: 1-50J/cm²

Plus 1-9ms (adjustable)

Fast model plus: 1-9ms (adjustable)

Frequency: 1-10Hz


Power: 3000W

Packing size: 54*74*111cm

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy, Vascular therapy, Acne therapy, tattoo remove, skin tightening



360 Magneto optic hair removal machine (FQAD5-3)

Feature:

Display: 8" colorful touch screen

Power: 2000W

Shoot: 400000

Size: 15*40mm²

Energy: 1-50J/cm²

Interval time: 1-4s

Pulse repetition rate: 1-6 P/s

Plus time: 1-9ms

Frequency: 1-6Hz

Packing size: 62*51*61CM

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair removal range: armpit hair, hand hair, leg hair, bikini, hair, hairline, back hair



Hair remove single mode



Hair remove fast mode



360 Magneto optic hair removal machine (FQA20-3)

Feature:

Display: 8" colorful touch screen

Power: 2000W

Shoot: 400000

Handle: 1(HR)

Size: 15*40mm²

Energy: 1-50J/cm²

Interval time: 1-4s

Pulse repetition rate: 1-6 P/s

Pulse time: 1-9ms

Frequency: 1-6Hz

Packing size: 115*56*56CM

Application:

Hair removal range: armpit hair, hand hair, leg hair, bikini hair, hair, hairline, back hair



Hair remove single mode



Hair remove fast mode



360 Magneto optic hair removal machine (FQAD5-1)

Feature:

Display: 8" colorful touch screen

Handle: 2(HR, SR)

Filter: no

Size: 15*40mm²

Energy: 1-50J/cm²

Interval time: 1-4s

Pulse repetition rate: 1-6 P/s

Plus time: 1-9ms

Fast model plus: 1-8ms

Frequency: 1-10Hz

Power: 3000W

Packing size: 77*65*113cm

Power supply: 220V, 50Hz or 110V, 60Hz

Application:

Hair remove, skin rejuvenation, Pigmentation therapy,

Vascular therapy, Acne therapy



HR handle

SR handle



DPL black gold superphoton(FQ058)

Feature:

Package dimensions:116*64*56cm

Purifying air volume:1120-600M3/h

Screen size:10.4-inch screen

Technology:In-motion sliding technology
energy size:10-40J/cm

Rated voltage:AC110V-220V/50-60Hz

Application:

freckle whitening,remove red blood,skin
rejuvenation,lighten stains,Shrink
pores,skin firming effect



10 inch HD display



DPL handle



Solve multiple skin problems



Picosecond laser machine (FQA32-6)

Feature:

Display: 10.4 inch touch screen

Laser wave length:

1064nm/755nm/532nm


Spot diameter: 1-10mm adjustable

Frequency: 1-10Hz

Energy: 0.4-2.0j/cm²

Application:

Tattoo and eyebrow removal, freckle, Chloasma, age pigment, chorioplague



Picosecond laser machine (FQA32-1)

Feature:

Display: 10.4 inch touch screen

Laser wave length:

1064nm/755nm/532nm

The pulse energy: 1-2000mj

Frequency: 1-10Hz

Cooling system: Water and air cooling


Warranty: 12 month

G.W: 80kg

Package size: 130*68*120cm

Application:

Tattoo and eyebrow removal, freckle, Chloasma, age pigment, chorioplague



Portable Picosecond laser machine(FA31-4)

Feature:

Display: 10 inch touch screen

Laser wave length:

1064nm/785nm/532nm

energy:0.45-1.0J/cm

Frequency:1-10Hz

Cooling system:Water and air cooling

Instrument size:167*53*78

Warranty:12 monthes

Application:


Tattoo and eyebrow removal, freckle, Chloasma, age pigment, chorioplague




Operating arm



1064*532cm



1064cm



Operation interface 1



Operation interface 2



Operation interface 3



Operation interface 4

Non-invasive eyebrow wash (FQ046)

Feature:

Display: LCD 10.4 inches touching screen

Voltage: AC 110V-220V/50-60Hz

Energy: 0.2-1.0mj

Frequency: 01-10Hz

Laser tip:

532nm; 1064nm; 1320nm

Application:

Wash eyebrows, Tattoo Removal age pigment birthmark



3 non-invasive eyebrow washer head



3 common eyebrow washer head



Portable ND yag Laser tattoo removal machine (FQ015-1)

Feature:

Display: LCD 4 inches touching screen

Laser wave length:

Laser tip: 3

1064nm/532nm/Black doll

The pulse energy: 500-2000mj

Pulse width: 1-10P/S

Gross weight: 8Kg

Size: 60*40*48cm

Warranty: 12 monthes

OEM: YES

Application:

Tattoo Removal age pigment birthmark



Laser eyebrow machine

can help you solve any problems



Effect comparison chart



80% off courses of 4
or more treatments



808 semiconductor hair removal|(FQ015-3)

Feature:

Screen size: 4.5 inches

Basic parameter model: FQ015-3

voltage ac: 110v-220v50-60hz

Wavelength: 532; 755; 1064


Energy.: 1-2500mf

frequency:1-10hz Machine weight: 14kg

cooling mode: fan cooling + water cooling

Application:

Remove all color tattoos, Skin regeneration,
Remove melani



Portable Picosecond laser machine (FQA32-4)

Feature:

Laser Type: Emerald Laser

Display: 10.4 inch touch screen

Laser wave length:

1064nm/755nm/532nm/1320nm

The pulse energy: 1-2000mj

Frequency: 1-10Hz

Cooling system: Water and air cooling



G.W: 30KGS

Packing size: 60*40*40cm

Warranty: 12 monthes

Application:

Tattoo and eyebrow removal, freckle, Chloasma, age pigment, chorioplague



Portable Picosecond laser machine (FQA32-5)

Feature:

Laser Type: Emerald Laser

Display: 10.4 inch touch screen

Laser wave length:

1064nm/755nm/532nm/1320nm

The pulse energy: 1-2000mj

Frequency: 1-10Hz

Cooling system: Water and air cooling


G.W: 30KGS

Packing size: 60*40*40cm

Warranty: 12 monthes

Application:

Tattoo and eyebrow removal, freckle, Chloasma, age pigment, chorioplague



Fractional co2 laser machine(FQ045)

Feature:


Product name: Co2 laser
Screen size: 10 inch HD display
Air outlet: Water cooling plus air cooling
language: Chinese, English
Voltage: AC110V-220V 50-60Hz
Instrument size: 137cm*37cm*38cm,
Handle: ABS POWER: 1-60W
model: Lattice Mode, PULSE MODE, PRIVATE MODE, VULVAE MODE

Application:

Vaginal tightening, Repair acne pits,
Whitening and rejuvenation, Wrinkle
Removal, Health and safety




Vaginal probe 1



Vaginal probe 2



Vaginal probe 3



Normal Co2 laser



THE FUNCTION OF THE HEAD

Co2 Fractional

Skin rejuvenation Acne
treatment Freckle spot
removal Stretch marks
removal



Cutting head

Mole Removal Stretch
marks treatment Medical
Cutting use



Vaginal Tighten

Skin rejuvenation Acne
treatment Freckle spot
removal Stretch marks
removal



Fractional co2 laser machine "Function use"



Screen 1



Screen 2



Screen 3



Screen 4

Vertical Vaginal tightening co2 laser machine (FQC02)

Feature:

Display: 10.4 inch touching screen

Laser wave length: 10600nm

Output power: 60w

Duration: 0.1~10ms

Distance: 0.1~2.6ms

Overlap: 1th~20th

Scan: Order, disorder, midsplit


Warranty: 12 month

Gross weight: 65Kg


Packing size: 87*57*114cm

Application:


Scar removal, Remove pigment lesions, birth mark removal, skin rejuvenation, wrinkle removal, vaginal tighten, acne scars treatment




1. Laser irradiation, the depth can be adjusted




2. Laser can irradiation rotate 360-degrees



3. It is tightened and restored laser irradiated sites



4. Vaginal wall contraction and firming



5 in 1 Cavitation rf machine (FQ082)

Feature:

Display: LCD 7 inches touching screen

Cavitation Frequency: 40Khz

Pressure for vacuum: 0.1 Mpa

Power: 220V/110V 50Hz/60Hz

RF Frequency: 1-3MHz

RF Energy: 1-50J

Size: 50x50x40cm

Warranty: 12 monthes

OEM: YES

G.W: 14kg

Application:

Body slimming, Face lifting, Weight loss



Portable Vacuum Cavitation RF machine(FQ083)

Feature:

Display: 7 "TFT color touch screen
 Cavitation frequency:40KHz
 RF frequency:5MHz, multipolar rf
 Vacuum handle:1
 Voltage : 220V/50Hz or 110V/50Hz
 G.W:30kg
 Packing size:56*45*40cm

Application:

Body slimming,Weight loss,tighten skin



New generation cavitation rf slimming machine(FQ043)

Feature:

Screen:10.4 inch touch Screen

Voltage:220V/50Hz or 110V/50Hz

Adsorption head:

Negative:01-10/s -----Release:01-10/s

Intensity:01-08/j -----Time:01-60/min

Lose weight head:

Intensity:01-15/j -----Time:01-60/min

Massage head:

Intensity:01-15/j -----Time:01-60/min

mode:A.B.C.D

Application:

Body slimming,Weight loss,tighten skin



HI-EMT EMsculpture slimming machine (FQEM01)

Feature:

Power:2000W

Screen:10.4 inch touch Screen

Power consumption:2.3KVA

Output intensity:0- 7 Tesla

Shape of stimulation pulse:HI-EMT

Style: Portable/Vertical

Handles:2/4 handles

Voltage : 220V/50Hz or 110V/50Hz

G.W:74kg

Packing size:570x510x1350 (mm)

Application:

Peach hip·Vest line·Beautiful legs and arm lines
19% Reduction in Fat ,16% Increase in Muscle
Mass,Improve confidence in appearance



Gluteus improvement



Abdominal improved



one treatment



FOUR-WAY FOUR-CONTROL+ BIG RF

RF on the whole face of paddles
30 minutes = 30,000 muscle contractions
Comfortable treatment: just lie down and no feelings
180° handle design: fits the curve of arms and thighs better

Cryo cavitation vacuum rf lipo laser slimming beauty machine (FQ039)

Feature:

Power:2000W

Screen:10.4 inch touch Screen

Handpiece:40K,cryo,RF,Lipo laser

Cryo handle size:20*7.5*9(big),16*7.5*7(small)

Temperature:-10~10

Cooling system:Air cooling+Water cooling

Language:Chinese,English,Spanish

Voltage : 220V/50Hz or 110V/50Hz

G.W:60kg

Packing size:132*88*65cm

Application:

fat dissolving,anti cellulite,fat freezing,body shaping, Tightening skin,strengthen skin elasticity,skin lift



R.E.T body slimming machine (FQ079-2)

Feature:

Voltage : 220V/50Hz or 110V/50Hz

LCD screen:8 inch Color touch screen

Package size:115*56*59cm

Intensity:01-15/j

Time:01-60/min

Mode:A.B.C.D

Application:

fat dissolving,anti cellulite,fat freezing,body shaping, Tightening skin,strengthen skin elasticity,skin lift



4 Cryo handle bodyslimming beauty machine (FQ040)

Feature:

Power:2000W

Screen:10.4 inch touch Screen

Handpiece:4 Cryo handle

Size:20*7.5*9(#200),15*6*6(#150)*2pcs,
10*5*6(#100)

Temperature:-10~45 ° C

Cooling system:Air cooling+Water cooling

Language:Chinese,English

Voltage : 220V/50Hz or 110V/50Hz


Packing size:140*70*60cm G.W:60kg

Application:

fat dissolving,anti cellulite,fat freezing,body
shaping, Tightening skin,strengthen skin
elasticity,skin lift



👍 4 handle work at same time or sepearte



360 Cryolipolysis Machine with lipo laser and Cavitation and RF(FQ074)

Feature:

Display: 10 inch color touch screen
 Cryo handles: 5 pcs(Different size)
 Cavitation frequency:40KHz
 RF frequency:1MHz, multipolar rf
 Lipo laser:4 Big, 2 small
 Cryolipolysis temperature:-10 ° C~15° C
 Vacuum pressure: 0-900Kpa
 Lipo laser:650nm
 Voltage : 220V/50Hz or 110V/50Hz
 NW:40KG GW:55kg
 Packing size:130*65*61cm

Application:

Body slimming,Weight loss,tighten skin



BEFORE AND AFTER



Portable 360 Cryolipolysis slimming machine

Feature:

Display: 8.4 inch touch screen

Cryo handle: 3 pcs

Cavitation frequency: 40KHz

Body RF: 5MHz

Lipo laser: 650nm

Pressure output: 0-90Kpa

Freezing temperature: $-10^{\circ}\text{C} \sim 10^{\circ}\text{C}$

Voltage: 220V/50Hz or 110V/50Hz

G.W: 38kg

Packing size: 58*48*79CM

Application:

Body slimming and shap, Weight loss, skin tighten



BEFORE AND AFTER



Cryolipolysis Machine with Cavitation and RF(FQ032)

Feature:

Display: 8 "TFT color touch screen

Cryo handles: 3 handles(big, middle, small)

Cavitation frequency:40KHz

RF frequency:5MHz, multipolar rf

Vacuum: 650mmHg

Pump flow rate: 60L/min

Vacuum pressure: 0-106Kpa

Freezing temperature: -2°C to -10°C

Voltage : 220V/50Hz or 110V/50Hz

G.W:50kg

Packing size: 134*63*62cm

Application:

Body slimming,Weight loss,tighten skin



New arrival 360 degree cryolipolysis slimming machine(FQA18)

Feature:

Display: 8 "TFT color touch screen

Cryo handles: 6 pcs

Cavitation frequency:40KHz

RF frequency:5MHz, multipolar rf

Lipo laser:4Big,2 Small

Vacuum: 650mmHg

Pump flow rate: 60L/min

Vacuum pressure: 0-106Kpa


Freezing temperature:-10~45 ° C

Voltage : 220V/50Hz or 110V/50Hz

G.W:30kg Packing size:61*62*44cm

Application:

Body slimming,Weight loss,tighten skin,double chin



Portable lipo laser Cavitation RF machine(FQ072)

Feature:

Display: color touch screen
 Cavitation frequency: 40KHz
 RF frequency: 5MHz, multipolar rf
 Vacuum: 650mmHg
 Pump flow rate: 60L/min
 Vacuum pressure: 0-106Kpa
 Voltage: 220V/50Hz or 110V/50Hz
 Lipo laser: 6 big, 2 small
 G.W: 30kg
 Packing size: 52cm × 50cm × 67cm

Application:

Body slimming, Weight loss, tighten skin



Portable lipo laser Cavitation RF machine(FQ073)

Feature:

Display: color touch screen
 Cavitation frequency: 40KHz
 RF frequency: 5MHz, multipolar rf
 Vacuum: 650mmHg
 Pump flow rate: 60L/min
 Vacuum pressure: 0-106Kpa
 Voltage: 220V/50Hz or 110V/50Hz
 Lipo laser: 6 big, 2 small
 G.W: 30kg
 Packing size: 52cm × 50cm × 67cm

Application:

Body slimming, Weight loss, tighten skin





Portable lipo laser Cavitation RF machine(FQ061-2)

Feature:

Display: color touch screen
 Cavitation frequency: 40KHz
 RF frequency: 5MHz, multipolar rf
 Vacuum: 650mmHg
 Pump flow rate: 60L/min
 Vacuum pressure: 0-106Kpa
 Lipo laser: 6 big, 2 small
 Voltage : 220V/50Hz or 110V/50Hz
 G.W: 18kg
 Packing size: 52cm × 50cm × 67cm

Application:

Body slimming, Weight loss, tighten skin




Portable lipo laser Cavitation RF machine(FQ062)

Feature:

Display: color touch screen
 Cavitation frequency: 40KHz
 RF frequency: 5MHz, multipolar rf
 Vacuum: 650mmHg
 Pump flow rate: 60L/min
 Vacuum pressure: 0-106Kpa
 Lipo laser: 6 big, 2 small
 Voltage: 220V/50Hz or 110V/50Hz
 G.W: 18kg
 Packing size: 43cm × 42cm × 40cm

Application:

Body slimming, Weight loss, tighten skin



RF face lifting machine (FQ090)

Feature:

1. Display: 8 inch
2. RF Handle: 3pcs (Neck, face, eye)
3. Cooling handle: 1pcs
4. Trolley: Optional
5. Aroma outlet
6. Size: 30*42*48CM
7. Weight: 22kg

Application:

Face lifting, wrinkle remove, skin rejuvenation,



BTL body slimming shaping beauty machine (FQA01)

Feature:

Display: 8 inch color touch screen
 Handles: Two, one for face and one for body
 Handles shots: No limited
 Body: 1-120 adjust (Power)
 10-30°C adjust (Temperature)
 Face: 1-50 adjust (Power)
 10-100 adjust (Duty factor)
 Voltage: 220V/50Hz or 110V/50Hz
 G.W: 50kg
 Packing size: 117*65*65cm

Application:

Fat Reduction, Body Shaping, Skin
 Tightening, Skin Lifting



RFO + Emsculpt(FQ052)

Feature:

Screen size: 10 inch HD display

Air outlet: air cooling

Language: Chinese, English, Japanese

Voltage: AC110V-220V 50-60Hz



Instrument size: 75cm*55cm*55cm

Handle: four/two

Material: sheet metal

Application:

Cellulite reduction Safe and painless Body shaping, Build a perfect figure



Arm



Abdomen



Buttocks



Legs



V



SQUARE



ROUND



HEART



Exercising the abdominal muscles to shape the vest line.



Exercising the buttocks muscles to shape the peach buttocks.



Exercising the oblique muscles to adomen to shape mermaid line.

RET RFO + Emsculpt(FQEM3)

Feature:

Screen size: 10 inch HD display

Air outlet: Water cooling plus air cooling

Language: Chinese, English, Japanese

Voltage: AC110V-220V 50-60Hz

Instrument size: 135cm*55cm*55cm

Handle: four

Material: sheet metal

Application:

1. Fat Reduction
2. Muscle Building
3. Tone Your Body



448k RF(FQ055)

Feature:

Screen size: 10 inch HD display
 Air outlet: Large fan cooling system
 Language: Chinese, English, Japanese
 Voltage: AC110V-220V 50-60Hz
 Instrument size: 104cm*70cm*35cm
 Handle: three
 Material: ABS

Application:

Principle: Lymphatic drainage, improve blood circulation, repair cells, nerves, bones, muscles and ligaments, regulate endocrine, improve sleep, brighten skin to



RES DEEP HEAT SOLUBLE
GREASE HEAD



CAP RF CERAMIC PROBE




FASCIA KNIFE



COMFORT PATCH



PRODUCTS ARE WIDELY USED



Therapeutic principle



Treatment to injury



Body management



Facial care

OUR PRODUCT APPEARANCE



10.4-inch large screen



Visual transparency



Unique handle design




Double handle



Unique handle design



360 four wheels



448k RF(FQ056)

Feature:

Screen size:7 inch HD display

Air outlet:Large fan cooling system

Language;Chinese, English, Japanese

Voltage:Ac110v-220v50-60hz

Instrument size:95cm*45cm*38cm

Handle:three,Material:ABS

Application:

Principle: Lymphatic drainage, improve blood circulation, repair cells, nerves, bones, muscles and ligaments, regulate endocrine , improve sleep, brighten skin to




Res handle



Cap handle



Comfort patch



Liposonic body slimming machine (FQA29/FQA29-1)

Feature:

Display: Touch screen
 CARTRIDGES: 1.3cm, 0.8cm
 SCREEN: 10.4 inch TFT touch screen
 CARTRIDGE TIMES: 525
 ULTRASONIC FREQUENCY: 4MHz
 INPUT: AC110V-220V 50-60Hz
 G.W: 40kg
 Packing size: 117*65*65cm

Application:

Body thinner, skin tightening, Fat reduction,
 Body slimming, body shaping



B-333 Microcurrent EMS Electric Muscle Stimulator (FQB333)

Feature:

Gross weight: 7.5KG

Dimension: 47*38*20cm

Voltage: AC 220V

Frequency: 50Hz

Power: 35W

Width of pulse: 100 μ S~500 μ S

Shrink time: 0.1S~5.0S

Loose time: 0.1S~5.0S

Application:

Losing weight, slimming the body, Firm
slack muscles and muscle mass



980nm Spider vein blood vessels Removal (FQ980)

Feature:

Laser wavelength:980nm

Output power:15w/20w/25w/30w(Optional)

Pulsewidth:15ms-100ms

Frequency:1-5HZ

Modl:Pulse mode

Net/Gross weight:5kg/13kg

Screen:8 inch touch screen

Packing size:54*40*36m

Warranty:12 monthes

Application:

Vascular removal: face, arms, legs and the whole body speckle, age spots, sunburn,pigmentation



HIFU Anti-aging Instrument (FQ041)

Feature:

Screen size: 10 inch HD display

Air outlet: Large fan cooling system

Language: Chinese, English

Voltage: AC110V-220V 50-60Hz

Instrument size: 130cm*60cm*50cm

Handle: four

Material: ABS

Application:

Accurately conduct energy, deeply cover the epidermis, dermis and fat layer, and the superficial fascia layer can carry out comprehensive and in-depth treatment and anti-aging



Lifting handle
(1.5mm, 3.0mm, 4.5mm)



Activation handle
(DC-Y/DC-V)



ice head



zoom probe



Super gun frozen age gun | multifunctional beauty instrument(FQ017-03)

Feature:

Screen size: 16 inch HD display

Air outlet: Large fan cooling system

Language: Chinese, English, Japanese

Voltage: AC110V-220V 50-60Hz

Instrument size: 52cm*47cm*55cm

Handle: two

Material: sheet metal

Application:


Accurately conduct energy, deeply cover the epidermis, dermis and fat layer, and the superficial fascia layer can carry out comprehensive and in-depth treatment and anti-aging



Lift Handle
(1.5mm, 3.0mm, 4.5mm)



sonic cannon
(1.5mm, 2.0mm, 3.0mm, 4.5mm)



Mini home use hifu beauty equipment (FQA28-2)

Feature:

Display: 8 inches touching screen

Treatment cartridge

DS10-1.5 DS7-3.0 DS4-4.5

(DS4-8.0 DS4-13.0 Optional)

Material:ABS

Frequency:7MHz 4MHz


Power:110V~240V 50ZH/60ZH

Packing size:42*34*35m

G.W:15kg

Application:

Firming,Face lifting,shaping



3D hifu machine (A28D)

Feature:

Display: 12 inches touching screen
 Treatment cartridge
 3.0/4.5(Standard)
 1.5/6.0/8.0/10.0/13.0/16.0(Optional)
 Row number:1-11
 Pitch:1.0-4.0mm
 Power:1.0-2.0J
 Lift:10000 shots
 Power:110V~240V 50ZH/60ZH
 Packing size:20*40*54cm
 G.W:15kg
Application:
 Firming,Face lifting,shaping



Face and vaginal hifu machine (FQA28-3)

Feature:

Display:15 inches touching screen

Vaginal cartridges:4.5mm,3.0mm

Hifu cartridges:4.5mm,3.0mm,1.5mm

Power:220V/110V 50Hz/60Hz Trolley:Optional

Hifu cartridge:

1.5mm:Energy:0.1-0.4J; Length:5-14mm

3.0mm:Energy:0.1-0.75J;Length:5-25mm

4.5mm:Energy:0.3-1.2J; Length:5-25mm

Vaginal cartridges:Energy:0.2-1.5J;

Pitch:1.0-5.0mm Length:5-25mm



Packing size:63*38*68cm

Warranty:12 month

OEM:YES G.W:28kg

Application:

Face Lift, Skin Rejuven, vaginal tightening



3D hifu machine (FQA28)

Feature:

Display: 15 inches touching screen

Vaginal cartridges: 4.5mm, 3.0mm

Hifu cartridges: 4.5mm, 3.0mm (11 lines)

Power: 220V/110V 50Hz/60Hz Trolley: Optional

Hifu cartridge:

1.5mm: Energy: 0.1-0.4J; Length: 5-14mm

3.0mm: Energy: 0.1-0.75J; Length: 5-25mm

4.5mm: Energy: 0.3-1.2J; Length: 5-25mm

Vaginal cartridges: Energy: 0.2-1.5J;

Pitch: 1.0-5.0mm Length: 5-25mm


Packing size: 57*43*66cm

Warranty: 12 month



OEM: YES G.W: 19kg

Application:

Skin Tightening, Wrinkle Remover, Anti-aging



Choose The Right Cartridge



7D/9D hifu machine

Feature:

Display: 10.4 inch touch screen
 cartridge: 1.5mm 2.0mm 3.0mm 4.5mm 6.0mm
 9.0mm 13.0mm Life: 20000 shots
 Dot distance: 5.0mm-25mm
 Energy: 0.1-2.0J Power: 220V/110V 50Hz/60Hz
 Package Size: 51*51*131 G.W: 35kg

Application:

1. Remove wrinkles around the forehead, eyes, mouth, etc.
2. Lifting and tightening both cheeks and skin.
3. Improving skin elasticity and shaping contour.
4. Improving jaw line, reducing "marionette lines".
5. Tightening the skin tissue on the forehead, and lifting the eyebrows lines.
6. Improving skin complexion, making the skin delicate and bright.



Face and vaginal hifu machine (FQA28-S)

Feature:

Display:10 inches touching screen

Vaginal cartridges:4.5mm,3.0mm

Hifu cartridges:4.5mm,3.0mm,1.5mm(Standard)
8.0mm,13.0mm(Optional)

Hifu cartridge:

1.5mm:Energy:0.1-0.4J; Length:5-14mm

3.0mm:Energy:0.1-0.75J;Length:5-25mm

4.5mm:Energy:0.3-1.2J; Length:5-25mm

Vaginal cartridges:Energy:0.2-1.5J;

Pitch:1.0-5.0mm Length:5-25mm

Packing size40*54*35cm

Power:220V/110V 50Hz/60Hz


Warranty:12 month OEM:YES G.W:15kg

Application:

Face Lift, Skin Rejuven, vaginal tightening



VAGINAL TIGHTENING PRINCIPLE



Radar anti wrinkle facial lifting line carving ultrasound (FQA28-7)

Feature:

Name: radar line carving

Voltage: 220V/60HZ

Technology principle: Ultrasonic Technology

Net weight: 5kg

Cartridges: 1.5/3.0/4.5nm(Standard)
8.0/13.0nm(Optional)

OEM: YES

Warranty: 12 monthes

Application:

Skin Tightening, Face Lifting, Skin Lifting
around Eyes, Face Wrinkle

Removal, Forehead Wrinkle Removal



Radar anti wrinkle facial lifting line carving ultrasound (FQA28-7A)

Feature:

Name: radar line carving

Voltage: 220V/60HZ

Technology principle: Ultrasonic Technology

Net weight: 5kg

Cartridges: 1.5/3.0/4.5nm(Standard)
8.0/13.0nm(Optional)

OEM: YES

Warranty: 12 monthes

Application:

Skin Tightening, Face Lifting, Skin Lifting
around Eyes, Face Wrinkle
Removal, Forehead Wrinkle Removal



Micro-Needle Fractional RF Anti-ageing machine (FQR01)

Feature:

Display:8" touch screen

Repetition:0.2 / 0.5 / 1 / 2sec / Single

Needle:10 ,25,64 electrode needle

Store temperature:-40~50° C

Wavelength: 650nm

RF frequency:3Mhz

G.W:60KG

Packing size:68*38*113cm

Application:

Wrinkle Reduction, Skin Tightening, Skin Rejuvenation ,Pore Reduction, Birth mark remove



Microneedle Fractional RF (MRF)



MRF(81 needles)



MRF(49 needles)



MRF(25 needles)



Hifu Micro-Needle Fractional RF machine (FQR02)

Feature:

Max. power: 20W $\pm$ 20%
 RF Frequency: 2Hz $\pm$ 20%
 Maximum RF Output Power: 50W
 Intensity: 0 ~ 100 Level (10 Step)
 Repetition: 0.2 / 0.5 / 1 / 2sec / Single
 Needle Thickness: 0.3mm
 Needle Figure: 10, 25, 64 electrode needle
 Treatment Area: 1cm*1cm
 Adjustable Depth: 0.5-3.5 (0.1step)
 Treatment Duration: 50ms-950ms
 G.W: 1KG
 Packing size: 48*35*40cm 7cm

Application:

Wrinkle Reduction, Skin Tightening, Skin
 Rejuvenation, Pore Reduction



No needle

64 pins for body

25 pins for face

10 pins for dead
angle part

Face Lifting

Anti Wrinkles

Marks Remove

Rejuvenation

Fading Scars

Acne Removal

EMRF Face

Feature:

EMRF is a non-invasive facial anti-aging care machine, which uses high-intensity facial electromagnetic stimulation (high-intensity facial electromagnetic stimulation) + Mono-polar Radio Frequency (RF)

G.W:15KG

Packing size:40*43*42cm

Application:

1.wrinkles reduction 2.face lifting & contouring 3.20 mins treatments 4.skin tightening 5.Hand free non invasive lifting effect



3 IN 1 EFFECT

23% MORE LIFT
Lifting Effect on average

30% MUSCLE
Muscle tone on average

30% SKIN
Reduction wrinkle on average

- ✓ Non-invasive and non-surgical plastic surgery
- ✓ Needle-free and non-invasive
- ✓ No fillers and toxins
- ✓ New experience and comfortable feeling
- ✓ Reduce wrinkles
- ✓ Lift loose skin wrinkles
- ✓ Thin face shaping
- ✓ Increase collagen and elastin
- ✓ Simple 20-minute care
- ✓ No manual operation required



Meso gun Beauty Equipment

Function:

Anti-Wrinkle
Extra BIO Whitening
Wrinkle Removal
Skin Firming Serum
Skin Regeneration
BIO Skin Purify
Skin Sensitivity Decreasing
Striae of Pregnancy Removal
BIO Freckle Dispelling
Body Slimming
G.W:1KG
Packing size:43*15*38cm

Package:

Main Machine X 1 Mesogun X 1
Tips X 2 needle and syringe X 1



Far infrared pressotherapy massager lymphatic drainage machine

Feature:

Product name:Lymphatic Drainage evicce

Used for Skin care center,Spa,Medical spa,Clinic,Hospital

Working principle infrared and air ressure
Infrared heater

Material for the suit:Waterproof , fireproof,
PVC polyester fiber Color:blue/gold/pink

Aipressure:100mmHg - 120mmHg

Temperature:35-80° C Power:680W

Voltage 110V-1240V ,50-60HZ

Size Of Machine: 1/2:85cm*46cm*15cm
2/2:45*36*23cm

Weight Of Machine:16kg

Application:

1. Lose weight, eliminate toxin
2. Accelerate metabolism
3. Strengthens nutrition absorbing
4. Dispel tire, remove liver's strain and unknown pain



Before



After



Before



After

Co2 bubble deep cleaning RF facial skin care machine(FQ0025)

Feature:

Principle: Collagen regeneration

Suitable skin types: All type of skin

Handle 1: RF radio frequency head

Handle 2: Beautiful white skin

Handle 3: Ultrasound (Optional)


Touch the screen: 10 inches

Input Power: AC110/220 $\pm$ 10%

50HZ/60HZ GW: 4.5KG

Application:

Whitening and moisturizing and
moisturizing the skin to activate
collagen



Space Oxygen Beauty Facial Spa Salon Machine(FQ078-5)

Feature:

Touch the screen:10 inches

Language:2 Solution:4

Function:Bio,Scrubber,Cooling,Ultrasound,
Camera,Hydra,Spray gun,LED Oxygen
mask,

Packing size:65*90*58cm

Input Power:AC110/220 $\pm$ 10% 50HZ/60HZ

Application:

skin oxygen filling and water injection,
nutrition supplement; hydrating and
whitening, wrinkle removal, skin smoothing,
elasticity and youth restorage

THE PURE OXYGEN FACIAL

An anti-aging facial that delivers
instant results giving your skin a
smooth, moisturized surface

- ✓ Softened expression lines
- ✓ Improved skin elasticity
- ✓ Deep infusion of oxygenated
vitamins and antioxidants




hydra facial(FQ057)

Feature:

Product name: hydra facial
 Screen size: 10 inch HD display
 Air outlet: Water cooling plus air cooling
 Spot size: Chinese, English
 Voltage: AC110V-220V 50-60Hz
 size: 120cm*55cm*43cm, Handle: twelve
 Material: sheet metal, Mode: A B
 Energy: 1-10J, Intensity: 1-10J/m
 Luminous intensity: 1-10J
 work time 1-60/min
 Oxygen injection bottle: S1 S2 S3 S4

Application:

Removal of wrinkles and blackheads, deep skin cleansing and replenishing moisture



Skin wrinkles



sagging skin



Large pores



Skin lack of water



Dead skin horny



Eye bags



Vertical hydral facial skin care machine(FQ078-2)

Feature:

Language:2

Work time:1-60minute

Intensity:1-10

Interval:1-10

Mode:3

Solution:4

Healing head:5

Packing size:123* 54*52cm

Touch the screen:8 inches

Input Power:AC110/220 $\pm$ 10% 50HZ/60HZ

Application:

Replenishment, dilute wrinkles, whitening rejuvenation, remove blackheads, freckle



Acne contrast



Blackhead



Shrink pores



Skin whitening



Wrinkle



Tighten lift



New hydal facial skin care machine(FQ078-4)

Feature:

Work time:1-60minute

Intensity:1-10

Solution:4

Healing head:8pieces

Weight:36kg

Packing size:127*56*59cm

Touch the screen:8 inches

Input Power:AC110/220 $\pm$ 10% 50HZ/60HZ

Application:

Acne treatment,increase skin elasticity,remove black heads,improve skin tone,skin mosturizing



11 In 1 Hydra Dermabrasion (FQ741)

Feature:

Product name: hydra facial

Screen size: 10 inch HD display

Air outlet: Water cooling plus air cooling

language: Chinese, English

Voltage: AC110V-220V 50-60Hz

Instrument size: 120cm*55cm*43cm,

Handle: eight Material: sheet metal,

Mode: A B C D

Intensity: 1-10/jm

work time 1-60/min

Oxygen injection bottle: S1 S2 S3 S4

Application:

Removal of wrinkles and blackheads, deep skin cleansing and replenishing moisture



1 11IN1 functions

Multifunctional skincare system water oxygen meter can satisfy the customer's different requirements.



Big Touch Screen

User-friendly Interface, easy to operation.



8 In 1 Hydra Dermabrasion (FQ039)

Feature:

Product name: hydra facial
 Screen size: 15.6 inch HD display
 Air outlet: Water cooling plus air cooling
 language: Chinese, English
 Voltage: AC110V-220V 50-60Hz
 Instrument size: 120cm*55cm*43cm,
 Handle: eight Material: sheet metal,
 Mode: A B C D
 Intensity: 1-10/jm
 work time 1-60/min
 Oxygen injection bottle: S1 S2 S3 S4

Application:

Removal of wrinkles and blackheads, deep skin
 cleansing and replenishing moistur



Suction pen



Oxygen spray gun



Soft light laser pen



Cold and hot head



Magic bubble pen



Ultrasonic head



Cleaning hand tools



Collagen spheres



Double whitening, save dull yellow skin




Improve facial skin moisturization



Effectively cleans fur excess for dead skin



Effectively remove blackheads



COLLAGEN SPHERES

Shrink pores, determine the ability to repair skinny muscles, and eliminate redness, swelling and sensitive heat effects. Make the skin cell concentration reach 455 degrees, when the time reaches 4550 meters or more, reduce the shrinkage rate of skin particles and achieve the purpose of lightening the skin.



Soft light laser pen

Chondro head therapy to dredge to improve sleep Dredge lymph and accelerate head circulation and metabolism



14 in 1 Hydra Dermabrasion (FQ014)

Feature:

Product name: hydra facial
 Screen size: 15.6 inch HD display
 Air outlet: Water cooling plus air cooling
 language: Chinese, English
 Voltage: AC110V-220V 50-60Hz
 Instrument size: 120cm*55cm*43cm,
 Handle: 14
 Mode: A B C D
 Intensity: 1-10/jm
 work time 1-60/min
 Oxygen injection bottle: S1 S2 S3 S4

Application:

Removal of wrinkles and blackheads, deep skin
 cleansing and replenishing moistur

14 in 1 hydra facial machine



MACHINE APPLICATION



7 In 1 Hydra Dermabrasion(FQ059)

Feature:

Screen size:HD dot touch screen

Air outlet::Large fan cooling system

Voltage:Ac110v-220v50-60hz

Instrument size:45cm*40cm*40cm

Handle:seven, Material:ABS

suction:S1 S2 S3 work time: 1-60 minutes


Ultrasonic handle: 1-10/jm

Lifting handle: 1-10/jm Golden eye:1-10/jm

blade:1-10/jm, Intensity:1-5/jm

Application:

Anti-wrinkle and anti-aging, deep moisturizing, exfoliating, lifting and firming, smoothing fine lines



PRODUCT STRUCTURE



6 In 1 Hydra Dermabrasion(FQ035)

Feature:

Screen size: 8 inch HD display

Air outlet: Large fan cooling system

Language: Chinese, English, Japanese

Voltage: AC110V-220V 50-60Hz

Instrument size: 34cm*44cm*50cm


Handle: six, Material: sheet metal

strength: 1-10J, work time 1-60/min


model: A B C D

Application:

Anti-wrinkle and anti-aging, deep moisturizing, exfoliating, lifting and firming



Screen size: 4.5-inch high-definition display



New hydal facial skin care machine(FQ078-3)

Feature:

Work time:1-60minute

Intensity:1-10

Solution:4

Healing head:6/7

Function:RF lifting,Spray gun,Hydra,
Ultrasound scrubber,Ultrasound,Cooling


Packing size:49*33*45cm

Touch the screen:7 inches

Input Power:AC110/220 $\pm$ 10% 50HZ/60HZ

Application:

Replenishment, dilute wrinkles, whitening
rejuvenation, remove blackheads, freckle



7 in 1 Hydra Dermabrasion Ultrasound RF Bio-lifting Cold Hammer Facial Machine(FQ077-1)

Function:

High-Frequency Skin Scrubber: Through the action of ultrasonic and high-frequency vibration

Ultrasonic Head: The macromolecule nutrients are converted into small molecules

Cold Hammer: Reduces pores, firms and shapes skin.

Ion Spray Pen: replenishes skin with moisture and nutrients.

Needling RF: Promotes the production of collagen and collagen fibers, plays a role in rejuvenating and reducing pores, and improving fine wrinkles.

Hydro Dermabrasion: cleans the skin surface and excess pores in the pores, aging cutin and dirt.

Camera Detection Pen: detects skin problems.

GW: 11.8KG Package size: 62*52*40cm



Facial dryness



Tingling, itching, redness



It was dull, yellowish, and shiny



Roughness of the skin

Intelligent Ice Blue Machine

microdermabrasion & skin scanner analyzer



7/8 in 1 Hydra Dermabrasion Ultrasound RF Bio-lifting Cold Hammer Facial Machine(FQ077-8)

Function:

- 1) Ultrasound head: promote microcirculation, skin rejuvenating
 - 2) Multipolar RF: compact shaping, face lifting
 - 3) Bio microcurrent: let jaw line tightened, Bio raise the overall tightening of facial skin
 - 4) Hydrodermabrasion: facial cleansing, peeling, pore cleansing
 - 5) cold hammer : cold down skin after treatment
 - 6) oxygen spray gun : skin rejuvenation
 - 7) LED face mask 8) Meso gun
- Voltage: 220v or 110v/50Hz-60Hz
Package size: 49*39*39cm GW: 8KG



治疗前

治疗2次后



New 8/10/11 in 1 hydra facial machine

Function:

1. Eye Ultrasound Probe
2. Facial ultrasound
3. Cold and Hot Hammer
4. RF Lifting Probe
5. Young Girl Line
6. Shovel Knife Skin scrubber
7. Hydra Water Facial Bubble probe
8. Oxygen Spray Gun
9. Massage Magic Bolls
10. Oxygen Bubble Pen
11. LED Mask

Voltage: 220v or 110v/50Hz-60Hz

GW: 9.6KG

Package size: 45*43*40.5cm

We have 3 different
Accessories combination

Please contact us!!!


8 in 1 hydrofacial



10 in 1 hydrofacial



11 in 1 hydrofacial



Before use



After use



Lifting
and firming



Clean pores



Hyperbaric
oxygen injection



Antioxidant
repair



Whitening



Moisturizing



Fade
fine lines



Deep
absorption



Acne



Dry and short of water



Dark deposit



Blackhead acne



Flabby skin



Coarse pores



Facial wrinkles



Makeup residue

Body Composition Analyzer machine (FQ108)

Feature:

- 1.Operation System: bioelectrical impedance analysis
- 2.Feature: body fat analyzer, body composition analyzer
- 3.Testing weight range: 10-200 Kg.
- 4.Test items: 25 items
- 5.Weight: 40kg
- 6.Size: 50.4x81.8x116.5cm
- 7.Controller:320*240 LCD touch screen
- 8 . Electrode method: 8 point contact electrode
9. Operation language: English/multi-language optional

Application:

Analyze body moisture, Analyze fat, Analyze body composition, Measure weight .



BODY ANALYZER
Body Composition Analysis



Radar anti wrinkle facial lifting line carving ultrasound (FQA28-5)

Feature:

Name: radar line carving
 Voltage: 220V/60HZ
 Technology principle: Ultrasonic Technology
 Net weight: 12kg
 Product size: 37*40*20 (CM)
 Package size: 60*57*32 (CM)
 OEM: YES
 Warranty: 12 monthes

Application:

Desalination of wrinkles, pulling tight skin, whitening and skin rejuvenation, freckle



VMax Wonder Lifting



Dr.pen derma roller beauty machine (FQDP)

Feature:

Needle Numbers:12pins

Needle length:From 0.25mm to 2.0mm

Adapter:5V –3000MA

Property:With vibrating ,interchangeable head

Voltage : 220V/50Hz or 110V/50Hz

G.W:3kg

Packing size:25*25*13cm

Application:

Anti-Puffiness.Acne Scar Removal&treatment.

Anti Wrinkle Treatment, Anti Aging.Surgical Scar Removal.Hyper pigmentation treatment.Shrink pores.



Reference Chart

REGION	THIN SKIN	THICK SKIN
Forehead	0.25-0.5mm	0.5-0.75mm
Between Eyebrow	0.25-0.5mm	0.5-1.0mm
Nose	0.25mm	0.5mm
Eye Area	0.25mm	0.25-0.5mm
Cheek Bone	0.5mm	0.5-1.0mm
Cheek and Chin	0.5-1.0mm	1.0-2.0mm
Upper Lip	0.25mm	0.25-0.75mm
Facial Scarring	1.0-1.25mm	1.5-2.0mm
Stretch Mark	1.5-2.0mm	1.5-2.0mm
Surgical & Bum Scars	2.0-2.5mm	



Omega light beauty machine (FQW01)

Feature:

Power:36W

Lamp:3EA

Bulb:3528

Bulbs:1014pcs

Cooling system:Air fans

Color:4 color,red,blue,green,yellow

Voltage : 220V/50Hz or 110V/50Hz


Machine Size:55*34*32CM

G.W:5kg

Packing size:57*35*33cm

Application:

aging, sagging, fine lines, wrinkles, enlarged pores, pigmentation freckle, sunburn, acne inflammation, gloomy



LED face Mask (FQW02)

Feature:

Output voltage:DC12V

Output power:10W (Max)

Security: 1-Type BF

Color:7 colors(red,blue,green,yellow,purple,clear blue,laser white) to be adjustable

Voltage : 220V/50Hz or 110V/50Hz


Machine Size:26.5*21.5*13cm

G.W:2kg

Packing size:33*28*20cm

Application:

- Kill bacteria that cause spots and pimples
- Great alternative to oral medications
- Prevent breakouts,Minimize scarring
- Effective for kinds of skin care treatments.



Omega light beauty machine (FQW03)

Feature:

Output Power:36W

Light Source:LED photo light pulsation

Light Wavelength:

Red light: 640nm, Blue light: 423nm, IR light: 830nm

Led Tubes:290 pcs

Material of shell:PC

Operation Modulation:Red , Blue , IR Color

Time Setting:5-60 minutes

Control Mode:Digital control

Work Voltage:DC 12V 3A

Input Voltage :AC100-240V

Dimensions:473(W)*484 (D)* 343 (H) mm

Weight:G.W: 4.5Kg / N.W:2.29Kg

Application:

Skin Rejuvenation , Wrinkle Removal , Acne treatment ,
Remove Age spots and Dark Pigmentation . Face and
hand spa treatment to be beautiful and younger.



Photosynthesis SPA beauty equipment(FQ033)

Feature:

Screen: touch screen

Oxygen: 1-10 / J (adjustable)

Atomization: Energy: 1-15 / (adjustable)

Cold handle: 1-10 (adjustable)

WD intensity: 1-12 (adjustable)

time: 1-60 (minute)

Application:

 Photosynthetic spa instrument 1 Space pure oxygen mask has a pure oxygen rate of more than 90% and improves the lack of oxygen of the skin. 2 collagen ball uses high-tech alternating cold and heat technology to eliminate fatigue and restore skin elasticity. 3 water electrophoresis probe can introduce various active ingredients into mesoderm, reduce wrinkles, remove spots and whiten











Production room and showroom



Guangzhou Baiyun District Fuqiang Beauty Equipment Factory is a professional beauty equipment manufacturer, we have our own technical team, D&R department, Sale department. Our main product is hifu, co2 laser, diode laser, cryolipolysis, cavitation rf, liposonix, nd yag laser, shr hair removal and so on.



Visitor & Exhibition & Shipping





Certification of Conformity

EU- ELECTROMAGNETIC COMPLIANCE DIRECTIVE-2014/30/EU

Registration No.: ENC1809078GZ42-1

Applicant : Guangzhou Fujiang Beauty Equipment Co., Limited
 Applicant Address : No. 103, Huasheng South Road, Donghua Industrial Park, Renhe Town, Baiyun District, Guangzhou, Guangdong, China (Mainland)
 Product Designation : Beauty Machine
 Model Number : FQ808, FQD03, FQD05, FQAD5, FQAD5-1, FQAD5-2, FQZ03, FQA31-4, FQA20, FQA31, FQA32, FQ039, FQ077, FQ078, FQ082, FQ032, FQ075, FQA01, FQA19, FQ015, FQ016, FQA28, FQA28-3, FQG01, FQA28-6, FQDP, FQA32-3, FQA32-4, FQA17, FQAD5-3, FQM01
 Brand Name : FQ
 Manufacturer : Guangzhou Fujiang Beauty Equipment Co., Limited
 Manufacturer Address : No. 103, Huasheng South Road, Donghua Industrial Park, Renhe Town, Baiyun District, Guangzhou, Guangdong, China (Mainland)

The submitted products have been tested by us with the listed standard and found in compliance with the following European Standards:

Directive	Applied Standards	Test Report No.
EMC Directive 2014/30/EU	EN 55014-1:2017	ENC1809078GZ42E1
	EN 61000-3-2:2014	
	EN 61000-3-3:2013	
	EN 55014-2:2015	

This certificate is based on an evaluation of a sample of the above mentioned product. Technical report and documentation are at the licensee applicant's disposal. This certificate does not imply assessment of the series-production of the product. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation.



Ray Zhou / General Manager
 Date of Issue: Sep. 17, 2018

East Notice Certification Service Co., Ltd.

1/F, Haohui Commercial Building, Zhuji Street, Dongpu Town, Tianhe District, Guangzhou City
 Tel: +86-020-8256 8534 Fax: +86-020-8256 8534 E-mail: enc@enc-lab.com Web: www.enc-lab.com



Certification of Conformity

EU-LOW VOLTAGE DIRECTIVE-2014/35/EU

Registration No.: ENC1809078GZ42-2

Applicant : Guangzhou Fujiang Beauty Equipment Co., Limited
 Applicant Address : No. 103, Huasheng South Road, Donghua Industrial Park, Renhe Town, Baiyun District, Guangzhou, Guangdong, China (Mainland)
 Product Designation : Beauty Machine
 Model Number : FQ808, FQD03, FQD05, FQAD5, FQAD5-1, FQAD5-2, FQZ03, FQA31-4, FQA20, FQA31, FQA32, FQ039, FQ077, FQ078, FQ082, FQ032, FQ075, FQA01, FQA19, FQ015, FQ016, FQA28, FQA28-3, FQG01, FQA28-6, FQDP, FQA32-3, FQA32-4, FQA17, FQAD5-3, FQM01
 Brand Name : FQ
 Manufacturer : Guangzhou Fujiang Beauty Equipment Co., Limited
 Manufacturer Address : No. 103, Huasheng South Road, Donghua Industrial Park, Renhe Town, Baiyun District, Guangzhou, Guangdong, China (Mainland)

The submitted products have been tested by us with the listed standard and found in compliance with the following European Standards:

Directive	Applied Standards	Test Report No.
LVD Directive 2014/35/EU	EN 60335-1:2012+A11:2014+AC:2014+A12:2017	ENC1809078GZ42L1
	EN 60335-2-23:2003+A1:2008+A2:2015	
	EN 62233:2008	

This certificate is based on an evaluation of a sample of the above mentioned product. Technical report and documentation are at the licensee applicant's disposal. This certificate does not imply assessment of the series-production of the product. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation.



Ray Zhou / General Manager
 Date of Issue: Sep. 17, 2018

East Notice Certification Service Co., Ltd.

1/F, Haohui Commercial Building, Zhuji Street, Dongpu Town, Tianhe District, Guangzhou City
 Tel: +86-020-8256 8534 Fax: +86-020-8256 8534 E-mail: enc@enc-lab.com Web: www.enc-lab.com

We are very old in this market,our quality is our culture and we believe we sell good service not only sell products.

Hope we can have chance to serve your company.

OEM/ODM Welcome
Exclusive design for
your company only