

Precision, Passion, Progress



EAGER Surgical

Reliable Surgical Instruments

Empowering dental professionals with precision-crafted instruments designed for accuracy, reliability, and performance

**DENTAL INSTRUMENTS CATALOGUE
2025 EDITION**

EAGER Surgical Instruments
Dallas, Texas, USA
Website: www.eagersurgical.com
Email: info@eagersurgical.com
+1 469 290 0254





 Precision, Passion, Progress

www.eagersurgical.com

About Us

At EAGER Surgical, headquartered in Dallas, Texas, we provide precision-crafted surgical, dental, and veterinary instruments trusted by professionals worldwide. Supported by our advanced manufacturing facility in Pakistan and guided by a commitment to uncompromising quality, we blend skilled craftsmanship with customer-focused service to equip clinics and hospitals with tools that enhance safety, reliability, and care.

info@eagersurgical.com

Dallas, Texas, USA

| DIAGNOSTIC

Diagnostic kits.....
Explorers, Single-end.....
Explorers.....
Probes, single-end.....
Handles for mouth mirrors.....
Mouth mirrors.....
Dressing pliers.....
Articulating papers forceps.....

| ENDODONTICS

Pliers for holding point with self-locking.....
Pliers for removing broken broaches.....
Spreaders, Pluggers, Condensers.....

| RESTORATIVE

Excavators
Amalgam plugger/condensers.....
Plastic filling instruments.....
Composite instruments (TIN)
Amalgam burnishes.....
Carver.....
Titanium nitride (TIN) carvers.....
Cement spatulas.....
Heidemann spatulas.....
Placements instruments.....
Crown remover.....

| PERIODONTAL

Sickle Scalers
Universal curettes.....
Gracey curettes.....
Gracy mini curettes.....
Sharpening accessories.
Sharpening stone.....
Sharpening stone.....



| SURGERY

Periodontal Knives.....
Surgical chisels.....
Scalpel handles.....
scalpel blade remover.....
Surgical mallets.....
Surgical curettes.....
Periosteals.....
Sinus-lift instruments.....
Curettes for sinus- lift.....
Spatole modellbaili-
Retractor.....
Retractors.....
Osteotomes.....
Bone scrapers.....
Compactors.....
Bone Well.....
Bone mill.....
Surgical aspirators.....
Gum scissors.
Surgical Scissors.....
Suture Scissors.....
Needle Holders.....
Needle holders with tungsten carbide inserts.....
Hemostats.....
Sterilizing Forceps.....
Tissue Pliers.....
Peritome.....
Extracting forceps.....
Bone Rongeurs.....
Root elevators.....
Root tip pick / teaser.....
Root fragment removing pliers.....
Luxators.....

| ORTHODONTICS


Pliers
Cutters / Distal and cutters with tungsten carbide inserts.....
D.B. Brackets holders.....
Orthodontic ligature pliers.....
Utility instruments.....

| LABORATORY

Wax carves.....
P.K. Thomas waxing instruments.....
Ceramic sculpturing instruments.....
Wax knives.....
Plaster knives.....
Mixing Spatulas.....
Crown Scissors.....
Gold Foil Scissors.....
Scalpel handles.....
Alcohol lamp.....
Pliers / Laboratory Pliers.....
Mallets-calipers.....
Insulated soldering tweezers.....
Crucible tongs.....



Atraumatic Root Extraction System



*Precision for Minimally
Invasive Extractions*



Handle Options



Scaler Single End 6mm Solid Handle



Scaler Single End 8mm Hollow Handle



Scaler Single End 10mm Hollow Handle



Scaler Double End 6mm Solid Handle



Scaler Double End 8mm Hollow Handle



Scaler Double End 10mm Hollow Handle



Dental Tweezer Plain Pattern



Dental Tweezer Learning Pattern






Handle Options



Elevator Pattern #G01



Elevator Pattern #G02



Extraction Forceps Pattern # P01



Extraction Forceps Pattern # S02



Extraction Forceps Pattern # SF03



Diagnostic



Diagnostic Kit



EG 001
Explorer single-end #2

EG 002
Handle for mouth mirrors Ø 6 mm. round #2
Ø 22 mm. Plain #4P
Ø 24 mm. Plain #5P
Mouth mirror

EG 003
Dressing pliers #2 College



EG 004
Explorer single-end #23

EG 005
Handle for mouth mirrors Ø 6 mm. round #2
Ø 22 mm. - Plain #4P
Ø 24 mm. - Plain #5P
Mouth mirror

EG 006
Dressing pliers #2 College



EG 007
Explorer single-end #9

EG 008
Handle for mouth mirrors Ø 6 mm. round #2
Ø 22 mm. - Plain #4P
Ø 24 mm. - Plain #5P
Mouth mirror

EG 009
Dressing pliers #2 College



EG 010
Explorer single-end #9

EG 011
Handle for mouth mirrors Ø 6 mm. round #2
Mouth mirror plain Ø 22 mm. #4P

EG 012
Dressing pliers #17 Flagg



EG 013
Explorer single-end #11



EG 014
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 015
Dressing pliers #2 College



EG 016
Explorer single-end #6A



EG 017
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 018
Dressing pliers #2 College



EG 019
Explorer single-end #23



EG 020
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 021
Dressing pliers #18 Meriam



EG 022
Explorer single-end #9



EG 023
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 024
Dressing pliers #18 Meriam



EG 025
Explorer double-end #23/17



EG 026
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 027
Dressing pliers #18 Meriam



EG 028
Explorer single-end #8



EG 029
Handle for mouth mirrors Ø 6 mm, round #2
Mouth mirror plain Ø 22 mm, #4P



EG 030
Dressing pliers #18 Meriam

| Explorers - single - endEG 031
2EG 032
23EG 033
8EG 034
9EG 035
12EG 036
6EG 037
17EG 038
11EG 039
17AEG 040
6A**sonde - punta doppia | explorers - double-end**EG 041
23/17EG 042
DG16EG 043
ODU 11/12

| Probes - Single - End



EG 044
W - Williams
1-2-3-5-7-8-9-10



EG 045
12S - WHO
Ball end
3-6-9-12



EG 046
11.5 - WHO
Punta sferica - Ball end
3,5-5,5-8,5-11,5



EG 047
11
3-6-8-11



EG 047
15
1-2-3-4-5-6-7-8-9-10-11-12-13-14-15



EG 049
12
3-6-9-12



EG 050
C
Ball end
3-6-9-12



| Probes - double - end



EG 051
 # GF/W - Goldman-Fox/Williams
 1-2-3-5-7-8-9-10



EG 052
 # 1N - Nabers



EG 053
 # 23/W - Williams
 1-2-3-5-7-8-9-10



EG 054
 # Q2N - Nabers
 3-6-9-12



EG 055
 # 23/15
 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15

| handles for mouth mirrors



EG 056
 # 1
 Octagonal



EG 057
 # 2
 Ø 6 mm.



EG 058
 # 2A
 Ø 7 mm.

| Mouth Mirrors



EG 059
4P
Ø mm. 22 - Plane




EG 060
5P
Ø mm. 24 - Plane




EG 061
4DS
Ø mm. 22 - Plane
Double side Rhodium surface.



EG 062
5DS
Ø mm. 24 - Plane
Double side Rhodium surface.



ALUMINIUM - BACK SURFACE



RHODIUM - FRONT SURFACE

- 1 - stainless steel frame
- 2 - sealed with waterproof seal
- 3 - reflective aluminum-coated back surface
- 4 - reflective rhodium-coated front surface

In mirrors with rhodium-coated front surface, the reflective surface is on the upper part of the mirror for precise and sharp reflected images. It is recommended to avoid contact between the mirror and abrasive surfaces. When cleaning, it also is recommended to avoid immersing the mirror in an ultrasonic tank.

| Dressing Pliers



EG 063
1 - Perry



EG 064
2 - College



EG 065
17 - Flagg



EG 066
18 - Meriam



EG 067
18L
Autobloccante - Self-locking

articulating paper forceps




EG 068
M1 - Miller
Straight



EG 069
M2 - Miller
Curved

Endodontics



| Pliers for holding points with self-locking



EG 070
EP1


| pliers for removing broken broaches



EG 071
1
Straight

EG 072
2
Curved

| spreaders, pluggers, condensers



EG 073
D11 Standard
23 mm. - Ø 0.40 mm.

1S Wakai
22 mm. - Ø 0.30 mm.

Ø 1.00
Ø 0.90 20
Ø 0.80 15
Ø 0.70 10
Ø 0.60 5

EG 074
60
Ø 0.60 mm.


Ø 1.20
Ø 1.10 20
Ø 1.00 15
Ø 0.90 10
Ø 0.80 5

EG 075
80
Ø 0.80 mm.

Ø 1.20
Ø 1.15 20
Ø 1.10 15
Ø 1.05 10
Ø 1.00 5

EG 076
100
Ø 1.00 mm

Restorative



| Excavators



EG 077
0



EG 078
1



EG 079
17W



EG 080
2



EG 081
18



EG 082
3



EG 083
18W



EG 084
4



EG 085
38/39



EG 086
6



EG 087
31L



EG 088
14



EG 089
32L



EG 090
17



EG 091
33L

| Amalgam Pluggers / Condensers



EG 092
1CL



EG 093
2CL



EG 094
3CL



EG 095
0 CL



EG 096
1 CLS
Serrated



EG 097
2 CLS
Serrated



EG 098
3 CLS
Serrated



| Plastic Filling Instruments



ES 099
10



ES 100
11



ES 101
3P



ES 102
12



ES 103
1PS
1.40 2.00
mm.



ES 104
1P



ES 105
2PS
1.90 2.50
mm.



ES 106
2P



ES 107
3P
3.00 3.00
mm.



TIN

TIN (titanium nitride), which has always been very fascinating for its shimmering gold coloring, was the first coating to be successfully used in dentistry due to its high biocompatibility. The TIN coating significantly increases the micro-hardness of the instrument. These features, together with a low coefficient of friction, render it the ideal coating for all the instruments that come in contact with any composites and cements.

| Titanium Nitride (TIN) Composite Instruments



ES 108
Ø6 mm. handle



ES 109
Ø8 mm. hollow handle



ES 110
4



ES 111
4



ES 112
5



ES 113
5



ES 114
6



ES 115
6



ES 116
9



ES 117
9



ES 118
F4 - Felt



ES 119
F4 - Felt



ES 120
G4 - Goldstein



ES 121
G4 - Goldstein



ES 122
3R



ES 123
3R



ES 124
10



ES 125
10

| Amalgam Burnishers




ES 126
 # 1PL
 1.00 1.20




ES 127
 # 2PL
 1.40 1.70




ES 128
 # 3PL
 1.90 2.20



ES 129
 # 1WT - Westcott



ES 130
 # 2WT - Westcott




ES 131
 # 27/29B - Ball
 1.20 2.20

| Carvers



ES 132
1/2HL - Hollenback



ES 133
3SHL - Hollenback



ES 134
2 IPOA - Interproximal Off-Angle



ES 135
3HL - Hollenback



ES 136
8



ES 137
1 IP - Interproximal

| cleoid-discoid



ES 138
CD 1/2



ES 139
CD 89/92

| titanium nitride (tin) carvers



ES 140
1 IP - Interproximal
Hollow handle
Ø 8 mm.



ES 141
2 IPOA - Interproximal Off-Angle
Hollow handle Ø 8 mm.

| Cement Spatulas



ES 142
00



ES 143
0



ES 144
1



ES 145
2



ES 146
3



ES 147
25

| Heidemann spatulas



ES 148
HD1 # HD1R
Flexible *Rigid*



ES 149
HD2 # HD2R
Flexible *Rigid*



ES 150
HD3 # HD3R
Flexible *Rigid*



| Placement Instruments



ES 151
CHP 1



ES 152
CHP 2



ES 153
CHP 3

| crown remover



ES 154
CR
3 tips included



ES 155
1




ES 156
2



ES 157
3



Periodontal



Titanium Blue Gracey-Curettes

- Hollow handle Ø 10 mm.



3...1	1...3
3...1	1...3

ES 158
1/2 GR - Gracey



5...1	1...5
5...1	1...5

ES 159
3/4 GR - Gracey
Lingual/palatal surfaces



8...4	4...8
8...4	4...8

ES 160
9/10 GR - Gracey
*Buccal and lingual/palatal surfaces -
More angled*



5...1	1...5
5...1	1...5

ES 161
5/6 GR - Gracey



8...4	4...8
8...4	4...8

ES 162
11/12 GR - Gracey
Mesial



8...4	4...8
8...4	4...8

ES 163
7/8 GR - Gracey



8...4	4...8
8...4	4...8

ES 164
13/14 GR - Gracey
Distal

Titanium Blue Scalers

Hollow handle Ø 10 mm.



5...1	1...5
5...1	1...5

ES 165
MIM817
H 5/33



5...1	1...5
5...1	1...5

ES 166
H 6/7



8...6	6...8
8...6	6...8

ES 167
204 S



8...1	1...8
8...1	1...8

ES 168
T2/T3 - Taylor



8...4	4...8
8...4	4...8

ES 169
204 SD

| Gracey mini-curettes

Hollow handle Ø 8 mm.



3...1	1...3
3...1	1...3

ES 170
1/2 GRL 3M - Gracey



5...1	1...5
5...1	1...5



ES 171
3/4 GRL 3M - Gracey



8...4	4...8
8...4	4...8



ES 172
11/12 GRL 3M - Gracey
Mesial



5...1	1...5
5...1	1...5



ES 173
5/6 GRL 3M - Gracey



8...4	4...8
8...4	4...8



ES 174
13/14 GRL 3M - Gracey
Distal



8...4	4...8
8...4	4...8



ES 175
7/8 GRL 3M - Gracey
palatali

| Gracey Curettes

- Hollow handle Ø 10 mm.



3...1	1...3
3...1	1...3

ES 176
1/2 GR - Gracey
Buccal surfaces



5...1	1...5
5...1	1...5

ES 177
3/4 GR - Gracey
Lingual/palatal surfaces



8...4	4...8
8...4	4...8

ES 179
11/12 GR - Gracey
Mesial



5...1	1...5
5...1	1...5

ES 178
5/6 GR - Gracey



8...4	4...8
8...4	4...8

ES 181
13/14 GR - Gracey
Distal



8...4	4...8
8...4	4...8

ES 180
7/8 GR - Gracey
palatal surfaces



8...4	4...8
8...4	4...8

ES 183
15/16 GR - Gracey
Mesial




8...4	4...8
8...4	4...8

ES 182
9/10 GR - Gracey
*palatal surfaces -
More angled*




8...4	4...8
8...4	4...8

ES 184
17/18 GR - Gracey
Distal



UPPER



LOWER

- # 1/2 - # 3/4 - # 5/6
- # 7/8 - # 9/10
- # 11/12 - # 15/16
- # 13/14 - # 17/18

Gracey curettes

- Hollow handle Ø 8 mm.



3...1	1...3
3...1	1...3

ES 185
1/2 GR - Gracey
Buccal surfaces



5...1	1...5
5...1	1...5

ES 186
3/4 GR - Gracey
palatal surfaces



5...1	1...5
5...1	1...5

ES 187
5/6 GR - Gracey



8...4	4...8
8...4	4...8

ES 188
7/8 GR - Gracey
palatal surfaces



8...4	4...8
8...4	4...8

ES 189
9/10 GR - Gracey
*palatal surfaces -
More angled*




8...4	4...8
8...4	4...8

ES 190
11/12 GR - Gracey
Mesial



8...4	4...8
8...4	4...8

ES 191
13/14 GR - Gracey
Distal



Universal Curettes

Hollow handle 08 mm.




3...1	1...3
3...1	1...3

ES 192
2R/2L C - Columbia



8...4	4...8
8...4	4...8

ES 193
4R/4L C - Columbia



8...4	4...8
-	-

ES 194
L 3/4 - Langer



8...4	4...8
8...4	4...8

ES 195
13/14 C - Columbia



3...1	1...3
3...1	1...3

ES 196
L 5/6 - Langer



8...6	6...8
8...6	6...8

ES 197
YG7/8 - Younger Good



8...1	1...8
8...1	1...8

ES 198
M 23A




-	-
8...4	4...8

ES 199
L 1/2 - Langer




8...1	1...8
8...1	1...8

ES 200
M 23B




Straight and curved sickle scalers feature:

- 2 cutting edges
- curved facial component
- tipped ends



Universal curettes are characterized by an angle of 90° between the terminal stem and the cutting component.



Gracey curettes are characterized by an angle of 70° between the terminal stem and the cutting component.

Sickle Scalers

Hollow handle 08 mm.



3...1	1...3	ES 201 # U15/33 - Towner-Jaquette
3...1	1...3	



3...1	1...3	ES 202 # U15/30 - Towner-Jaquette
3...1	1...3	



8...6	6...8	ES 203 # 11/12 MC - McCall <i>Mesial</i>
8...6	6...8	



5...1	1...5	ES 204 # H5/33
5...1	1...5	



8...6	6...8	ES 205 # 13/14 MC - McCall <i>Distal</i>
8...6	6...8	



5...1	1...5	ES 206 # H6/H7
5...1	1...5	



8...1	1...8	ES 207 # 13S/14S MC - McCall
8...1	1...8	



5...1	1...5	ES 208 # 30/33 - Jaquette
5...1	1...5	



8...6	6...8	ES 209 # 17S/18S MC - McCall <i>Distal</i>
8...6	6...8	



8...1	1...8	ES 210 # T2/T3 - Taylor
8...1	1...8	



8...1	1...8	ES 211 # M23
8...1	1...8	



8...4	4...8	ES 212 # 204
8...4	4...8	




8...4	4...8	ES 213 # 204SD
8...4	4...8	



8...6	6...8	ES 214 # 204S
8...6	6...8	

Surgery



| **Periodontal Knives**



ES 215
15K/16K - Kirkland



ES 216
1/2 O - Orban



ES 217
8 GF - Goldman-Fox



ES 218
11 GF - Goldman-Fox

| **surgical chisels**



ES 219
36/37 RH - Rhodes



ES 220
1 OB - Ochsenbein



ES 221
2 OB - Ochsenbein

| Scalpel Handles



ES 222
#3



ES 223
#3



ES 224
#4



ES 225
#3L



ES 226



ES 227
14.5cm



ES 228
16.2cm Scalpel ajustable
180 degree

| scalpel blade remover



| **Surgical Mallets**



ES 229
MM20 - Mead
cm. 18,5



ES 230
MM19 - Mead
cm. 18,5

| Surgical Currettes

ES 231
1 HY - Hemingway



ES 232
2 HY - Hemingway



ES 233
3 HY - Hemingway



ES 234
4 HY - Hemingway



ES 235
85 L - Lucas



ES 236
86 L - Lucas



ES 237
87 L - Lucas



ES 238
88 L - Lucas

| Periosteals



ES 239
M 9 - Molt



ES 240
PR 3 - Prichard



ES 241
F 15 - Freer
Curved



ES 242
F 16 - Freer
Strongly curved



ES 243
H 1 - Hopkins



ES 244
H 2 - Hopkins



ES 245
ME 3 - Mead



ES 246
M 8 - Molt



ES 247
GF 14 - Goldman-Fox

| Periosteals



ES 248
24 G



ES 249
SCS 02
Curved discoid

| sinus-lift instruments



ES 250
1



ES 251
2



ES 252
3



ES 253
4



ES 254
5



ES 255
6

| Currettes for sinus-left



| Tongue Depressors and malleable spatula / retractor



ES 262
A145
cm. 14,5




ES 263
A17
cm. 17

ES 264
SM 16
cm. 20



ES 265
SM 18
cm. 20



ES 266
SM 8
cm. 24



| Retractors



ES 267
S16 - Stenberg
cm. 16



ES 268
C14 - Columbia
cm. 14



ES 269
F 13 - Farabeuf
cm. 13



ES 270
F 15 - Farabeuf
cm. 15

| Retractors



ES 271
K1 - Kocher
cm. 21 (mm. 35x11)

ES 272
K2 - Kocher
cm. 21 (mm. 40x11)

ES 273
A1
Atraumatic
cm. 21 (mm. 42x11)

ES 274
cm. 21 (mm. 35 x 19)

ES 275
cm. 21 (mm. 40 x 19)

ES 276
cm. 21 (mm. 55 x 19)

| Compactors



ES 277
1




ES 278
2



ES 279
BB

| bone well



ES 280
Bone Crusher



ES 281
Bone Crusher

| Aspirators

With removable wire stylet for cleaning of the tip.

Ø 3 mm.



Ø 2 mm.



Ø 2,5 mm.



Ø 3,5 mm.



Ø 3 mm.



Ø 4 mm.



Ø 5 mm.



| Gum Scissors



ES 290
17 - Iris Delicate
Straight
cm. 11,5

ES 291
17 TC - Iris
Straight
cm. 11,5

ES 292
17 SC - Iris
Straight
cm. 11,5

ES 293
17 - Iris
Straight
cm. 11,5



ES 294
18 - Iris Delicate
Curved
cm. 11,5

ES 295
18 TC - Iris
Curved
cm. 11,5

ES 296
18 SC - Iris
Curved
cm. 11,5

ES 297
18 - Iris
Curved
cm. 11,5



1 blade serrated, double curved.



ES 298
14 - La Grange
cm. 11,5

ES 299
14 SC - La Grange
1 blade serrated
cm. 11,5

ES 300
14 - La Grange
1 blade serrated
cm. 11,5



see pages 10-13

| Gum Scissors



ES 301
 # 1 SF - Super Fine
Straight
 cm. 11,5

ES 302
 # 2 SF - Super Fine
Curved
 cm. 11,5



ES 303
 # 1 - Kelly
Straight
 cm. 16

ES 304
 # 2 - Kelly
Curved
 cm. 16



ES 305
 # 31 - Noyes
 cm. 11



ES 306
 # 16 RL - Goldman-Fox
Straight
 cm. 13

ES 307
 # 16 RTC - Goldman-Fox
Straight
 cm. 13



ES 308
 # 16 L - Goldman-Fox
Curved
 cm. 13

ES 309
 # 16 TC - Goldman-Fox
Curved
 cm. 13



| Surgical Scissors

**ES 310**# 1MD - Metzenbaum
Delicate
Retta - Straight
cm. 14,5**ES 311**# 2MD - Metzenbaum
Delicate
Curva - Curved
cm. 14,5**ES 312**# 1MFTC
Metzenbaum Fino
cm. 14,5**ES 313**# 1J TC - Joseph
Curved sharp
cm. 14**ES 314**# 2J TC - Joseph
Curved blunt
cm. 14

| suture scissors

**ES 315**# S 11 - Spencer
For sutures
cm. 11

| Needle Holders

ES 316
MH16 - Mayo-Hegar
cm. 16



ES 317
CW15 - Crile-Wood
cm. 15




ES 318
M14 - Mathieu
cm. 14



ES 319
CV14R - Castroviejo
Straight
cm. 14



ES 320
CV14C - Castroviejo
Curved
cm. 14



| Needle Holders with tungsten carbide inserts

ES 321
MH16T - Mayo-Hegar
cm. 16



ES 322
CW15T - Crile-Wood
cm. 15



ES 323
MR15T - Mini Ryder
cm. 15



ES 324
OH14T - Olsen-Hegar
cm. 14



ES 325
M14T - Mathieu
cm. 14




*Tungsten carbide inserts with
ground pyramidal profile:
firm and positive grip of needle.*

| Needle holders with tungsten carbide inserts

ES 326
J18RT - Jacobson
Straight
cm. 18



ES 327
J18CT - Jacobson
Curved
cm. 18



ES 328
MC15TC - Micro-Vascular
cm. 15



ES 329
MC18TC - Micro-Vascular
cm. 18



ES 330
CV14RT - Castroviejo
Straight
cm. 14



ES 331
CV14CT - Castroviejo
Curved
cm. 14



ES 332
CR14RT - Castroviejo
Straight
cm. 14



ES 333
CR14CT - Castroviejo
Curved
cm. 14

| Hemostats



ES 334
HMR - Mosquito
Straight
cm. 12



ES 335
HMC - Mosquito
Curved
cm. 12



ES 336
KR - Kocher
Straight
cm. 14




ES 337
KC - Kocher
Curved
cm. 14

| sterilizing forceps



ES 338
PS 20
cm. 20



| Tissue Pliers



ES 339
 # 41 - Adson
 cm. 12 - mm.1,5



ES 340
 # 41D - Adson Delicate
 cm. 12 - mm.1,2



ES 341
 # 42 - Adson 1x2
 cm. 12 - mm.1,5




ES 342
 # 42D - Adson Delicate 1x2
 cm. 12 - mm.1,2



ES 343
 # 41L - Adson
 cm. 15 - mm.1,5



ES 344
 # 42L - Adson 1x2
 cm. 15 - mm.1,5



Surgical tweezers

The wide range of tweezers, made of a special steel alloy that gives them a proper elasticity, sensitivity and optimal grip, minimizes the risk of excessive trauma in periodontal tissues.

*Detail of the retentive profile of
 tweezers 1x2.*

| Tissue Pliers



ES 345
 # 31 - Semken-Taylor
Straight
 cm. 12,5



ES 346
 # 32 - Semken-Taylor
Curved
 cm. 12,5



ES 347
 # 33 - Semken-Taylor 1x2
Straight
 cm. 12,5



ES 348
 # 31L - Semken-Taylor
Straight
 cm. 15



ES 349
 # 32L - Semken-Taylor
Curved
 cm. 15



ES 350
 # 34 - Semken-Taylor 1x2
Curved
 cm. 12,5



ES 351
 # 90 - De Bakey-Vascular
Straight
 cm. 16



ES 352
 # 90C - De Bakey-Vascular
Angulated
 cm. 16



ES 353
 # 62 - Taylor
Bayonet

| Periostomes



ES 354
 # PM
Handle



ES 355
 # PD1
Straight
 1,6 mm.



ES 356
 # PC1
Curved
 1,6 mm.



ES 357
 # PA1
Angulated
 1,6 mm.



ES 358
 # PD2
Straight
 2,5 mm.



ES 359
 # PC2
Curved
 2,5 mm.



ES 360
 # PA2
Angulated
 2,6 mm.

| syndesmotomes



ES 361
Handle



ES 362
 # SYB 0



ES 363
 # SYB 2



ES 364
 # SYB 3



ES 365
 # SYB 4



ES 366
 # SYB 5



ES 367
 # 2 - Chompret
 cm. 16

Extracting Forceps, Meissner

With only three pliers you can extract any tooth in the upper or lower arch, one the left or right side.

A new frontier has been achieved, Lascod has created the series of Meissner pliers, keeping in mind the importance of easy cleaning and sanitation. The smooth handle and, above all, the fact that they can be totally and easily dismantled ensures maximum hygiene even in those parts that are usually inaccessible.




Dimensions 1:1

Extracting Forceps, Meissner



ES 368
Meissner

For root fragment



ES 369
Meissner

Universal, medium

5...1	1...5
5...1	1...5

The special ergonomic design makes these pliers extremely handy. The inclination and the tips are designed to maximize the grip of any tooth in the two dental arches.



ES 370
Meissner

Universal, large

8...6	6...8
8...6	6...8

Meissner pliers.
Smooth handle.
May be dismantled.
Three models
to meet all
requirements
of extraction.
Exclusive design.



Dimensions 1:1

Extracting Forceps

*Upper incisors
and canines.*

3...1	1...3
-	-



*Lower incisors
and canines.*

-	-
3...1	1...3



*Upper
premolars.*

5...4	4...5
-	-



Dimensions 1:1

Extracting Forceps



ES 381
13

*Lower
premolars.*

-	-
5...4	4...5



ES 382
17

*Upper molars,
right.*

-	6...7
-	-



89



ES 383
18

*Upper molars,
left.*

7...6	-
-	-



90

Dimensions 1:1

Extracting Forceps

*Lower
molars.*

-	-
7...6	6...7



ES 384
73



ES 385
22



ES 386
321

*Upper third
molars.*

8	8
-	-



ES 387
67A

*Lower third
molars.*

-	-
8	8



ES 388
79A

Dimensions 1:1

| Bone Rongeurs



ES 389
 # B45 - Blumenthal 45°
 cm. 14,5



ES 390
 # B90 - Blumenthal 90°
 cm. 14,5



ES 391
 # MF - Mini Friedman
 cm. 15



ES 392
 # FL - Friedman
 cm. 16,5



ES 393
 # MF90 - Mini Friedman 90°
 cm. 14



ES 394
 # FS - Friedman
 cm. 14



ES 395
 # CH - Cohen
 - Tissue nipper
 cm. 10



ES 396
 BONE RONGEURS BEYER
 DOUBEL ACTION 18CM



ES 397
 BONE RONGEURS BEYER
 DOUBEL ACTION
 RUSKIN-MINI 15.5CM

| Extracting Forceps of Children - Klein Style



III...I	I...III
-	-

ES 398
2

Upper Incisors and Canines



V...IV	IV...V
-	-

ES 399
7

Upper Premolars



-	-
III...I	I...III

ES 400
33

Lower Incisors, Canines and Roots



-	-
V...IV	IV...V

ES 401
13

Lower Premolars



-	-
V...IV	IV...V

ES 402
22

Lower Molars



V...IV	IV...V
-	-

ES 403
51

Upper Roots



V...IV	IV...V
-	-

ES 404
52

Upper Molars

| Root Elevators



ES 405
 # 2B - Bein
 mm. 2



ES 406
 # 3B - Bein
 mm. 3



ES 407
 # 4B - Bein
 mm. 4



ES 408
 # 2BFT - Bein
Flat tip
 mm. 2



ES 409
 # 3BFT - Bein
Flat tip
 mm. 3



ES 410
 # 4BFT - Bein
Flat tip
 mm. 4



ES 411
 # F1 - Flohr



ES 412
 # F2 - Flohr



ES 413
 # F3 - Flohr



ES 414
 # 12A - Apical



ES 415
 # 13A - Apical



ES 416
 # 1C - Coupland
 mm. 2



ES 417
 # 2C - Coupland
 mm. 3



ES 418
 # 3C - Coupland
 mm. 4



| Root Elevators



ES 419
 # 20
 Curved tip
 mm. 2



ES 420
 # 30
 Curved tip
 mm. 3



ES 421
 # 40
 Curved tip
 mm. 4



ES 422
 # LL - Lindo Levien
 mm. 3



ES 423
 # 27CR - Cryer



ES 424
 # 28CR - Cryer



ES 425
 # 71MA - Miller-Apexo



ES 426
 # 72MA - Miller-Apexo



ES 427
 # 11B - Bernard



ES 428
 # 6B - Bernard



ES 429
 # 7B - Bernard

| Root Tip Pick/Teaser



ES 430
1HD - Heidbrink



ES 431
2HD - Heidbrink



ES 432
3HD - Heidbrink

root fragments removing pliers



ES 433
1 - Stieglitz



ES 434
2 - Stieglitz

| luxators



ES 435
L2ST
Straight
mm. 2



ES 436
L3ST
Straight
mm. 3



ES 437
L4ST
Straight
mm. 4



ES 438
L2C
Curved
mm. 2



ES 439
L3C
Curved
mm. 3



ES 440
L4C
Curved
mm. 4

| **Pliers**



ES 441

Round/concave plier

For wires up to 0.028" (0.71 mm)



ES 442

Tweed loop forming plier

For wires up to 0.060" (1.52 mm); 0.090" (2.28 mm); 0.120" (3.04 mm)



ES 443

Loop forming plier Tweed style

For round and rectangular wires. Replaceable cylindrical beak with three sizes: 0.045" (1.14 mm); 0.060" (1.52 mm); 0.075" (1.90 mm)



ES 444

Loop closing plier Nance style

For round and rectangular wires. Wire max 0.028" (0.71 mm)



ES 445

Rectangular arch/ ribbon arch plier

For square/rectangular section wire of 0.022" (0.56 mm); 0.028" (0.71 mm)



ES 446

Lingual arch plier

To form double or triple back bends in wires of 0.030" (0.76 mm); 0.036" (0.91 mm)



ES 447

Arch forming and contouring plier

For wires of 0.016" (0.41 mm); 0.018" (0.46 mm); 0.022" (0.56 mm)



ES 448

- Arch forming plier, without groves,



ES 449

Accessories for MIH504



ES 450

Accessories MIH504 - Accessories for MIH504



| Pliers



ES 451
*Three
medium jaw plier Aderer bending*
For wires up to 0.036"
(0.91 mm)



ES 452
Hollow-Chop plier
To form and contour arches, designed for
NICKEL TITANIUM



ES 453
Rounded jaw wire bending pliers, short beak
To bend wires up to 0.029" (0.74 mm)



ES 454
Rounded jaw wire bending pliers, long beak
To bend wires up to 0.020" (0.51 mm)



ES 455
*Bird beak
Angle plier*
To bend round wires up to 0.045" (1.14 mm)



ES 456
Bird beak long plier
For small
loop forming. For light wires of 0.016" (0.41mm); 0.018"
(0.46 mm)



ES 457
Jarabak plier
To curve thin wires and loop form-
ing; 3 flute up to 0.020" (0.51 mm)



ES 458
Adams plier
Designed for bending wires up to 0.036" (0.91 mm)




ES 459
Laboratory bird beak long plier
Designed for loop forming and
bending wires up to 0.018" (0.46 mm)

| Pliers



ES 460
Steel and ceramic D.B. bracket
removing plier
Replaceable ends



ES 461
Direct bond bracket removing angled plier.



ES 462
Posterior band removing long plier
Plastic pad



ES 463
Posterior band removing short plier
Plastic pad



ES 464
Band seating plier
For positioning of
D.B brackets



ES 465
How plier-straight
For modeling and positioning of wires and
bands, with serrated tips.



ES 466
How plier-curved
For modeling and positioning of wires and
bands, with serrated tips.



ES 467
Accessories for MIH530



ES 468
Accessories for ZLH530



ES 470
Utility plier
Weingart TC
Tungsten
carbide insert.



ES 469
Accessory for MIH532 and MIH533

| Cutters/distal and Cutters With Tungsten Carbide Inserts



ES 471
 Hard wire
 cutter TC, straight.
 Hard wire, max 0.022" x 0.028" (0.56 mm x
 0.71 mm)



ES 472
 Mini
 pin and ligature cutter TC, straight.
 Max 0,015" (0,38 mm).



ES 473
 - Mini pin and ligature cutter TC, angled 15°.
 Max 0,015" (0,38 mm).



ES 474
 filo - Mini distal and cutter TC with safety hold.
 Max 0,022" x 0,028" (0,56 mm x 0,71 mm).



ES 475
 - Distal and cutter TC with safety hold.
 Max 0,022" x 0,028" (0,56 mm x 0,71 mm).



ES 476
 Flush cut distal and
 cutter TC.
 Max 0,022" x 0,028" (0,56 mm x 0,71 mm).



ES 477
 Mini distal and cutter TC
 with safety hold, long handle.
 Max 0,022" x 0,028" (0,56 mm x 0,71 mm).

| D.B. Brackets Holder



ES 478
Straight



ES 479
Angulated

| orthodontic ligature pliers



ES 480
MD14
Mathieu style
cm. 14



ES 481
OLP
Smaha style
cm. 14



| Utility Instruments



ES 482
*Ligature
tying plier*
Flared-tip metal plier for tying ligatures



ES 483
*Force
modular separator*
For positioning separating elastics




ES 484
*# AS
Band setting and
removing instrument*



ES 485
*# DC
Ligature instrument, curved.*



ES 486
*# DR
Ligature instrument, straight.*



ES 487
*# AL
Elastic setting instrument*



ES 488
Band pusher-Mershon
cm. 14

| Wax Carvers



ES 489
HL - Hylin



ES 490
HL-LC - Hylin-Lecron



ES 491
HL-C - Hylin-Ceramic



ES 492
HL-SP - Hylin-Spoon



ES 493
LC 1 - Lecron



ES 494
ZH 1 - Zahle



ES 495
ZH 2 - Zahle



ES 496
ZH 3 - Zahle



ES 497
GR - Gritman



ES 498
BA - Beale

| P.K. Thomas Waxing Instruments

- Aluminium handle Ø 8 mm.



ES 489
PKT 1V



ES 490
PKT 2G




ES 491
PKT 3V



ES 492
PKT 4R



ES 493
PKT 5A



PK THOMAS spatulas:

The complete set of PK THOMAS spatulas used for waxing is offered with a colored handle (green, gold, blue, red, steel) made of anodized aluminum to facilitate identification in the dental laboratory. The aluminum handle with a larger diameter (8 mm) makes the instrument lightweight, sensitive and precise, all important features to reproduce the correct anatomical morphology of natural teeth.

| Ceramic Sculpturing Instruments

- Aluminium handle
- 2 tips included



ES 494
1 G



ES 495
2 R



ES 496
3 B



ES 497
GR 1G - Greenstein



ES 498
GR 2R - Greenstein



ES 499
GR 3B - Greenstein



Superflexible tips, thickness 1/10



ES 500
1 A SF



ES 501
2 A SF



ES 502
3 A SF

Flexible tips, thickness 2/10



ES 503
1 BF



ES 504
2 BF



ES 505
3 BF

| Wax Knives



ES 506
1 F - Fahnenstock
cm. 12,5



ES 507
2 F - Fahnenstock
cm. 17



ES 508
1 L - Lessman
cm. 12,5



ES 509
2 L - Lessman
cm. 17

| plaster knives



ES 510
PK 1
cm. 17



ES 511
GR - Gritman
cm. 17



ES 512
3

| **Mixing Spatulas**



ES 513
1 SA
Rigid
cm. 18,5



ES 514
3 SA
Flexible
cm. 21,5



ES 515
1 SS
Flexible
cm. 18,5



ES 516
2 SS
Flexible
cm. 20,5

| **Crown Scissors**



ES 517
1 L 100
Retta - Straight
cm. 10

ES 518
2 L 100
Curva - Curved
cm. 10



| **gold foil scissors**



ES 519
1 G
Straight
cm. 12

ES 520
2G
Curved
cm. 12



| **scalpel handles**



ES 521
4



ES 522
8

| **alcohol lamp**



ES 523
LA



ES 524
Spare wicks for alcohol lamp

| Laboratory Pliers



ES 526
U 18 - Universal
cm. 18



ES 527
U25 - Universal
cm. 25



ES 525
U - Universal
cm. 15
For wires up to 0.036" (0.91 mm)



ES 528
U3 - Universal
cm. 16
For wires up to 0.036" (0.91 mm)



ES 529
PCP
Plaster cutting plier
cm. 19

| Mallets - Calipers



ES 531
MP
Plate tip
cm. 22



ES 532
BO - Boley
0 - 100 mm.



ES 530
IWW - Iwanson
For wax



ES 533
IWM - Iwanson
For metal

| Insulated Soldering Tweezers



ES 534
1
Straight
cm. 16,5



ES 535
2
Curved
cm. 16,5

| crucible tongs



ES 536
P 51
cm. 51



ES 537
P 43
cm. 43

Impression Trays

The absence of solder and retention edges facilitate the removal of residual material which is a breeding ground for blood, spores and bacteria in the mouth.



EXPLOIT has been designed bearing in mind the oral anatomy and an optimal distribution of the impression material, thus offering greater comfort to the patient.



Perforated and non-perforated EXPLOIT trays are compatible with all impression materials and are available with a micro-sand blasted satin finish for even more retention.

Impression Trays



XS

S

M

L



XL

SHAPES	Width mm.	Length mm.	Perforated	Non-perforated
UPPER				
XS	62	50	EXF112	EXL312
S	66	53	EXF113	EXL313
M	69	57	EXF114	EXL314
L	72	60	EXF115	EXL315
XL	80	63	EXF115	EXL316

	mm.	mm.		
LOWER				
XS	60	49	EXF212	EXL412
S	72	53	EXF213	EXL413
M	75	57	EXF214	EXL414
L	78	60	EXF215	EXL415
XL	80	63	EXF216	EXL416

"We manufacture custom design impression trays, if you are looking for the impression trays that are not listed here, you can let us know their image, reference number of any catalog, or drawing. We will manufacture for you cost-effectively with top quality standard."

Precision, Passion, Progress



For Inquiries and orders, reach us at:

Phone: +1 469-463-5877

Email: info@eagersurgical.com

www.eagersurgical.com