



SWIFTMED

Manufacturer for Injection Devices



SWIFTMED



SWIFTPEN I-PRO

SPRING MECHANISM (INWARD ROTATION)

Usage type	Multidose use. The injection pen shall be discarded once the cartridge is empty.
Compatible Cartridge	1.5mL, 3mL (ISO 11608-3, Others)
Compatible Needle	ISO 11608-2, Others
Minimum Injection Dose	0.01mL
Maximum Injection Dose	0.8mL
Injection Increment	0.01mL, 1U
Injection Accuracy	±5% (ISO 11608-1)
Pen Type	Plastic Disposable Multidose Single Cartridge Pen
Medicine Compatibility	This Pen can be fully customized to meet variety of medicine dosing requirements such as GLP-1s, Hormones, Biologicals, etc.
Other	Pen can be fully customized for different Colors, Imprints, Dosing marks, etc.

SWIFTPEN O-PRO

DIAL MECHANISM (OUTWARD ROTATION)

Usage type	Multidose use. Disposable: The injection pen shall be discarded once the cartridge is empty. Reusable: The injection pen can be reloaded with a new cartridge.
Compatible Cartridge	1.5mL, 3mL (ISO 11608-3, Others)
Compatible Needle	ISO 11608-2, Others
Minimum Injection Dose	0.01mL
Maximum Injection Dose	0.8mL
Injection Increment	0.01mL, 1U
Injection Accuracy	±5% (ISO 11608-1)
Pen Type	Plastic Disposable/Reusable Multidose Pen
Medicine Compatibility	This Pen can be fully customized to meet variety of medicine dosing requirements such as GLP-1s, Hormones, Biologicals, etc.
Other	Pen can be fully customized for different Colors, Imprints, Dosing marks, etc.

SWIFTPEN AUTO

SINGLE DOSE AUTOINJECTOR

Usage type	Single dose use. The injection pen shall be discarded after first use. Easy 3 step injection procedure: Open Safety Cap, Inject, Discard.
Compatible PFS	1mL, 2.25mL (ISO 11608-3, Others)
Viscosity	≤20 cP
Needle Diameter	27G to 32G
Pen Type	Plastic Disposable Single Dose Pen
Medicine Compatibility	This Pen can be fully customized to meet variety of medicine dosing requirements such as GLP-1s, Hormones, Biologicals, etc.
Other	Pen can be fully customized for different Colors, Imprints, etc.

SWIFTMED MEDICAL DEVICES PRIVATE LIMITED, Changodar, Ahmedabad, Gujarat (India)

Email: corporate@swiftmed.ca