

WIDEX MOMENT APP

THE PERFECTION OF NATURAL SOUND IS WHEN I FEEL LIKE I'M IN CONTROL

With Widex Moment™ hearing aids, you have the option to fully control your hearing aid using the Widex Moment app on your Apple or Android smartphone. It's never been easier to fine-tune your hearing, so you get exactly the sound you prefer - in exactly the moments that are important to you.

FEATURES AT A GLANCE



Personalise your sound

- Equalizer
- Sound Mixer
- Personal programs
- Made for you
- Made by you
- Wellbeing and relaxation programs
- Zen programs



Remote control

Use the app as a remote control for your hearing aids

- Hearing program selection
- Volume control
- Mute / unmute hearing aids
- Directional focus
- Right / left volume adjustments



Audio streaming


Stream audio directly from your iPhone and many models of Androids phones to your hearing aids (similar to earbuds)

- Control and adjust the sound to your liking
- Control and adjust audio streamed from TV Play device



Monitor your hearing aids

- View detailed battery status
- View connection status to app
- Find my hearing aids
- Turn on/off Bluetooth in hearing aids (flightmode)



Help on app and hearing aids

- Help section
- FAQs on handling hearing aids
- How-to videos



My Guide

- Learn to use and wear hearing aids
- Troubleshooting
- Courses
- Hearing aids reminder
- Find my hearing aids




Hearing aid update

- Firmware update
- How -to



Appearance

- Switch between dark and light mode



Pairing and connection

- Connecting guide
- In-app pairing

MYSOUND SOUND MADE FOR YOU

What if a hearing aid could give you access to the listening experiences of thousands of real users just like you, so they could give you a personally recommended sound in an instant?

With Widex MySound™, you can. MySound combines real human experiences with the immense power of artificial intelligence and condenses it into a simple and easy-to-use smartphone app.

WIDEX TV PLAY

GREAT SOUND FROM YOUR TV NEVER LOOKED SO GOOD



TV PLAY WITH NEW DESIGN

Enjoy stable, high-quality sound streamed directly from your TV to your hearing aids in style.

The upgraded Widex TV Play features the same elegant design accents as the rest of the Moment range, and the cooperation with Danish textile brand Kvadrat gives TV Play an unprecedented touch of tactility to an already award-winning design. Kvadrat fabrics have been used in countless pieces of classic Danish design and have sustainable principles at their core.

BIG TV MOMENTS DESERVE THE BEST SOUND

Whether you want to get updated on the latest TV news or just relax with a movie, TV Play helps you hear the sound from the TV better.

TV Play lets you fully enjoy what you want to see and hear through your hearing aids. With the dual antenna technology, TV Play offers maximum streaming stability, so you will never miss out on important TV moments.





HEAR IT YOUR WAY

Would you like to control the sound through your smartphone? Or do you just need a remote control?

With TV Play you can change programs and volume on your hearing aid, with an [RC-DEX control](#), or via an app on your smartphone. With the smartphone app you can also balance TV volume with ambient sounds or conversation. It's up to you what you want to hear.

WIDEX SOUND ASSIST

MAKES IT EVEN EASIER TO HEAR THE SOUNDS THAT MATTER TO YOU

Widex Sound Assist™ connects to your hearing aids to help you hear your life better.



THE BENEFITS OF WIDEX SOUND ASSIST

- **Partner microphone mode** – One-to-one conversations with friends and family sound clear wherever you are.
- **Table microphone mode** – Hear important meetings, lively dinners, and group conversations more easily.
- **Hands-free phone calls** – Great sounding calls on the go.
- **Streaming from any Bluetooth device** – Send music, film, or TV audio to your hearing aids.
- **Remote control for Widex Moment Bluetooth hearing aids** – Control your sound, easily.
- **Telecoil mode** – Telecoil functionality even if your hearing aids don't have telecoil.

WIDEX SOUND ASSIST CAN HELP IN ALMOST ALL SITUATIONS

Whether it's a one-to-one conversation, a lively dinner, or an important call, Sound Assist can help you feel more connected to the sounds that matter most to you.



A LOT OF POWER IN A TINY PACKAGE

- Battery runtime: 10 hrs (mixed use)
- Full charge: 3 hrs
- Bluetooth protocol: version 4.2
- Li-ion battery
- Micro USB power supply
- Small size: 43x43x16 mm
- Lightweight: 32 grams



WIDEX RC-DEX (remote control)

YOUR REMOTE CONTROL

RC-DEX is designed for basic wireless remote control of your hearing aids. It offers easy and discreet adjustment of your hearing aid volume and switching between programs.



SIMPLE AND INTUITIVE

RC-DEX is a compact and user-friendly hearing aid remote control. Its simple and intuitive design gives you easy control of basic hearing aid features like volume control and switching programs – also when you're watching TV with WIDEX TV PLAY™.

- **Stylish and user-friendly** - Looks good and is easy to use.
- **Discreet** - You won't need to touch your hearing aids to adjust them.
- **Function Control** of hearing aid functions like volume and programs.
- **Feature activation** - Verbal or tone cues confirm your selection.
- **Small and intuitive** design that will fit on your keychain.
- **Lock switch** - Prevents accidental activation.



SMALL ENOUGH TO FIT ON A KEYCHAIN

RC-DEX is intelligently laid out with three simple keys and is small enough to fit on a keychain. That means that you can quickly and easily start using it. And as features are activated, verbal or tone cues in your hearing aids let you know where you are. An LED indicator on the remote control also lights up when it is activated.

WIDEX TV-DEX

WATCH TV WITH NO ECHO

Enjoy distortion-free wireless transmission of TV-sound directly to your hearing aids.



EXPERIENCE SOUND EXACTLY AS IT WAS INTENDED

The TV-DEX is your user-friendly wireless assistive listening device specifically designed for enjoying TV and audio. The main advantage of the TV-DEX is its real-time high-quality stereo sound. Experience TV or audio sound exactly as it was intended - without annoying distortion or more importantly - echoes.

NO BACKGROUND NOISE

Another unique feature of the TV-DEX is the 'Room Off' function. It allows you to temporarily switch off the hearing aid's microphone and hear the TV sound only. This way you can enjoy your favourite TV programmes without unnecessary background noise. Naturally, Room off can easily be switched on and off.