



ZydroMed
INNOVATING HEALTHCARE
SINCE 2021


Pressure Relief

PATIENT POSITIONER


www.zydromed.com

Donut/ Open Head Ring

Multifunctional positioner



01 piece/pack



Item ¹	Size	Dimensions (mm) X* Y*Z
OHR-N	Neonatal	90 x 45 x 20
OHR-N-F	Neonatal	
OHR-M	Medium	150 x 80 x 30
OHR-M-F	Medium	
OHR-L	Large	200 x 100 x 40
OHR-L-F	Large	

¹Item code with suffix 'F' indicate Gel encapsulated in Film


- Supports patient's head during prolonged surgeries (air and gel cushioning)
- Protect skull, neck and ears in multiple positions by conforming to head's contour

Closed Head Ring

Multifunctional positioner



01 piece/pack



Item ¹	Size	Dimensions (mm) X* Y*Z
CHR-N	Neonatal	90 x 45 x 20
CHR-N-F	Neonatal	
CHR-M	Medium	140 x 70 x 35
CHR-M-F	Medium	
CHR-L	Large	200 x 100 x 45
CHR-L-F	Large	

¹Item code with suffix 'F' indicate Gel encapsulated in Film


- Supports patient's head with complete coverage to back of head during prolonged surgeries
- Protect head, neck and ears in various positions by complete pressure redistribution

C-Shape/U-Shape Head Ring

Multifunctional positioner



01 piece/pack




Item ¹	Size	Dimensions (mm) X* Y*Z
UHR-M	Medium	140 x 60 x 45
UHR-M-F	Medium	
UHR-L	Large	200 x 100 x 60
UHR-L-F	Large	

- Supports patient's face during surgery allows access to intubation due to its contoured features.
- No pressure on eyes especially in prone position due to ergonomic design and elastomeric gel technology.




Supine/Contoured Head Rest

Unipositional Protector



01 piece/pack



Item	Size	Dimensions (mm) L*W*H [d*h]
SHR-SD	Standard	175 x 180 x 70 [145x35]
SHR-SD-F	Standard	


- Relieves skull compression from head / calf in supine position (like eye surgeries) by providing cross-sectional stability
- Contoured to match the anatomical shape of head by pressure redistribution in all directions

Ophthalmic Cradle Head Rest

Uni-functional positioner



01 piece/pack



Item	Size	Dimensions (mm) L*W*H [l*w*h*s]
O-CHR	Standard	240 x 200 x 40 [170 x 140 x 20 x 50]
O-CHR-F	Standard	


- Contoured cradling and extra support/ stability to patient's head in supine or other positions preventing risk of alopecia in ophthalmic long tenured surgeries.
- Extensive pressure redistribution for complete head rest.

Horse Shoe Ring -Eye Pressure Relief

Uni-functional protector



01 piece/pack




Item	Size	Dimensions (mm) D*H [l1*w1] [l2*w2]
HORSE-SD	Standard	200 x 75 x [50 x 140] [90 x 75]
HORSE-SD-F	Standard	

- Supports patient's Face during surgery allow access to intubation due to its contoured features
- No pressure on eyes especially in prone position due to ergonomic design and elastomeric gel technology



Prone Head Rest

Prone position protector




Type-I
01 piece/pack



Type-II
01 piece/pack



Type-III
01 piece/pack



Template
(Indicative)

Item	Size	Dimensions (mm) L*W*H
PHR-I-SD	Standard	260 x 215 x 140
PHR-I-SD-F	Standard	
PHR-II-M	Adolscnt	230 x 190 x 130
PHR-II-M-F	Adolscnt	
PHR-II-L	Adult	270 x 230 x 140
PHR-II-L-F	Adult	
PHR-III-SD	Standard	260 x 240 x 145
PHR-III-SD-F	Standard	


- Supports patient's head in prone position during prolonged surgeries
- Protect head, neck and ears in prone position by complete pressure redistribution and ergonomic design allows ease of intubation.
- PHR - III provides see through area at its bottom

Gel Base for Metallic Neurohead Frame

Critical functional positioner



01 set of 02 half pieces/pack
(without Neurohead Frame Device)



Item	Size	Dimensions (mm) L* W*H [l * w]
NEURO-L	Large	200 x 55 x 65 [145 x 85]
NEURO-L-F	Large	
NEURO-M	Medium	175 x 50 x 60 [125 x 80]
NEURO-M-F	Medium	

- Provide cushioning interface to metallic neuro-head frame during prolonged brain surgeries.
- Supports patient's facial contours with complete pressure redistribution.



Organic



Single Phase



Non Film



Elastomeric




No Leak




Non Allergenic

Dome Flat Bottom Chest Roll

Multi-utility bolster



01 piece/pack




Item	Size	Dimensions (mm) L*W*H
DOME-S	Small	350 x 145 x 75
DOME-S-F	Small	
DOME-M	Medium	405 x 150 x 100
DOME-M-F	Medium	
DOME-L	Large	510 x 155 x 125
DOME-L-F	Large	


- Provides excellent support and stability to chest areas & other body parts by complete pressure redistribution
- Elevate the chest area and also useful in multiple sites

Flat Bottom Chest Roll

Standard positioner



01 piece/pack




Item	Size	Dimensions (mm) L*W*H
FBCR-M	Medium	405 x 140 x 100
FBCR-M-F	Medium	
FBCR-L	Large	510 x 150 x 100
FBCR-L-F	Large	


- Provides excellent support and stability to chest areas during prolonged surgeries
- Pressure redistribution in chest area preventing bed sores

Universal Positioner

All-site positioner



01 piece/pack




Item	Size	Dimensions (mm) L* W*H
UP-SD	Standard	420 x 145 x 60
UP-SD-F	Standard	

- It has unique and smooth curve structural line allows it to position at every body part.
- It has single phase continuous surface with complete pressure redistribution & support even twisted body parts




Patient Positioner Device

Multi-site positioner



01 piece/pack




Item	Size	Dimensions (mm) L* W*H [(l*w)]
PPD-M	Medium	420 x 130 x 55 [200 x 95]
PPD-M-F	Medium	


- Supports patient's twisted body parts during prolonged surgeries
- Protect body parts in various uncomfortable positions against torque by complete pressure redistribution

Wedge Positioner (15°)

Body Elevator



01 piece/pack




Item	Size	Dimensions (mm) L* W*H
WP-15	15°	280 x 280 x 75
WP-15-F	15°	


- It is a positioning device that rests body part at an angle of 15° slant
- Elevate body part while distribute pressure equally in all directions

Lateral Positioner

Side Positioner



01 piece/pack



Item	Size	Dimensions (mm) L* W*H [(h*d)]
LP-SD	Standard	710 x 480 x 150 (60 x 25)
LP-SD-F	Standard	


- It provides platform that rests the body part like upper part of arm while the patient needs to extend its arm in contoured arm board
- Elevates upper arm provides similar positioning in various positions



Prone Gel Bed Upper Body Protector



01 piece/pack




Item	Size	Dimensions (mm) L* W*H (t _i *t _w)(i*g*n)
PGB-SD	Standard	780 X.... X 125 (600 X 530) (360 X 280 X 210)
PGB-SD-F	Standard	

- Provides extreme cushioning to upper abdominal area while using it with Prone Head Rest
- Cut-outs for positioning Head Rests & Arm Boards to stabilise upper body in prone position

Elevation Pillow Upper Body Elevator



01 piece/pack



Item	Size	Dimensions (mm) L*W*H1*H2(RP*RP _w)
EP-SD	Standard	750 X 510 X 137 X 37 (175 X 180)
EP-SD-F	Standard	

- Elevates upper/lower body part based on positioning requirement
- Slant neck support aligns body with optimal blood circulation in stabilised body parts



Organic



Single Phase



Non Film



Elastomeric




No Leak




Non Allergenic

Contoured Arm Board Pad Multifunctional extensor



01 piece/pack



Item	Size	Dimensions (mm) L*W*H
CONT-ABP	Standard	500 x 130 x 35
CONT-ABP-F	Standard	


- Anatomically shaped armboard pad to provide support by pressure redistribution
- Slightly grooved to hold the arm with release of pressure on either direction

Sacral Protector

Specific-site padded device



01 piece/pack




Item	Size	Dimensions (mm) L* W*H [b*d]
SACP-SD	Standard	445 x 440 x 30 [130 x 100]
SACP-SD-F	Standard	


- Provides protection to unprotected sacral area in lithotomy and supine position, by relieving pressure and transverting it to peripheral areas

Hip Pad Position Device

Critical functional positioner



01 piece/pack



Item	Size	Dimensions (mm) L*W*H
HPPD-S	Small	400 x 400 x 12
HPPD-S-F	Small	
HPPD-L	Large	520 x 520 x 12
HPPD-L-F	Large	


- Supports patient's hip during prolonged surgeries
- Hip protection with centrifugal equidistant pressure distribution

Flat Arm Board Pad

Multi-utility Extensor



01 piece/pack




Item	Size	Dimensions (mm) L* W*H
FABP-SD	Standard	600 x 125 x 12
FABP-SD-F	Standard	

- FABP provides stabilized resting state to the extended arm during prolonged surgeries
- Standard OT Table Extensor useful in various body positions




Ulnar/Brachial Nerve Protector

Bi-functional protector



01 piece/pack




Item	Size	Dimensions (mm) L* W*H
UBNP-M	Medium	350 x 150 x 10
UBNP-M-F	Medium	
UBNP-L	Large	450 x 200 x 10
UBNP-L-F	Large	


- Providing protection to the sensitive areas of elbows and arms during prolonged surgeries
- Comforting sensitive nerve areas by extreme cushioning

Contoured Knee Pad

Multifunctional device



01 piece/pack




Item	Size	Dimensions (mm) L*W*H
CKP-SD	Standard	500 x 130 x 35
CKP-SD-F	Standard	


- Anatomically shaped knee pad to provide added support and stability by pressure redistribution
- Standard OT Table Extensor useful in multiple positions in lower body parts .

Knee Crutch Pad

Knee support device



01 piece/pack



Item	Size	Dimensions (mm) L* W*H
KCP-S	Small	230 x 300 x 12
KCP-S-F	Small	
KCP-M	Medium	280 x 400 x 12
KCP-M-F	Medium	
KCP-L	Large	300 x 500 x 12
KCP-L-F	Large	

- Support to base of knee and leg, stabilize the position while knee and leg is aligned at multiple angles, resting during prolonged surgeries
- Protect nerve endings and complete pressure redistribution allows blood flow




Contoured Leg Rest

Leg resting positioner



01 pair/pack




Item	Size	Dimensions (mm) L* W*H
CLP-SD	Standard	500 x 170 x 40
CLP-SD-F	Standard	


- Ergonomic designed support to leg and heel area protects popliteal nerve by aligning with knee groove till anterior end of ankle region.
- Standard OT Table Resting device useful in providing stability to multiple body positions.

Contoured Calf Rest

Calf Resting Device



01 pair/pack



Item	Size	Dimensions (mm) L* W*H
CCR-SD	Standard	300 x 130 x 35
CCR-SD-F	Standard	


- It has grooved structure with slight curved edge as resting surface for calf area.
- Single phase continuous surface with enhanced pressure distribution on either side supporting calf as padded support device

Crutch Stirrup Pad

Knee & leg resting device



01 pair/pack
(without Stirrup Device)



Item	Size	Dimensions (mm) L*W*H
CSP-S	Small	320 x 220 x 10
CSP-S-F	Small	320 x 220 x 10
CSP-M	Medium	400 x 250 x 10
CSP-M-F	Medium	400 x 250 x 10
CSP-L	Large	520 x 300 x 10
CSP-L-F	Large	520 x 300 x 10

- Designed to fit both Lloyd Davies & various types of stirrups / leg-holders to relieve pressure from nerves and allows blood flow especially during lithotomy procedures.
- Protect against frictional impact, shear force and nerve injury.



Organic



Single Phase



Non Film



Elastomeric




No Leak




Non Allergenic

Heel Protector

Lower Extremity Protector



01 pair/pack



Item	Size	Dimensions (mm) L* W*H [(b*d)(l1*w1*w2)]
HPR-SD	Standard	175 x 180 x 40 [20 x 60] [170 x 60 x 10]
HPR-SD-F	Standard	


- It looks like heart shaped with elevated areas to protect heel on either side. Inner curvature provide padded platform for the heels during prolonged surgeries
- Single phase continuous surface, ergonomically designed to relieve pressure, enhanced blood flow to extremities with padded support till ankle region

Wilson Frame Pad

Spine aligner -cum-cushioner



01 pair/pack
(without Wilson Frame Device)



Item	Size	Dimensions (mm) L*W*H
WFP-M	Medium	660 x 100 x 12
WFP-M-F	Medium	
WFP-L	Large	750 x 255 x 12
WFP-L-F	Large	750 x 255 x 12

- Supports patient's spine by aligning it on wilson frame and WFP provide extreme elastomeric cushioning with added protection from shear & pressure. WFP is provided alone and not with Wilson frame.
- Prefixed on frame with hoop & loop attachment used frequently in neurological & orthopaedic prolonged surgical procedures



Organic



Single Phase



Non Film



Elastomeric



No Leak




Non Allergenic

Heel Cup (Without Strap)

Lower Extremity Resting Device



01 pair/pack




Item	Size	Dimensions (mm) L*W*H [l*w*a]
HCWS-SD	Standard	130x 180 x 65 [60 x 80 x 40]
HCWS-SD-F	Standard	


- Designed to provide maximum protection to heel and ankle area during prolonged surgeries.
- Elastomeric rebound at body extremity by extreme resting state and adequate blood flow preventing foot ulceration.

Heel Support Pad

Lower Extremity Protector



01 pair/pack



Item	Size	Dimensions (mm) L* W*H [*h]
HSP-M	Medium	150 x 80 x 55 [50 x 30]
HSP-M-F	Medium	
HSP-L	Large	200 x 100 x 75 [75 x 35]
HSP-L-F	Large	


- Relieves pressure & enhanced blood flow to lower extremity by elastomeric ankle support and air cushioning to heels.
- Especially used in foot sole surgeries by exposing foot area.

Boot Stirrup Pad

Foot & ankle resting device



01 pair/pack




Item	Size	Dimensions (mm) L*W*H
BSP-SD	Standard	600 x 150 x 10
BSP-SD	Standard	600 x 150 x 10


- Designed to protect patient's foot & ankle area in universal stirrups.
- Pads are used in lithotomy position for gynecological & urological procedures; colon/rectal cases.

Face/ Fore Head Protector

Uni-site protector



01 piece/pack



Item	Size	Dimensions (mm) L* W*H
FHP-SD	Standard	165 x 80 x 8
FHP-SD-F	Standard	

- Protect the patient's forehead / facial contours from pressure caused by endotracheal tubes
- Universal protector and it can be used at any similar body site



Organic



Single Phase



Non Film



Elastomeric




No Leak



Non Allergenic

Operation Theatre Table Pads



OT Table Pad (Flat Type) - 3Pc set
Whole body cushioning prevents bed sores

Item	Operation Theatre Table Pads*	Dimensions (mm) L*W*H
UOTP-F	Universal Operating Table Pad, Flat Type	500 X 500 X 12
UOTP-F	Universal Operating Table Pad, Flat Type	
HLOTP-F	Half Length Operating Table Pad, Flat Type	900 X 500 X 12
HLOTP-F	Half Length Operating Table Pad, Flat Type	
¾ LOTP-F	3/4 Length Operating Table Pad, Flat Type	1150 X 500 X 12
¾ LOTP-F	3/4 Length Operating Table Pad, Flat Type	
¾ LOTP-WCO	3/4 Length Operating Table Pad, Flat Type, With Cut-Out	1150 X 500 X 12
¾ LOTP-WCO	3/4 Length Operating Table Pad, Flat Type, With Cut-Out	
PTP-F	Perineal Table Pad, Flat Type, With Cut-Out	500 X 500 X 12
PTP-F	Perineal Table Pad, Flat Type, With Cut-Out	
OTTP-3-F	Operation Theatre Table Pad 3pc. set Flat Type	
	Head Section	225 X 500 X 12
	Head Section-F	
	Torso Section	990 X 500 X 12
	Torso Section-F	
	Foot Section	500 X 500 X 12
	Foot Section-F	

* All OT Table Pads except OTTP-3-F are single unit packs

Elastomeric Cushioning
(~2 times stretchability)

Natural Skin Touch
(Feel it yourself)

Organic - Environmental Friendly
(non toxic)

Non leaky
(no release of any fluid)

Stick On
(Self attached Device)

Moisture Resistant
(Dip water test)

Resilience Factor
(~2000 pounds, 900Kgs)

Nature Colour Variation
(Blue, Pink & Green)

Dealer/ Distributor

H.R. ENTERPRISES

Call: +91 9318466233

E-mail: info@zydromed.com

Website: www.zydromed.com

Address: E 1/117. Sonia Vihar Delhi-110094