


VISTA[®] 300 ELECTRIC OPERATING TABLE



VISTA® 300

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External Size (LxWxH):	2070x550x680-1000mm
▪ Built-in Kidney Bridge:	0-120 mm
▪ Back Plate Up/Down:	75°/35°
▪ Head Plate Up/Down:	45°/90°
▪ Leg Plate Up/Down:	30°/90°
▪ Trendelenburg/Reverse Trendelenburg:	30°/30°
▪ Lateral Tilt Left/Right:	22°/22°
▪ Horizontal Sliding:	350 mm
▪ X-rays translucent:	YES
▪ C-arm compatible:	YES
▪ Multi-functional orthopedic frame:	OPTIONAL
▪ Supported Weight	250kg
▪ Battery	3 hours

Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Support with Cushion	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pcs
<input checked="" type="checkbox"/> Column Controller	1pcs

Remote Panel

All function controlled by hand controller, and can be locked.
Optional: Double control system on column



C-arm

Extra-long horizontal sliding, available for c-arm operation.




Kidney Bridge

Built-in kidney bridge, more convenient and nice appearance.



Wash Basin Hole

Five section table top including Head, Back, Shoulder, Hip and Leg, all section applicable for X ray use and under each section with X ray cassette layer.



Head Plate

Detachable double head plate and leg plate meet operation position needs and use Germany cylinder assist system.



Power Supply

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.



Shoulder Support with Cushion



Body Support with Cushion



Arm Support with Cushion



Leg Support with Cushion



Lowest Position of Table Top
680 mm



Highest Position of Table Top
1000 mm



Longitudinal Adjustment
0-350mm



Longitudinal Adjustment
0-350mm



Reversed Trendelenburg
0°-30°



Trendelenburg
0°-30°



Lateral Tilt Left
0°-22°



Lateral Tilt Right
0°-22°



Kidney Bridge Elevation
120 mm



Anorectum Operation
Position



E.N.T Operation Position



(Optional) Head Rest



(Optional) Lateral Arm
Support



(Optional) Waste Basin For
Gynecology



(Optional) Orthopedics Traction
Frame





 **CONSTANT
INNOVATION**



Design



Research



Manufacturing



Maintenance



Service

 **WORLDWIDE
BRANCHES**

United Kingdom
Russia
East Africa
Middle East



ISO13485

Visor Healthcare Solutions
North London Business Park, Building 5, Office 1
Oakleigh Rd South, New Southgate,
London N11 1GN, UK
www.visorhealthcare.co.uk