

# Software release notes



**Version:** 1.4.14

**Released:** 27th March 2024

**Download link:** [ [Link](#) ]

- More STL exports added.

**Version:** 1.4.06

**Released:** 18th Jan 2024

**Download link:** [ [Link](#) ]

- Faster and more reliable
  - Software opens 4x faster.
  - Scanner data transfer 16x faster AND more reliable.
    - Less wifi drop outs during scanning.
  - Calibration
    - Faster image capture.
    - Faster calibration processing.
    - Calibration progress bar now also on the LCD.
- UI simplifications.
- Text easier to read on scanner LCD.
- More STL exports added.
  
- Bug fixes
  - Software uses almost half the RAM memory compared to previous versions.
    - Easier on your laptops resources, and more stable.
  - Main window closes without freezing.
  - Occasional crash bugs fixed.

**Version:** 1.3.50 Initial Release

**Released:** 15<sup>th</sup> December 2023

**Download link:** [ [link](#) ]

Initial Release

- Calibration now lasts for 3 hours (not 1 hour)
- Foreign language wifi bug fixed.

## How to install new software

Download and run the installation file using the links in this document.

### Notes:

- The links will make a fresh install or update an existing version.
- Your scan data will not be overwritten.

