

ADULT AIRWAY MANAGEMENT TRAINER

Realistic practice is the key to developing proficiency in airway management skills. The Laerdal Airway Management Trainer realistically simulates an adult airway, and the lifelike upper torso and head can be used to demonstrate and practice intubation, ventilation, and suction techniques. In addition, it can be used to demonstrate upper airway bronchoscopy and bronchoscopy guided endotracheal intubation.

Bag-Valve Ventilation

- Provides visual inspection of lung expansion
- Provides auscultation of breath sounds
- Sellick Maneuver (Cricoid pressure)
- Laryngospasm



Intubation

- Tracheal (oral and nasal)
- Pharyngeal (oral and nasal)
- Esophageal
- Bronchial
- LMA insertion
- Endotracheal intubation
- Placement of supraglottic airways
- Correct tube placement
- Feedback for excessive laryngoscope pressure

- Bronchoscopic evaluation of tip position

Suction and Clearing

- Oral cavity
- Oro- or nasopharynx
- Oro- or nasotracheal, via endotracheal tube
- Gastric drainage

This hand-held teaching tool helps instructors explain the anatomy of the airway and show how the life size anatomy piece some medical devices interact with the airway in a visual manner. It allows for learners to familiarize themselves with the anatomy of adult airways in a more hands-on manner.

